## Multicalibrador CFAS Humano MonlabTest®

IVD

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## Human Multicalibrator CFAS MonlabTest®

IVD

### Serum calibrator for clinical chemistry assays

Only for professional *in vitro* diagnostic use.  
Store at 2-8°C.

#### PRODUCT CHARACTERISTICS

Human Multicalibrator CFAS is a human lyophilized serum. The concentrations and activities have been selected so as to ensure optimum calibration for use in manual and automatic analyzers.

#### REAGENTS

Human serum. Biological additives. Bacteriostatic agents. The concentration / activities of the components are lot-specific. The exact values and ranges valid for reagents are given informational purpose in the value sheet.

#### PRECAUTIONS

The individual units comprising lot or lots of components from human origin, have been tested for HBsAg and antibodies to HIV and found negative / non reactive by an FDA approved test. However, handle cautiously as potentially infectious.

#### CALIBRATION

Information available on request.

#### PREPARATION

- Reconstitute (→) with 3.0 mL of distilled water.
  - Mix thoroughly, avoiding foam forming.
  - Bring to room temperature for about 30 minutes before use.
- Improper handling and/or storage can affect results. Inaccurate reconstitution and errors in assay technique can cause erroneous results.

#### STORAGE AND STABILITY

The calibrator is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date or if there is visible evidence of microbial growth. Store calibrators tightly capped when not use.

- After reconstitution is stable for:
  - At 15°C to 25°C..... 8 hours
  - At 2°C to 8°C ..... 2 days
  - At -25°C to -15°C ..... 4 weeks
- Bilirubin (stored protected from light):
  - At 15°C to 25°C..... 6 - 3 hours (total-direct)
  - At 2°C to 8°C ..... 24 - 8 hours (total-direct)
  - At -25°C to -15°C ..... 2 weeks
- Acid phosphatase:
  - At 15°C to 25°C..... 4 hours
  - At 2°C to 8°C ..... 1 day
  - At -25°C to -15°C ..... 2 weeks

#### PRESENTATION

MO-165106

10 x 3 mL