Ferritina MonlabTest®

Turbidimetría Látex

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Ferritin MonlabTest®

Latex Turbidimetry

Quantitative determination of Ferritin

Only for professional *in vitro* diagnostic use. Store at 2-8°C.

INTENDED USE

Turbidimetric immunoassay for the quantitative determination of ferritin in human serum or plasma.

PRINCIPLE OF THE METHOD

Ferritin MonlabTest is a quantitative turbidimetric test for the measurement of ferritin in human serum or plasma. Latex particles coated with specific anti-human ferritin are agglutinated when mixed with samples containing ferritin. The agglutination causes an absorbance change, dependent upon the ferritin contents of the sample that can be quantified by comparison from a calibrator of known ferritin concentration.

CLINICAL SIGNIFICANCE

Serum ferritin concentration usually reflects body iron stores and is considered one of the most reliable indicators of iron status of patients. Whereas low serum concentrations of ferritin are always indicative of an iron deficiency, elevated concentrations can occur for variety of reasons. Thus, although elevated concentrations often indicate an excessive iron intake, they are also caused by liver disease, chronic inflammation and malignancies. Pregnant women, blood donors, hemodialysis patients, adolescents and children are groups particularly at risk.

REAGENTS

Diluent (R1)	Tris Buffer 20 mmol/L, pH 8.2. Preservative.
Latex (R2)	Latex particles coated with rabbit IgG anti-human ferritin, pH 8.2. Preservative.
Ferritin Calibrator	Calibrator. Ferritin concentration is stated on the vial.
Optional	MO-165056 Ferritin Control.

PRECAUTIONS

R1 and R2: H317 - May cause an allergic skin reaction. Contains 2-Methylisothiazol-3(2H)-one (Proclin 950). Follow the precautionary statements given in MSDS and label of the product. Components from human origin have been tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However, handle cautiously as potentially infectious.

PREPARATION

Ferritin Calibrator: Reconstitute (→) with 3.0 mL of distilled water. Mix gently and incubate at room temperature for about 10 minutes before testing.

CALIBRATION

Use Ferritin Calibrator MonlabTest (MO-165055). The sensitivity of the assay and the target value of the calibrator have been standardized against the 3rd International Standard of Ferritin (94/572, 2008 WHO). Recalibrate when control results are out of specified values; when using a different lot of reagent and when the instrument is adjusted. **Calibration curve:** Prepare the following dilutions of the Ferritin Calibrator using NaCl 9 g/L. To obtain the concentration of each dilution, multiply using the dilution factor shown in the next table:

Calibrator dilution	1	2	3	4	5	6
Calibrator Ferritin (µL)	--	25	50	100	200	400
NaCl 9 g/L (µL)	400	375	350	300	200	--
Dilution Factor	0	1/16	1/8	1/4	1/2	1.0

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Reagents should not be left inside the analyzer after use, they must be stored refrigerated at 2-8°C. Latex may sediment. Mix reagents gently before use. Do not use reagents over the expiration date. **Do not freeze.** Frozen Latex or Diluent could change the functionality of the test. **Reagent deterioration:** Presence of particles (R1, R2) and turbidity (R1).

ADDITIONAL EQUIPMENT

- Thermostatic bath at 37°C.
- Spectrophotometer or photometer thermostatable at 37°C with a 540 nm filter.

SAMPLES

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C. The samples with presence of fibrin should be centrifuged before testing. Do not use highly hemolyzed or lipemic samples.

PROCEDURE

1. Bring the reagents and the photometer (cuvette holder) to 37°C.
2. Assay conditions:
Wavelength: 540 nm (530-550)
Temperature: 37°C
Cuvette light path: 1 cm
3. Adjust the instrument to zero with distilled water.
4. Pipette into a cuvette:

Diluent (R1)	800 µL
Latex (R2)	200 µL
Calibrator or sample	90 µL

5. Mix and read the absorbance immediately (A₁) and after 5 minutes (A₂) of the sample addition.

MONLAB has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.

CALCULATIONS

Calculate the absorbance difference (A₂-A₁) of each point of the calibration curve and plot the values obtained against the Ferritin concentration of each calibrator dilution. Ferritin concentration in the sample is calculated by interpolation of its (A₂-A₁) in the calibration curve.

QUALITY CONTROL

Control Sera are recommended to monitor the performance of manual and automated assay procedures. It should be used the Ferritin Control MonlabTest (MO-165056).

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES

- Men: 30 – 220 µg/L
- Women: 20 – 110 µg/L

Each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Measuring range: Up to 600 µg/L. Samples with higher values should be diluted 1/5 in NaCl 9 g/L and retested. The upper linearity limit increases as the sample volume and the sensitivity decrease.

Detection limit: 5.04 µg/L.

Quantification limit: Values under 6.6 µg/L may give non-reproducible results.

Prozone effect: No prozone effect was detected at least up to 9000 µg/L.

Precision: According to the EP5-A2 standards (CLSI), the reagent has been tested for 20 days, measuring each level per duplicate twice a day (n=80):

Mean (µg/L)	Intra-assay (n= 80)			Total (n= 80)		
	33.4	114.5	289.8	33.4	114.5	289.8
SD	1.7	1.4	2.4	2.1	3.4	7.5
CV (%)	5.1	1.2	0.8	6.3	2.9	2.6

Method comparison: The reagent was compared to another commercially available Ferritin reagent by testing 144 samples (male and female), with concentrations between 6.97 and 730 µg/L. The coefficient of correlation (r) was 0.988, and the equation y = 0.96x + 1,15.

Performance characteristics depend on the analyzer used.

INTERFERENCES

Bilirubin (40 mg/dL), hemoglobin (5 g/L), and rheumatoid factor (750 UI/mL), do not interfere. Lipids (≥ 2,5 g/L) do interfere. Other substances may interfere⁵.

NOTES

Clinical diagnosis should not be based on findings of a single test result but should integrate both clinical and laboratory data.

BIBLIOGRAPHY

1. Knovich MA et al., Blood Rev. 2009 23(3):95-104.
2. Mazza J et al. Can Med Assoc J 1978; 119: 884-886
3. Rodríguez Perez J et al. Revista Clínica Española 1980; 156 (1): 39-43
4. Milman N et al. Eur J Haematol 1994; 53: 16-20.
5. Young DS. Effects of drugs on clinical laboratory test, 5th ed. AACCC Press, 1999.

PACKAGING

MO-165032 R1: 1 x 40 mL
R2: 1 x 10 mL
Ferritin CAL: 1 x 3 mL

SYMBOLS FOR IVD COMPONENTS AND REAGENTS

	Manufacturer		For <i>in vitro</i> diagnostic use only
	Don't re-use		Consult instructions for use
	Contains sufficient for <n> tests		Keep dry
	Catalogue Code		Temperature limitation
	Lot Number		Use by