Strep A MonlabTest®

DIM-100001 25 TESTS

One Step Streptococcus Group A Test Device

IVD

A rapid, one step test for the qualitative detection of Group A *Streptococcus* (GAS) from throat swab.

For professional *in vitro* diagnostic use only.

INTENDED USE

The Strep A MonlabTest® is a rapid chromatographic immunoassay for the qualitative detection of Group A *Streptococcus* from throat swab to aid in the diagnosis of GAS infection ("strep throat").

Strep A MonlabTest® used to detect Group A *Streptococcus* colonies from culture.

SYNTHESIS

Group A *Streptococcus* is a bacterium often found in the throat and on the skin. People may carry group A streptococci in the throat or on the skin and have no symptoms of illness. Most GAS infections are relatively mild illnesses such as "strep throat," or impetigo. On rare occasions, these bacteria can cause other severe and even life-threatening diseases.

Strep throat is an infection caused by group A *Streptococcus* bacteria, and it's very common among kids and teens. The symptoms of "strep throat" include fever, stomach pain, and red, swollen tonsils. Strep throat may produce mild or severe symptoms.

PRINCIPLE

Strep A MonlabTest® is a qualitative lateral flow immunoassay for the detection of Group A *Streptococcus* antigen in human throat swab samples. The membrane is pre-coated with monoclonal antibodies against Group A *Streptococcus* antigen on the test line region. During testing, the sample reacts with the particle coated with anti-Group A antibodies which was pre-dried on the test strip. The mixture moves upward on the membrane by capillary action. In the case of a positive result the specific antibodies present on the membrane will react with the mixture conjugate and generate coloured lines. A blue coloured band always appears in the control line and serves as verification that sufficient volume was added, that proper flow was obtained and as an internal control for the reagents.

PRECAUTIONS

- For professional *in vitro* diagnostic use only.
- Do not use after expiration date.
- The test should remain in the sealed pouch until use.
- Do not use the test if pouch is damaged.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, do not eat, drink or smoke in the area.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The test should be discarded in a proper biohazard container after testing.
- The presence of yellow lines in the results window (control and test line zone) that are visible before using the test are completely normal. That does not mean failure on test functionality.
- Excepting from the Diluent A and the Diluent B, all the components of the kit are for single use. Do not reuse.
- Diluent A (hazardous ingredient: Sodium Nitrite): Avoid contact with eyes and skin. Do not ingest or inhale. Harmful if swallowed. Hazard statement(s): H302.

STORAGE AND STABILITY

Store as packaged in the sealed pouch either at refrigerated or room temperature (2-30°C/36-86°F). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

MATERIALS PROVIDED

- 25 Tests
- Instructions for use
- 1 Diluent B (Sample diluent)
- 1 Diluent A (Sample diluent)
- 25 Sterile swabs
- 1 Strep A control +
- 25 Testing tubes or vials
- 25 Plastic pipettes

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer

SPECIMEN COLLECTION AND PREPARATION

Collect the throat swab sample with a sterile swab, from the posterior pharynx, tonsils and other inflamed areas. Avoid touching the tongue, cheeks and teeth with the swab.

Send specimen to lab immediately (testing sensitivity decrease over time)

Swab sample may be stored and transport in a clean and dry container, samples should be process as soon as possible after collection. If this is not possible, the samples can be stored in the refrigerator (2-8°C) until 8 hours since the samples was taken.

Strep A MonlabTest® is compatible with the use of AMIES transport medium. The proportion between AMIES medium and diluent A plus diluent B should be equal or minor that 50%, taking note that the volume of diluent A plus diluent B should be the indicated in the test procedure.

Culture confirmation

Strep A MonlabTest® used to detect Group A *Streptococcus* colonies on Blood Agar plates (anaerobic incubation 72 hours/37°C).

PROCEDURES


Allow the tests, samples and diluents to reach to room temperature (15-30°C/59-86°F) prior to testing. Do not open the tube until ready to perform the assay.

To process the collected throat swab sample (see illustration 1):

Use a separate testing tube or vial for each sample (swab). Add the diluent A (5 drops) into the testing tube or vial. Add the diluent B (5 drops) and mix, the colour of the solution changes from light pink to light yellow (colourless). Put the throat swab, mix and extract as much liquid possible from the swab.


If the sample is from AMIES medium: add 5 drops of diluent A, 5 drops of diluent B and 10 drops of AMIES medium into a testing tube, then follow the procedure.

Illustration 1


Remove the Strep A MonlabTest® from its sealed pouch and use it as soon as possible. Use a separate device for each sample. Dispense exactly 3 drops into the specimen well (S). Start the timer. Read the result at 10 minutes after dispensing the sample.

Illustration 2

INTERPRETATION OF RESULTS

Illustration 3

POSITIVE: Two lines appear across the central window, in the result line region (**red** test line marked with the letter T) and in the control line region (**blue** control line marked with the letter C).

NEGATIVE: Only one **blue** band appears across the control line region marked with the letter C (control line).

INVALID: A total absence of the **blue** control coloured band regardless the appearance or not of the **red** test line. Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the result line region (T) will vary depending on the concentration of antigens in the specimen. However, neither the quantitative value, nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural controls are included in the test:
 - A blue line appearing in the control line region (C). It confirms sufficient specimen volume and correct procedural technique.

External Quality Control


Each kit contains a positive control material. Use the control swab to check that the extraction reagents and the test are working properly. Also use the controls to test that you are able to correctly perform the test procedure.

Quality Control Procedure:

Strep A Positive Control: Remove the Strep A positive control from its sealed pouch. Add the diluent A (5 drops) in a testing tube. Add the diluent B (5 drops) and mix, the colour of the solution changes from light pink to light yellow (colourless). Put the Strep A, positive control swab, mix 15 seconds and extract as much liquid possible from the swab. Discard the swab. Remove the test from its sealed pouch and dispense 3 drops of the positive control liquid into the specimen well (S).

Result: Strep A positive (see interpretation of results).

Illustration 4

LIMITATIONS

1. Strep A MonlabTest® will only indicate the presence of GAS in the specimen (qualitative detection) and should be used for the detection of GAS antigens in throat swab specimens only. Neither the quantitative value nor the rate of increase in GAS antigens concentration can be determined by this test.
2. The test must be carried out within 15 minutes after opening the sealed bag.
3. This test does not differentiate between carriers and acute infection. Pharyngitis may be caused by organisms other than Group A *Streptococcus*.
4. Strep A MonlabTest® should be used only with throat swabs or colonies taken directly from a plate. The use of swab specimens taken from other sites or the use of other samples such as saliva, sputum or urine has not been established. The quality of the test depends on the quality of the sample; proper throat swab specimens must be obtained.
5. The intensity of the test line may vary from very strong at high antigen concentration to faint when the antigens concentrations is close to the detection limit value of the test.
6. This test provides a presumptive diagnosis of GAS infections. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated must be based in the correlation of the results with further clinical observations.
7. Positive results determine the presence of Group A *Streptococcus* in throat samples; nevertheless, a positive result should be followed up with additional laboratory techniques to confirm the results.
8. Depending on the quantity of antigen in the sample and of the quantity of AMIES medium, final sensitivity could be affected.
9. A negative result is not meaningful because of it is possible the antigen concentration in the sample was too small (lower than the detection limit). If clinical symptoms persist, a Group A *Streptococcus* determination should be carried out, on a sample from an enrichment culture.
10. Mucous and/or bloody samples could cause non-specific reactions in the test. Mucous and/or bloody samples whose result is positive should be followed up with other techniques to confirm the result.

EXPECTED VALUES

There are several million cases of "Strep throat" each year. About 9,400 cases of invasive GAS disease occurred in the United States in 1999.

PERFORMANCE CHARACTERISTICS

Detection limit

The detection limit value for **Strep A MonlabTest®** is $9.8 \cdot 10^3$ CFU/mL.

Sensitivity and specificity

A retrospective evaluation with throat specimens from patients with pharyngitis symptoms was performed using **Strep A MonlabTest®** versus bacterial culture method (gold-standard). This is the evaluation criteria. The results were as follows:

	Evaluation criteria			Strep A MonlabTest® vs Evaluation criteria		
		+	-	Total	Mean Value	95% confidence interval
IC test: Strep A MonlabTest®	+	97	4	101	Sens. 83.6%	75.6-89.8%
	-	19	94	113	Esp. 95.9%	89.9-98.9%
	Total	116	98	214	VPP 96.0%	90.2-98.9%
					VPN 83.2%	75.0-89.6%

The results showed a high sensitivity and specificity to detect Group A Streptococcus using **Strep A MonlabTest®**.

N°	TP	FP	FN	TN	N
1	97	4	19	94	214

Sens. (%) CI (%)	Spec. (%) CI (%)	VPP (%) CI (%)	VPN (%) CI (%)
83.6 (75.6-89.8)	95.9 (89.9-98.9)	96.0 (90.2-98.9)	83.2 (75.0-89.6)

True positive (TP), false positive (FP), false negative (FN), true negative values (TN), and total samples (N). Sensitivity (Sens), specificity (Spec), Predictive Positive Values (PPV), Predictive Negative Values (NPV), 95% interval of confidence (CI) for **Strep A MonlabTest®**.

Cross-Reactivity

It was performed an evaluation to determine the cross reactivity of **Strep A MonlabTest®**. There is not cross reactivity with the following pathogens/substances tested:

<i>Adenovirus</i>	<i>Entamoeba histolytica</i>	<i>Norovirus GI/Norovirus GII</i>
<i>Astrovirus</i>	<i>Escherichia coli O:111; O:026; O157:H7</i>	<i>Respiratory Syncytial Virus</i>
<i>Campylobacter coli/jejuni</i>	<i>Giardia lamblia</i>	<i>Rotavirus</i>
<i>Clostridium difficile antigen GDH</i>	<i>Helicobacter pylori</i>	<i>Salmonella enteritidis/paratyphi A/typhi/typhimurium</i>
<i>Clostridium difficile Toxin A/Toxin B</i>	Hemoglobin (bovine and pig)	<i>Shigella boydii/dysenteriae/flexneri/sonnei</i>
<i>Clostridium perfringens</i>	Influenza A and B	<i>Streptococcus pneumococcal</i>
<i>Coronavirus</i>	<i>Legionella pneumophila</i>	<i>Yersinia enterocolitica O:3/O:9</i>
<i>Cryptosporidium parvum</i>	<i>Listeria monocytogenes</i>	
<i>Pig haemoglobin</i>	<i>Bovine haemoglobin</i>	

Interferences

It was performed an evaluation to determine the possible interferences of **Strep A MonlabTest®**. There are no interferences against the following substances tested:

Exogenous interferences			
Potential interfering substance	Concentration	Potential interfering substance	Concentration
Metronidazole	3.0 mg/mL	Loperamide hydrochloride (<i>Fortasec</i>)	0.15 mg/mL
Ampicillin	3.0 mg/mL	Heparin (<i>Hibor</i>)	350.0 IU/mL
Oseltamivir	0.3 mg/mL	Almagato (<i>Almax</i>)	3.0 mg/mL
Amantadine	0.3 mg/mL	Fosfamicin (<i>Monurol</i>)	3.0 · 10 ⁻³ mg/mL
Ribavirin	3.0 mg/mL	Acetylcysteine (<i>Fluimucil</i>)	3.0 mg/mL
Codeine (<i>Toseina</i>)	0.2 mg/mL	Dexketoprofen trometamol (<i>Enantyum</i>)	0.3 mg/mL
Benzocaine (<i>Angileptol</i>)	0.3 mg/mL	Levofloxacin	3.0 mg/mL
Cloperastine (<i>Flutox</i>)	0.3 mg/mL	Ciprofloxacin	0.3 mg/mL
Carbocisteine (<i>Iniston mucolítico</i>)	3.0 mg/mL	Rifampicin (<i>Rifaldin</i>)	0.3 mg/mL
Loratadine	0.3 mg/mL	Phenoxyethylpenicillin potassium	3.0 mg/mL
Dexchlorpheniramine (<i>Polaramine</i>)	0.3 mg/mL	Ambroxol hydrochloride (<i>Mucosan</i>)	0.3 mg/mL
Ebastine (<i>Ebastel</i>)	3.0 mg/mL	Macrogol 3350 (<i>Movicol</i>)	3.0 mg/mL
Acetyl Salicylic (<i>Adiro</i>)	0.3 mg/mL	Lysine Carbocysteinate (<i>Pectox</i>)	3.0 · 10 ⁻² mg/mL
Ibuprofen (<i>Espidifen</i>)	3.0 mg/mL	Hydroxyzine dihydrochloride	0.3 mg/mL
Paracetamol (<i>Dolocatil</i>)	5.0 mg/mL	Lorazepam	3.0 · 10 ⁻³ mg/mL
Metamizole (<i>Nolotil</i>)	5.0 mg/mL	Amoxicillin	3.0 mg/mL

Exogenous interferences			
Potential interfering substance	Concentration	Potential interfering substance	Concentration
Prednisone	0.3 mg/mL	Mercaptopurine	0.3 mg/mL
Omeprazole	2.0 · 10 ⁻³ mg/mL	Biotine	0.1 mg/mL
Naso GEL	0.9 mg/mL	Sore Throat Phenol spray	0.5 mg/mL
CVS Nasal Spray (<i>Cromolyn</i>)	4 mg/mL	Tobramycin	0.3 mg/mL
Afrin (<i>Oxymetazoline</i>)	5.0 · 10 ⁻² mg/mL	Mupirocin	2.5 · 10 ⁻² mg/mL
CVS Nasal Drops (<i>Phenylephrine</i>)	10.0 mg/mL	Fluticasone Propionate	5.0 · 10 ⁻² mg/mL
ZICAM	1.0 mg/mL	Homeopathic	DIL 1/10
Phenylpropanolamine	0.3 mg/mL		
Endogenous interferences			
Human haemoglobin	5.0 µg/mL	Human lactoferrin	5.0 µg/mL
Human transferrin	5.0 µg/mL	Mucine	0.5 %
Human calprotectin	5.0 µg/mL	Human blood	5.0 %

Repeatability and Reproducibility

It was performed a study of repeatability and reproducibility using different internal samples, negative and positive. There are no differences observed within the evaluations.

REFERENCES

- Weiss KA LM. Group A Streptococcus invasive infections: a review. *Can J Chem.* 1997;40(1):18-25.
- Walker MJ, Barnett TC, McArthur JD, Cole JN, Gillen CM, Henningham A, et al. Disease manifestations and pathogenic mechanisms of group A Streptococcus. *Clin Microbiol Rev.* 2014;27(2):264-301.
- Thompson TZ, McMullen AR. Group a streptococcus testing in pediatrics: The move to point-of-care molecular testing. *J Clin Microbiol.* 2020;58(6).
- Cohen JF, Bertille N, Cohen R, Chalumeau M. Rapid antigen detection test for group A streptococcus in children with pharyngitis. *Cochrane Database Syst Rev.* 2016;7(7).
- Bessen DE, McShan WM, Nguyen SV, Shetty A, Agrawal S TH. Molecular Epidemiology and Genomics of Group A Streptococcus. *Infect Genet Evol.* 2015; 33:393-418.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS

	Manufacturer		For <i>in vitro</i> diagnostic use only
	Don't re-use		Consult instructions for use
	Contains sufficient for <n> tests		Keep dry
	Catalogue Code		Temperature limitation
	Lot Number		Use by
	Date of manufacture		Sample diluent
	Positive Control		

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- Weiss KA LM. Group A Streptococcus invasive infections: a review. Can J Chem. 1997;40(1):18-25.
- Walker MJ, Barnett TC, McArthur JD, Cole JN, Gillen CM, Henningham A, et al. Disease manifestations and pathogenic mechanisms of group A Streptococcus. Clin Microbiol Rev. 2014;27(2):264-301.
- Thompson TZ, McMullen AR. Group a streptococcus testing in pediatrics: The move to point-of-care molecular testing. J Clin Microbiol. 2020;58(6).
- Cohen JF, Bertille N, Cohen R, Chalumeau M. Rapid antigen detection test for group A streptococcus in children with pharyngitis. Cochrane Database Syst Rev. 2016;7(7).
- Bessen DE, McShan WM, Nguyen SV, Shetty A, Agrawal S TH. Molecular Epidemiology and Genomics of Group A Streptococcus. Infect Genet Evol. 2015; 33:393-418.

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