## C. DIFFICILE TOXINS A+B MonlabTest®

MO-076011 20 TESTS

One-step immunochromatographic test for the differential detection of toxin A and toxin B from *C. difficile* in faeces

For *in vitro* use only. Store at 2 - 30°C.

### INTENDED USE

*C. difficile* toxins A+B MonlabTest chromatographic immunoassay provides a procedure for a qualitative detection of Clostridium difficile toxin A (TcdA) and toxin B (TcdB) in two separate bands. A positive band of either toxin is a sign of an underlying Clostridium difficile infection (CDI), which should be taken into consideration by the clinician.

Unlike other immunoassay systems (ELISA and rapid tests) that only allow detection of TcdA and TcdB in combination and without any differentiation between the two, the *C. difficile* toxins A+B MonlabTest allows you to run a single assay with a single strip that differentiates between both of the toxins using two separated bands - a red band below the blue control band when TcdA is present, and another red band above the control band when TcdB is in the sample (see Fig. 1).

The test is based on the immunological capture of coloured microparticles during their passage along a membrane on which specific monoclonal antibodies against TcdA and TcdB have been immobilised at two separated locations.

### SUMMARY

#### Presumptive diagnosis of Clostridium difficile infection (CDI)

*C. difficile* produces two different toxins that constitute the essential virulence factors for CDI induction. Recent investigations<sup>1,2</sup> have proven that each of the two alone will induce disease in hamsters.

*C. difficile* infection is considered responsible for approximately 25% of the diarrhoea incidents related to the consumption of antibiotics such as clindamycin, second and third generation cephalosporins, gyrase-inhibitors, ampicillin or amoxicillin. In addition to the diarrhoea symptoms, the disease can lead to pseudo-membranous colitis (PMC), requiring urgent treatment with antibiotics effective against *C. difficile* (metronidazole / vancomycin) and which, without treatment, may severely compromise the life of the patient. CDI mortality can be as high as 6% to 30%, particularly when the patient suffers from PMC. Patients suffering from CDI induced by previous treatment can lead to an increased hospital stay of 6-10 days at an additional cost of 6000-8000 euros<sup>3</sup>.

#### C. difficile description

*Clostridium difficile* is an anaerobic, gram-positive, spore-forming bacteria that is carried by approximately 5% of the healthy population. Following hospitalisation, the carrier rate climbs to approximately 30%.

Children are colonised by *C. difficile* early after birth, but they do not appear to suffer any clinical symptoms despite intensive investigations. The current discussion suggests a lack of enteric receptors for the toxins or differences in faecal pH in children that prevent toxin activity.

As previously mentioned, *C. difficile* can produce different toxins<sup>4</sup>:

- **Toxin A (TcdA)** (308 kDa) is called an enterotoxin given that it can induce full symptoms in the hamster animal model. TcdA also displays high cytotoxicity but only on specific TcdA sensitive cells such as HT-29<sup>5</sup>.
- **Toxin B (TcdB)** (279 kDa) is classified as a cytotoxin. In the majority of cells cultured in laboratories (e.g. Vero, CHO or HeLa) it is about 1,000 times more potent than toxin A. The TcdB amino acid sequence varies between different strains<sup>6</sup>.

#### TcdA/TcdB detection patterns

The following *C. difficile* strains can be identified based on the production of toxins:

- **Non-toxicogenic strains** are non-pathogenic and lack TcdA and TcdB production as well as production of the binary toxin.
- **TcdA+ / TcdB+ strains** are the most common CDI inducing strains of which ribotypes 001, 014 and 078 are the most prevalent in Europe<sup>7</sup>.
- **TcdA- / TcdB+ strains** were first identified by Depitre et al. in Belgium<sup>8</sup>. They are considered pathogenic despite not producing toxin A<sup>9</sup>. Ribotype 017 strains belong to this group and were responsible for various endemic outbreaks in North America.
- **TcdA+ / TcdB- strains** can be identified directly in stool samples with the *C. difficile* toxins A+B MonlabTest due to provides a simple, fast and separate detection of toxin A and toxin B in a single test. Very few isolates of these strains have been found until now. Recent results show that TcdB mutant *C. difficile* strains still induce CDI in hamsters owing to the production of toxin A<sup>2</sup>. These results appear to indicate the presence of these strains in clinical samples.

### BASIC PRINCIPLES OF THE TEST

The *C. difficile* toxins A+B MonlabTest employs a combination of:

- 1) Red latex particles conjugated to a specific toxin A antibody that cooperates with another antibody specific for toxin A and that is located on the membrane, below the control band,
- 2) Other red latex particles conjugated to a specific toxin B antibody that cooperates with another specific toxin B antibody that is located on the membrane, above the control band,
- 3) Blue latex particles conjugated to an antigen recognised by an antibody specific for that antigen and bound to the membrane, serving as the control band test.

To run the test, the sample is first treated with a sample diluent buffer (provided in the kit) that extracts the toxins from the stool matrix. Following extraction, an aliquot of the supernatant needs to be added to the test strip followed by a 15 minutes wait.

When the extract flows through the test membrane, the coloured particles begin to migrate. In the event of a positive sample, the specific antibodies on the membrane will capture antigen-covered coloured particles.

The pattern of lines obtained after 15 minutes of incubation at room temperature are used to interpret the result (see Fig.1).

### MATERIAL PROVIDED

- 20 cassettes
- 20 vials with diluent (1.5mL)
- Disposable plastic pipettes
- Instructions for use

### MATERIALS NOT INCLUDED IN THE KIT

- Vortex
- Timer

### PRECAUTIONS

1. Patient samples (faeces) should be handled with care as they may contain infectious agents. All the necessary protections should be used throughout handling (disposable gloves, goggles, laboratory coat and others).
2. The sample diluent buffer contains Sodium Azide as an antimicrobial agent. Avoid direct contact with the skin and mucous membranes. Dispose of appropriately. The buffer should not be used if there are signs of contamination or precipitation.
3. Do not eat, drink, smoke, store or prepare food in areas where the reagents and the samples are handled.
4. Once the work has been concluded, remove the gloves and first disinfect your hands with alcoholic disinfectant. Secondly, wash your hands with soap. Finally, the sink that has been used needs to be decontaminated with sporocidal disinfectants, as the *C. difficile* spores are not eliminated with alcohol.
5. Do not exchange components between kits with different lot numbers.
6. If the test is stored refrigerated, allow kit components and stool samples to reach room temperature before use, as cold reagents and/or samples may reduce test performances. About 20-30 minutes are usually sufficient for reaching room temperature.
7. All reagents are for *in vitro* use only.
8. Do not use kit components beyond their expiry date.
9. In case the primary packaging is damaged (aluminium pouch or diluent buffer vial) the product should not be used even though none of the components have been damaged.
10. It is very important to add the correct volume of extracted sample to the reactive device. If the volume is lower than indicated, chromatography may not occur because the sample may not reach the reaction area. If higher volumes are used, brown lines may appear instead of red or blue ones.
11. All products are for single use only and should be discarded according to current legislation.
12. Do not use the test if any coloured lines appear in the result area prior to performing the test.
13. It is critical to collect the correct sample quantity: approximately **110 mg** of a **solid sample** (a **small ball of 5 mm in diameter**). If the **sample is semi-liquid** (unable to take it with a pipette), take a sample amount so that it completely **covers the grooves of the stick** attached to the vial cap and **110 µl** of a liquid sample (**4 drops** if the disposable pipettes included in the kit are used). These quantities are extracted in the sample diluent supplied with the vials provided in the kit. An excess of sample in relation to the amount of buffer added prevents the chromatography from running correctly; this is especially critical in the case of solid samples, since it is harder to obtain an appropriate quantity.
14. Do not discard the external box of the kit until its content has been totally used. The external box contains essential information regarding the EC marking and lotification.

### STORAGE AND STABILITY

*C. difficile* toxins A+B MonlabTest can be stored at any temperature between 2 and 30°C.

Its expiry date is printed on the tube or on the aluminium wrap.

### SAMPLES

- This test is designed to analyse liquid or semi-liquid stool samples. Solid samples may be analysed, but it is unusual given that the leading symptom of *C. difficile* infection is diarrhoea.
- Do not use samples that have been collected in means of conveyance or to which enrichment media / preserving agents have been added (e.g. formalin, SAF, PVA or similar) as they may interfere with the test.
- The analysis of untreated fresh samples is recommended. If preservation is needed, they should not be stored in the refrigerator (+2-8°C) for longer than 1 to 2 days. For longer storage, samples must be frozen at -20°C, however please bear in mind that some samples turn negative when they have been frozen.
- Pay special attention when analyzing hemorrhagic samples as they often give false positive results when the blood content is high. An indicator of this destabilization of the test is usually the alteration of the colour in the control band (tends to show a purple or dark blue colour).
- Ensure that frozen samples have completely defrosted and reached room temperature prior to proceeding with their analysis.
- Avoid repeatedly freezing and thawing the stool samples as the integrity of the toxins may suffer.

### STOOL SAMPLES PREPARATION

**General remark:** all the necessary protections should be used throughout the test procedure due the handling of infectious samples. Once finished, do not forget to comply with the hygiene procedures detailed in point 4 of the "Precautions" section.

The sample diluent buffer is the same as for the MO-076012 *C. DIFFICILE* Ag (GDH) MonlabTest.

The protocol for the preparation of stool samples is as follows:

1. Homogenize previously the sample to get an aliquot as much representative as possible.

- Unscrew the cap from the vial with caution in order not to spill the sample diluent buffer. With **solid samples**, take with the stick attached to the vial cap an approximate amount of **110 mg** of faeces (a small portion of **5 mm in diameter**). If the sample is **semi-liquid** (unable to take it with a pipette), take a sample amount so that it completely **covers the grooves of the stick** attached to the vial cap. With **liquid samples**, take a volume of **110 µl (4 drops)** if the disposable pipettes included in the kit are used).
- Carefully add the sample into the vial containing the dilution buffer. Screw the cap well and shake vigorously to ensure a homogeneous mixture

## PROCEDURE

- Take the reaction device out of its aluminium pouch. Discard the desiccant as it only functions as to preserve the test from any excess of humidity.
- Put the vial upside down and add **3 drops** to the sample area of the reaction device (round window marked with an arrow).
- Wait for **15 minutes** and read the results.

## INTERPRETATION OF THE RESULTS

The five strips shown in Fig. 1 exemplify the various results that can be obtained using the *C. difficile* toxins A+B MonlabTest. There are three different coloured bands:

- Blue band:** the control band that indicates a correct performance of the test.
- Upper red band:** TcdB positive sample.
- Lower red band:** TcdA positive sample.

The blue band (control) should always appear. The additional appearance of any red band in the strip indicates the presence of *C. difficile* that produces TcdA and/or TcdB in the analysed sample.

- Strip 1:** NEGATIVE result: the sample does not contain *C. difficile* or it contains a strain that does not produce TcdA/TcdB. A single **BLUE** horizontal band appears within the central area of the reactive device (in the Simple format it is aligned with the letter "C" marked on the cassette). This is the control band and it should always appear as an indication of the chromatography running smoothly.
- Strips 2-4:** POSITIVES results:
  - Strip 2:** Detection of TcdA: a **BLUE** (control band) and a **RED** band appear just below the control band aligned with the letter "T1" marked on the cassette. The intensity depends on the concentration of toxin A in the sample.
  - Strip 3:** Detection of TcdB: a **BLUE** (control band) and a **RED** band appear just above the control band aligned with the letter "T2" marked on the cassette. The intensity depends on the concentration of toxin B in the sample.
  - Strip 4:** Detection of both TcdA and TcdB: a **BLUE** (control band) and two **RED** bands (one above [TcdB] and one below [TcdA] the control band).
- Strip 5:** INVALID result: the blue control band does not appear, or the blue colour of this band is clearly altered (very dark blue or a purple colour), also nonspecific colours in the positive bands (different from red) are invalid results. In general, any colour combination different from those indicated in 1-4 strips indicates an anomalous test performance. Possible reasons are:
  - some reagents have got damaged or the test has expired.
  - the sample was not prepared according to the instructions of use.
  - high blood content in the sample.

In the event of an invalid result it is recommended that another test is run, strictly following the protocol described in this manual. For blood samples, it is advisable the use of an alternative technique because the problem of destabilization does not usually depend on the strip but of the sample matrix. Other commercial rapid tests gave similar results with these bloody samples.

Any line appearing after the standard 15 minutes reaction time is of NO diagnostic value.

NOTE: The final and definitive diagnosis CDI/PMC is established by the clinician. This test only detects TcdA/TcdB in a sample but does not constitute a case to confirm whether a person has CDI.

Strip 1	Strip 2	Strip 3	Strip 4	Strip 5
				
TcdB- TcdA-	TcdB- TcdA+	TcdB+ TcdA-	TcdB+ TcdA+	Invalid test

Fig. 1: Pattern of possible results

## LIMITATIONS OF THE PROCEDURE

- C. difficile* toxins A+B MonlabTest analyses liquid or semi-liquid human stool samples; solid samples may be used; however, the test has not been optimised for their use since, in rare occasions, toxins sequestration phenomena have been observed with such solid matrices.
- This test is qualitative and not quantitative, although the intensity of the positive bands is associated with the quantity of toxins that are detected in the stool sample.
- Over 200 stool samples were evaluated to ensure the correct performance of the test. The correlation of the results with other techniques (ELISA and Cytotoxicity) was good. However, this study does not exclude interferences in the performance of the tests with other stool samples.
- Weak signals may be due to excessively low amounts of sample. In the event of this occurring, the test should be repeated with a greater amount of sample whilst maintaining the recommended sample diluent ratio (see "Sample preparation" section).
- An excess of sample can significantly slow down the development of the test or even prevent the test from running (control band remains invisible). In the event of this occurring, the test should be repeated with a reduced sample amount. This is particularly relevant in the analysis of solid samples.
- A negative result does not fully exclude the possibility of infection with a *C. difficile* (CDI). The test result must be interpreted in relation to the clinical symptoms of the patient. In addition, it is important to keep in mind that toxins are fragile molecules that can be easily degraded, for example owing to inappropriate storage of the sample or the presence of inhibitors. Under these conditions, toxin concentration may be below the test detection limit (see "Analytical Sensitivity" section).
- A positive result obtained from solid samples must be interpreted with a great amount of caution. In principle, diarrhoea is the leading symptom of *C. difficile* (CDI) infection CDI, and a solid stool implies an absence of diarrhoea. The person performing the test must provide the clinician with information on the nature of the sample. This, in combination with the patient's medical history, will allow the clinician to establish the most accurate diagnosis possible.
- A certain degree of cross-reaction has been observed with stool samples that are strongly positive for *Entamoeba histolytica*. The reaction was negative when using a preparation of the isolated parasite (without stool matrix). Other ELISA and rapid tests that are available on the market yielded similar results for these samples.
- It is observed that faecal samples with a high blood content may interfere negatively with the test. In this case, specificity problems may appear with *C. difficile* negative samples. This destabilization of the test is often accompanied by an alteration in the colour of the control band; a purple or dark blue colour appears instead of the expected light blue (see images in "Interpretation of the results" section).
- C. difficile* colonisation rates of up to 50% have been reported in infants. A high rate has also been reported in cystic fibrosis patients. Normally, these two groups of patients remain asymptomatic and do not require specific treatment. These positive results that are of no clinical significance should be treated with caution as affected patients do not require treatment for *C. difficile* infection.

## ANALYTICAL SENSITIVITY

In order to determine the sensitivity of the *C. difficile* toxins A+B MonlabTest toxins A and B from tgc BIOMICs were used diluted in the sample diluent buffer of this test. Most manufactured batches detect a concentration of 0.75 ng/mL of both toxins.

## DIAGNOSTIC SENSITIVITY AND SPECIFICITY

The average values of sensitivity and specificity for the *C. difficile* toxins A+B MonlabTest obtained from all the evaluations made with the test are:

Sensitivity: 97.4 %  
Specificity: 96.1 %

*C. difficile* toxins A+B MonlabTest shows an excellent performance with sensitivity and specificity values above 95%.

Two of the evaluations performed with the test are showed below:

### A. External Evaluation:

*C. difficile* toxins A+B MonlabTest was evaluated at a Spanish hospital by measuring a total of 150 negative samples and 44 positive samples according to the reference technique, the Cytotoxicity. All samples were fresh.

The results obtained were:

Sensitivity: 92.5 %  
Specificity: 95.5 %

### B. Internal Evaluation

*C. difficile* toxins A+B MonlabTest was evaluated internally by measuring 242 negative samples and 105 positive samples according to the Cytotoxicity developed at the place of origin of the samples (different Spanish hospitals). All samples analysed were frozen. The results obtained were the follows:

Sensitivity: 99.0 %  
Specificity: 96.3 %

## REPEATABILITY

Purified toxins A and B were used to design a sensitivity curve to measure the sensitivity of the *C. difficile* toxins A+B MonlabTest in different conditions. Dilutions were two-fold and the samples were assayed by the same person in triplicate during a single session. The results differed by less than a factor of 2 indicating a high accuracy of the test.

## REPRODUCIBILITY

**INTER-DAY PRECISION:** A sensitivity curve was measured with the same lot of *C. difficile* toxins A+B MonlabTest on four different days. The results were very reproducible (the same level of sensitivity is obtained for both TcdA and TcdB over the four measurement days).

**INTER-OPERATOR PRECISION:** Five people with no prior training measured a sensitivity curve in duplicate. Differences were observed in the stronger curve dilutions (weaker signals), never exceeded a factor of 2.

**INTER-LOT PRECISION:** Three different lots were used to measure the sensitivity curve in duplicate. The analysis was performed by a single person on the same day. Differences of a dilution factor of 2 were appreciated, which are acceptable and tolerable for the test that was carried out.

The differences found in the "Reproducibility" sections are acceptable for this qualitative immunochromatographic test due to the inherent variability associated to this technique.

## HOOK EFFECT

Several publications indicate that in the case of severe *C. difficile* infections, the highest concentration of toxins found in faeces is around 112 ng/mL<sup>10</sup>. To evaluate the effect of a very high concentration of toxins, values well above those that can be found among the population, affected by CDI, concentrations of toxin A and toxin B up to 5000 ng/mL were measured with the test. This concentration is about 400 times the LPC test for toxin A (12.5 ng/mL), about 1500 times the LPC test for toxin B (3 ng/mL) and about 40 times the highest concentration of toxins found in faeces from patients. No decrease in the intensity of the positive signals was observed

## INTERFERING SUBSTANCES

The substances indicated in the below table, at the specified concentration, did not interfere with the results of the test when they were added to stool samples (positive and negative ones):

Racecadotril	5% (p/v)	Ibuprofen	20% (p/v)
Cimetidine	10% (p/v)	Acetylsalicylic acid	30% (p/v)
Loperamide	5% (p/v)	Edulcorant	5% (p/v)
Metronidazole	10% (p/v)	Palmitic acid	40% (p/v)
Omeprazole	3% (p/v)	Barium Sulfate	5% (p/v)
Ampicillin	15% (p/v)	Mucin	5% (p/v)

## CROSS-REACTIVITY

This study was developed in two different places:

- MONLAB:** *C. difficile* toxins A+B MonlabTest was evaluated against strongly positive samples for the following microorganisms:

*Adenovirus - Rotavirus - Norovirus - Astrovirus - Helicobacter pylori - Entamoeba histolytica - Giardia lamblia - Cryptosporidium parvum.*

A certain degree of cross-reaction was observed with strongly positive samples for *Entamoeba histolytica* (see point 8 in the "Limitations of the Procedure" section).

- NATIONAL HOSPITAL CENTER (España):** *C. difficile* toxins A+B MonlabTest was evaluated against different micro-organisms likely to be present in the intestinal tract at any time in a sufficiently high concentration. Tests were carried out on bacterial suspensions at a concentration of 10<sup>8</sup> cfu/mL. *C. difficile* toxins A+B MonlabTest showed no cross reactivity with any of the micro-organism listed below:




*Aeromonas caviae, Aeromonas hydrophila, Bacillus cereus, Bacillus spp, Bacteroides nordii, Campylobacter jejuni, Candida albicans, Clostridium butyricum, Clostridium cadaveris, Clostridium perfringens, Clostridium sordellii, Enterobacter aerogenes, Enterobacter cloacae, Enterococcus faecalis, Enterococcus faecium, Escherichia coli 1, Escherichia coli 2, Klebsiella pneumoniae, Lactobacillus goseri, Listeria monocitogenes, Plesiomonas shigelloides, Proteus mirabilis, Pseudomonas aeruginosa, Salmonella enteritidis, Salmonella typhi, Salmonella typhimurium, Serratia marcescens, Shigella dysenteriae, Shigella flexnerii, Shigella sonnei, Staphylococcus aureus, Staphylococcus epidermidis, Vibrio cholerae, Vibrio parahaemolyticus, Yersinia enterocolitica.*

## REFERENCES

- Lyras O, O'Connor JR, Howarth PM, Sambol SP, Carter GP, Phumoonna T, Poon R, Adams V, Vedantam G, Johnson S, Gerding DN and Rood JI. Toxin B is essential for virulence of *Clostridium difficile*. Nature. Apr 2009, Vol: 458(7242) pp: 1176-1179.
- Kuehne SA, Cartman ST, Heap JT, Kelly ML, Cockayne A and Minton NP. The role of toxin A and toxin B in *Clostridium difficile* infection. Nature. Oct 2010, Vol: 467(7316) pp: 711-713.
- Vonberg RP, Reichardt C, Behnke M, Schwab F, Zindler S and Gastmeier P. Costs of nosocomial *Clostridium difficile* associated diarrhoea. J Hosp Infect. Sep 2008, Vol: 70(1) pp: 15-20.
- Rupnik M, Dupuy B, Fairweather NF, Gerding DN, Johnson S, Just I, Lyster DM, Popoff MR, Rood JI, Sonenshein AL, Thelestam M, Wren BW, Wilkins TD and von Eichel-Streiber C. Revised nomenclature of *Clostridium difficile* toxins and associated genes. J Med Microbiol. Feb 2005, Vol: 54(Pt 2) pp: 113-117.
- Tucker KD, Carrig PE and Wilkins TD. Toxin A of *Clostridium difficile* is a potent cytotoxin. J Clin Microbiol. May 1990, Vol: 28(5) pp: 869-871.
- Soehn F, Wagenknecht-Wiesner A, Leukel P, Kohl M, Weidmann M, van Eichel-Streiber C and Braun V. Genetic rearrangements in the pathogenicity locus of *Clostridium difficile* strain 8864-implications for transcription, expression and enzymatic activity of toxins A and B. Mol Gen Genet. May 1998, Vol: 258(3) pp: 222-232.
- Barbut F, Mastrantonio P, Delmée M, Brazier J, Kuijper E and Poxton I. European Study Group on *Clostridium difficile* (ESGCD). Prospective study of *Clostridium difficile* infections in Europe with phenotypic and genotypic characterisation of the isolates. Clin Microbiol Infect. Nov 2007, Vol: 13(11) pp: 1048-1057.
- Depitre C, Delmee M, Avesani V, L'Haridon R, Roels A, Popoff M and Corthier G. Serogroup F strains of *Clostridium difficile* produce toxin B but not toxin A J Med Microbiol. Jun 1993, Vol: 38(6) pp: 434-441.

- Drudy D, Harnedy N, Fanning S, Hannan M and Kyne L. Emergence and control of fluoroquinolone-resistant, toxin A-negative, toxin B-positive *Clostridium difficile*. Infect Control Hosp Epidemiol. Aug 2007, Vol: 28(8) pp: 932-940.
- Alex B. Ryder, Ying Huang, Haijing Li, Min Zheng, Xiaobo Wang, Charles W. Stratton, Xiao Xu and Yi-Wei Tang. Assesment of *Clostridium difficile* infections by quantitative detection of TcdB toxin by use of a real-time cell analysis system. Journal of Clinical Microbiology, Nov 2010, pp: 4129-4134.

## PRESENTATION SYMBOLS USED FOR IVD COMPONENTS AND REAGENTS

	Manufactured by		For <i>in vitro</i> diagnostic use
	Do not re-use		Please read pack insert
	Contains sufficient for <n> tests		Dry storage
	Catalogue number		Store at
	Lot number		Expiry date
	This product fulfils the requirements of Directive 98/79/EC on <i>in vitro</i> diagnostic medical devices		Dilution buffer

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1. Lyras O, O'Connor JR, Howarth PM, Sambol SP, Carter GP, Phumoonna T, Poon R, Adams V, Vedantam G, Johnson S, Gerding DN and Rood JI. Toxin B is essential for virulence of *Clostridium difficile*. Nature. Apr 2009, Vol: 458(7242) pp: 1176-1179.
2. Kuehne SA, Cartman ST, Heap JT, Kelly ML, Cockayne A and Minton NP. The role of toxin A and toxin B in *Clostridium difficile* infection. Nature. Oct 2010, Vol: 467(7316) pp: 711-713.
3. Vonberg RP, Reichardt C, Behne M, Schwab F, Zindler S and Gastmeier P. Costs of nosocomial *Clostridium difficile* associated diarrhoea. J Hosp Infect. Sep 2008, Vol: 70(1) pp: 15-20.
4. Rupnik M, Dupuy B, Fairweather NF, Gerding DN, Johnson S, Just I, Lyerly DM, Popoff MR, Rood JI, Sonenshein AL, Thelestam M, Wren BW, Wilkins TD and von Eichel-Streiber C. Revised nomenclature of *Clostridium difficile* toxins and associated genes. J Med Microbiol. Feb 2005, Vol: 54(Pt 2) pp: 113-117.
5. Tucker KD, Carrig PE and Wilkins TD. Toxin A of *Clostridium difficile* is a potent cytotoxin. J Clin Microbiol. May 1990, Vol: 28(5) pp: 869-871.
6. Soehn F, Wagenknecht-Wiesner A, Leukel P, Kohl M, Weidmann M, van Eichel-Streiber C and Braun V. Genetic rearrangements in the pathogenicity locus of *Clostridium difficile* strain 8864-implications for transcription, expression and enzymatic activity of toxins A and B. Mol Gen Genet. May 1998, Vol: 258(3) pp: 222-232.
7. Barbut F, Mastrantonio P, Delmée M, Brazier J, Kuijper E and Poxton I. European Study Group on *Clostridium difficile* (ESGCD). Prospective study of *Clostridium difficile* infections in Europe with phenotypic and genotypic characterisation of the isolates. Clin Microbiol Infect. Nov 2007, Vol: 13(11) pp: 1048-1057.
8. Depitre C, Delmee M, Avesani V, L'Haridon R, Roels A, Popoff M and Corthier G. Serogroup F strains of *Clostridium difficile* produce toxin B but not toxin A J Med Microbiol. Jun 1993, Vol: 38(6) pp: 434-441.
9. Drudy D, Harnedy N, Fanning S, Hannan M and Kyne L. Emergence and control of fluoroquinolone-resistant, toxin A-negative, toxin B-positive *Clostridium difficile*. Infect Control Hosp Epidemiol. Aug 2007, Vol: 28(8) pp: 932-940.
10. Alex B. Ryder, Ying Huang, Haijing Li, Min Zheng, Xiaobo Wang, Charles W. Stratton, Xiao Xu and Yi-Wei Tang. Assessment of *Clostridium difficile* infections by quantitative detection of TcdB toxin by use of a real-time cell analysis system. Journal of Clinical Microbiology, Nov 2010, pp: 4129-4134.

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