## One-Step Barbiturates MonlabTest®

Only for professional *in vitro* diagnostic use. Store at 2 - 30°C.

### INTENDED USE

The One-Step Barbiturates MonlabTest® is a lateral flow, one-step immunoassay for the qualitative detection of barbiturates, primarily secobarbital, in human urine at a cut-off of 300ng/mL. This product is used to obtain a visual, qualitative result and is intended for professional use. The assay should not be used without proper supervision and is not intended for over the counter sale to lay persons.

This assay provides only a preliminary analytical test result. A more specific alternative chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) has been established as the preferred confirmatory method by the National Institute on Drug Abuse (NIDA). Clinical consideration and professional judgment should be applied to any drugs of abuse test result, particularly when preliminary positive results are indicated.

### SUMMARY

Barbiturates are a class of central nervous system depressants. Phenoarbital is a long acting barbiturate derivative that has been used as a daytime sedative and very extensively as an anticonvulsant. Pentobarbital and secobarbital are two examples of a short acting barbiturate sedative. Abuse of barbiturates can lead not only to impaired motor coordination and mental disorder, but also to respiratory collapse, coma and even death. Barbiturates are taken orally, rectally, or by intravenous and intramuscular injections. Short-acting barbiturates will generally be excreted in urine as metabolites, while the long-acting barbiturates will primarily appear unchanged.

Urine based screening tests for drugs of abuse range from simple immunoassay tests to complex analytical procedures. The speed and sensitivity of immunoassays have made them the most widely accepted method for screening urine for drugs of abuse. The One-Step Barbiturates MonlabTest® is based on the principle of the highly specific immunochemical reactions of antigens and antibodies which are used for the analysis of specific compounds in biological fluids. This test is a rapid, visual, competitive immunoassay that can be used for the qualitative detection of barbiturates in human urine at 300ng/mL cut-off concentration.

### PRINCIPLE

The One-Step Barbiturates MonlabTest® is a one-step immunoassay in which a chemically labeled drug (drug-protein conjugate) competes with the drug which may be present in urine for limited antibody binding sites. The test device contains a membrane strip which has been pre-coated with drug-protein conjugate on the test band region. A colored anti-barbiturate antibody-colloidal gold conjugate pad is placed at the end of the membrane. The colored antibody-colloidal gold conjugate moves along with urine, chromatographically by capillary action, across the membrane.

In the absence of drug in the urine, the colored antibody colloidal gold conjugate attaches to the drug-protein conjugate on the test band region to form a visible line as the antibody/drug-protein conjugate complexes. Therefore, the formation of a visible

precipitant in the test band region occurs when the test urine is **negative** for the drug. When the drug is present in the urine, the drug/metabolite antigen competes with the drug-protein conjugate on the test band region for the limited antibody sites. When a sufficient amount of drug is present, it will fill the limited antibody binding sites. This will prevent attachment of the colored antibody-colloidal gold conjugate to the drug-protein conjugate zone on the test band region. Therefore, absence for the color band on the test region indicates a **positive** result.

A control or reference band with a different antigen/antibody reaction is also added to the immunochromatographic membrane strip to indicate that the test is performed properly. This control line should always appear regardless of the presence of drug of metabolite. This means that **negative** urine will produce **two** colored bands, and **positive** urine will produce only **one** band. The presence of this colored band in the control region also serves as verification that 1) sufficient volume has been added, and 2) that proper flow was obtained.

### MATERIALS SUPPLIED

- 40 individually wrapped test which include one disposable pipette each. Each test contains a membrane coated with drug-protein conjugate and a colloidal gold conjugate pad coated with anti-barbiturate monoclonal antibody.
- One instruction sheet.

### MATERIAL REQUIRED BUT NOT PROVIDED

- Specimen collection container.
- Timer.

### STORAGE AND STABILITY

The test kit should be stored refrigerated or at temperature 2-30°C. Each device should remain in its sealed pouch for the duration of the shelf-life.

### PRECAUTIONS

- FOR IN-VITRO DIAGNOSTIC USE.
- For professional use only.
- Urine specimens may be potentially infectious. Proper handling and disposal methods should be established.
- Avoid cross-contamination of urine samples by using a new specimen collection container and specimen pipette for each urine sample.
- Do not use it after expiration date.

### SPECIMEN COLLECTION AND HANDLING


The One-Step Barbiturates MonlabTest® is formulated for use with urine specimens. Fresh urine does not require any special handling or pretreatment. Urine samples should be collected such that testing may be performed as possible after the specimen collection, preferably during the same day. The specimen may be refrigerated at 2-8°C for 2 days or frozen at -20°C for a longer period of time. Specimens that have been refrigerated must be equilibrated to room temperature prior to testing. Specimens previously frozen must be thawed, equilibrated to room temperature, and mixed thoroughly prior to testing.

**Note:** Urine specimens and all materials coming in contact with them should be handled and disposed of as if capable of transmitting infection. Avoid contact with skin by wearing gloves and proper laboratory attire.

## TEST PROCEDURE


- Review "Specimen collection" instructions. Test device, patient's samples, and controls should be brought to room temperature (10-30°C) prior to testing. Do not open pouches until ready to perform the assay.
- Remove the test device from its protective pouch (bring the device to room temperature before opening the pouch to avoid condensation of moisture on the membrane). Label the device with patient or control number.
- Hold the dropper vertically and transfer 3 full drops of urine (approx. 0.15ml) to the specimen well (S) of the test device, and then start the timer. Use a separate pipette and device for each sample or control.
- Read result between 3 to 8 minutes after the addition of samples.** Do not read result after 8 minutes.

## INTERPRETATION OF RESULTS


### Positive

The control line appears in the viewing window, but the test line is not visible.

### Negative

**Two pink-rose lines (bands) are visible in the control (C) and test (T) region of the viewing window.** The intensity of the test line may be less than that of the control line; this still means negative result.

### Invalid

The test is invalid if the control line is not visible at five minutes. The test failed, or the test procedure was not followed properly. Verify the test procedure and repeat the test with a new test device.

**Note:** A very faint line on the test region indicates that the barbiturate in the sample is near to the cut-off level for the test. These samples should be re-tested or confirmed with a more specific method before a positive determination is made.

## LIMITATIONS OF PROCEDURE

- The assay is designed for use with human urine only.
- A positive result with any of the tests indicates the presence of a drug/metabolite only and does not indicate or measure intoxication.
- There is a possibility that technical or procedural errors as well other substances or factors not listed may interfere with the test and cause false results. See **SPECIFICITY** for lists of substances that will produce positive results, or that do not interfere with test performance.
- If it is suspected that the samples have been mislabeled or tampered with, a new specimen should be collected and the test should be repeated.

## QUALITY CONTROL

Good laboratory practice recommends the use of control materials to ensure proper kit performance. Before using a new kit with patient specimens, positive and negative controls should be tested. Quality control specimens are available from commercial sources. When testing the positive and negative controls, use the same assay procedure as with a urine specimen.

## PERFORMANCE CHARACTERISTICS

### Accuracy

The accuracy of the One-Step Barbiturates MonlabTest® was evaluated in comparison to a commercially available immunoassay at a cut-off of 300ng/mL for barbiturates. One hundred (100) urine samples, collected from presumed non-user volunteers, have been tested by both procedures with 100% agreement.

Fifty two (52) urine samples, (obtained from a clinical laboratory where they were confirmed as positives by GC/MS) were tested with the One-Step Barbiturates MonlabTest® and a commercially available one-step immunoassay. Of the forty five (45) samples with barbiturates concentration  $\geq 300$ ng/mL, all were found to be positives by both methods (100% agreement). Of the seven (7) samples with concentration from 150 to 250ng/mL, all were determined negatives by both methods (100% agreement).

### Reproducibility

The reproducibility of the One-Step Barbiturates MonlabTest® was evaluated at four different sites using blind controls. Of the sixty (60) samples without secobarbital all were determined negative. Of the sixty (60) samples with secobarbital concentration of 600ng/mL, all were determined positive.

### Precision

The precision of the One-Step Barbiturates MonlabTest® was determined by conducting the test with spiked controls. The control at 150ng/mL should give negative result and the control at 600ng/mL should give positive result.

Concentration (ng/mL)	Number Tested	Correct Result	Correct Result%
150	50	50	100
600	50	50	100

### Specificity

The specificity for the One-Step Barbiturates MonlabTest® was tested by adding various drugs, drug metabolites, and other compounds that are likely to be present in urine. All compounds were prepared in drug-free normal human urine.

The following structurally related compounds produced positive results when tested at levels equal to or greater than the concentrations listed below.

Compound	Concentration (ng/mL)
Alphenal	300
Aprobarbital	300
Barbital	300
Butabarbital	300
Butethal	300
Pentobarbital	300
Phenobarbital	300
Allobarbital	1,000
Amobarbital	1,000
Butalbital	2,000

The following compounds were found not to cross-react when tested at concentrations up to 100µg/mL.

Acetaminophen	Benzoylcegonine
Acetone	Bilirubin
Albumin	(+)-Brompheniramine
Amitriptyline	Caffeine
D-Amphetamine	Chloroquine
L-Amphetamine	(+/-)-Chlorpheniramine
Ampicillin	Chlorpromazine
Aspartame	Cocaine
Aspirin	Codeine
Atropine	Creatine
Benzocaine	(-)-Desoxyephedrine
Dextromethorphan	Natoxetine
Diazepam	Naltrexone
4-Dimethylaminoantipyrine	(+/-)-Naproxen
Dopamine	(+/-)-Norephedrine
Doxylamine	Nortriptyline
Ecgonine	Oxalic Acid
Ecgonine Methyl Ester	Oxazepam
(+/-)-Epinephrine	Oxycodone
(+)-Epinephrine	Penicillin-G
Erythromycin	Pentermine
Ethanol	Phencyclidine
Foursemide	Pheniramine
Glucose	Phenothiazine
Guaiacol Glyceryl Ether	L-Phenylephrine
Hemoglobin	β-Phenylethylamine
Hydrocodone	Procaine
Hydromorphone	Promethazine
Imipramine	d-Propoxyphene
(+/-)-Isoproterenol	Quinidine
Lidocaine	Ranitidine
Meperidine	Sodium Chloride
Methadone	Sulindac
Methamphetamine	Thioridazine
Methaqualone	Trifluoperazine
(1R, 2S)-(-)-N-Methyl-Ephedrine	Trimethobenzamide
Methylphenidate	Tyramine
Morphine	Vitamine C
11-nor-Δ <sup>9</sup> -Tetrahydrocannabinol-9-carboxylic acid	

## REFERENCES

- Baselt, R, C. Disposition of Toxic Drugs and Chemicals in Man. Biomedical Publications. Davis. CA 1982.
- Urine Testing for Drugs of Abuse. National Institute on Drug Abuse (NIDA). Research Monograph 73, 1986.
- Ellenhorn, M. J. And Barceloux, D. G. Medical Toxicology. Elsevier Science Publishing Company. Inc New York. 1988
- Fed. Register. Department of Health and Human Services. Mandatory Guidelines for Federal Workplace Drug Testing Programs, 53, 69, 11970-11979, 1988

- Gilman, A, G, and Goodman, L, S, The Pharmacological Basis of Therapeutics. eds. MacMillan Publishing. New York NY, 1980.
- Gorodetzky, C, W, Detection of Drugs of Abuse in Biological Fluids. In Martin WR(ed): Drug Addiction 1, New York, Spring-Verlag, 1977.
- Harvey, R, A, Champe, P, C, Lippincott's illustrated Reviews. Pharmacology, 91-95,1992.
- Hofmann F, E, A, Handbook on Drug and Alcohol Abuse: The Biomedical Aspects. New York, Oxford University Press, 1983.
- McBay, A, J, Clin. Chem, 22, 33B-40B, 1987

## PACKAGING

CHEM-7008C-40

40 Barbiturate Tests

## GRAPHICAL SYMBOLS USED

	Manufacturer		For <i>in vitro</i> diagnostic use only
	Don't re-use		Consult instructions for use
	Contains sufficient for <n> tests		Keep dry
	Catalogue Code		Temperature limitation
	Lot Number		Use by

**Manufacturer:** Shanghai Chemtron Biotech Co., Ltd.

Address: No. 518, Qingdai Rd., International

Medical Park, Pudong, 201318

Shanghai, PEOPLE'S REPUBLIC OF CHINA

**European Representative:**

Shanghai International Holding Corp. GmbH (Europe)

Address: Eiffestrasse 80, 20537 Hamburg Germany

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## REFERENCIAS

1. Baselt. R, C. Disposition of Toxic Drugs and Chemicals in Man. Biomedical Publications. Davis. CA 1982.
2. Urine Testing for Drugs of Abuse. National Institute on Drug Abuse (NIDA). Research Monograph 73, 1986.
3. Ellenhorn. M. J. And Barce1oux, D. G. Medical Toxicology. Elsevier Science Publishing Company. 1nc New York. 1988
4. Fed. Register. Department of Health and Human Services. Mandatory Guidelines for Federal Workplace Drug Testing Programs, 53, 69, 11970-11979, 1988
5. Gilman, A, G, and Goodman, L, S, The Pharmacological Basis of Therapeutics. eds. MacMillan Publishing. New York NY, 1980.
6. Gorodetzky. C, W, Detection of Drugs of Abuse in Biological Fluids. In Martin WR (ed): Drug Addiction 1, New York, Spring-Verlag, 1977.
7. Harvey, R, A, Champe, P, C, Lippincotts illustrated Reviews. Pharmacology, 91-95,1992.
8. Hofmann F. E, A, Handbook on Drug and Alcohol Abuse: The Biomedical Aspects. New York, Oxford University Press, 1983.
9. McBay, A, J, Clin. Chem, 22, 33B-40B, 1987

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**Fabricante:** Shanghai Chemtron Biotech Co., Ltd.