One-Step Phencyclidine MonlabTest®

Only for professional *in vitro* diagnostic use. Store at 2 - 30°C.

INTENDED USE

The One-Step Phencyclidine MonlabTest® is a lateral flow, one-step immunoassay for the qualitative detection of phencyclidine in human urine at a cut-off of 25ng/mL. This product is used to obtain a visual, qualitative result and is intended for use by drug testing professionals in drug testing programs. The assay should not be used without proper supervision and is not intended for over the counter sale to lay persons.

This assay provides only a preliminary analytical test result. A more specific alternative chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) has been established as the preferred confirmatory method by the National Institute on Drug Abuse (NIDA). Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are indicated.

SUMMARY

Phencyclidine is an arylcyclohexylamine that was originally used as an anesthetic agent and a veterinary tranquilizer. Phencyclidine can produce hallucinations, lethargy, disorientation, loss of coordination, trance-like ecstatic states, a sense of euphoria and visual distortions. It has many street names, such as "angel dust" and "crystal cyclone", etc. phencyclidine can be administered orally, by nasal ingestion, smoking, or intravenous injection. It is metabolized in the liver and excreted through the kidneys. The half-life of phencyclidine in the body is around three days. Suction and urinary acidification in the treatment of overdose typically reduces its half-life from three days to one day.

Urine based screening tests for drugs of abuse range from simple immunoassay tests to complex analytical procedures. The speed and sensitivity of immunoassays have made them the most widely accepted method for screening urine for drugs of abuse. The One-Step Phencyclidine MonlabTest® is based on the principle of the highly specific immunochemical reactions of antigens and antibodies which are used for the analysis of specific compounds in biological fluids. This test is a rapid, visual, competitive immunoassay that can be used for the qualitative detection of phencyclidine in human urine at a 25ng/mL cut-off concentration.

PRINCIPLE

The One-Step Phencyclidine MonlabTest® is a one-step immunoassay in which a chemically labeled drug (drug conjugate) competes with the drug which may be present in urine for limited antibody binding sites. The test device contains a membrane strip which was pre-coated with drug conjugate on the test band region. A colored anti-phencyclidine monoclonal antibody-colloidal gold conjugate pad is placed at the end of the membrane.

The colored antibody-colloidal gold conjugate moves along with urine, chromatographically by the capillary action, across the membrane, in the absence of drug in the urine, the colored

antibody colloidal gold conjugate attaches to the drug conjugate on the test band region to form a visible line as the antibody/drug conjugate complexes. Therefore, the formation of a visible precipitant in the test band region occurs when the test urine is negative for the drug. When the drug is present in the urine, the drug/metabolite antigen competes with the drug conjugate on the test band region for the limited antibody sites. When a sufficient amount of drug is present, it will fill the limited antibody binding sites. This will prevent attachment of the colored antibody-colloidal gold conjugate to the drug conjugate zone on the test band region. Therefore, absence of the color band on the test region indicates a positive result.

A control or reference band with a different antigen/antibody reaction is also added to the immunochromatographic membrane strip to indicate that the test is performed properly. This control line should always appear regardless of the presence of drug or metabolite. This means that **negative** urine will produce **two** colored bands, and **positive** urine will produce only **one** band. The presence of this colored band in the control region also serves as verification that 1) sufficient volume has been added, and 2) that proper flow was obtained.

MATERIALS SUPPLIED

- 40 individually wrapped test devices which include one disposable pipette each.
- One instruction sheet.

MATERIAL REQUIRED BUT NOT PROVIDED

- Specimen collection container.
- Timer.

STORAGE AND STABILITY

- The test kit should be stored refrigerated or at temperature 2-30°C.
- Each device should remain in its sealed pouch for the duration of the shelf-life.

PRECAUTIONS

- FOR IN VITRO DIAGNOSTIC USE.
- For professional use only.
- Urine specimens may be potentially infectious. Proper handling and disposal methods should be established.
- Avoid cross-contamination of urine samples by using a new specimen collection container and specimen pipette for each urine sample.
- Do not use it after expiration date.

SPECIMEN COLLECTION AND HANDLING



The One-Step Phencyclidine MonlabTest® is formulated for use with urine specimens. Fresh urine does not require any special handling or pretreatment. Urine samples should be collected such that testing can be performed as soon as possible after the specimen collection, preferably during the same day. The specimen may be refrigerated at 2-8°C for 2 days or frozen at -20°C for a longer period of time. Specimens that have been refrigerated must be equilibrated to room temperature prior to testing. Specimens previously frozen must be thawed, equilibrated to room temperature, and mixed thoroughly prior to testing.

Note: Urine specimens and all materials coming in contact with them should be handled and disposed of as if capable of transmitting infection. Avoid contact with skin by wearing gloves and proper laboratory attire.

TEST PROCEDURE


- Review "Specimen collection" instructions. Test device, patient's samples, and controls should be brought to room temperature (10-30°C) prior to testing. Do not open pouches until ready to perform the assay.
- Remove the test device from its protective pouch (bring the device to room temperature before opening the pouch to avoid condensation of moisture on the membrane). Label the device with patient or control number.
- Hold the dropper vertically and transfer 3 full drops of urine (approx. 0.15mL) to the specimen well (S) of the test device, and then start the timer. Use a separate pipette and device for each sample or control.
- Read result between 3 to 8 minutes after the addition of samples.** Do not read result after 8 minutes.

INTERPRETATION OF RESULTS


Positive

The control line appears in the test window, but the test line is not visible.

Negative

Two pink-rose lines (bands) are visible in the control (C) and test (T) areas of the test window. The intensity of the test line may be less than that of the control line; this still means negative result.

Invalid

The test is invalid if the control line is not visible at five minutes. The test failed, or the test procedure was not followed properly. Verify the test procedure and repeat the test with a new testing device.

Note: A very faint line on the test region indicates that the phencyclidine in the sample is near the cut-off level for the test. These samples should be re-tested or confirmed with a more specific method before a positive determination is made.

LIMITATIONS OF PROCEDURE

- The assay is designed for use with human urine only.
- A positive result with any of the tests indicates the presence of a drug/metabolite only and does not indicate or measure intoxication.
- There is a possibility that technical or procedural errors as well other substances or factors not listed may interfere with the test and cause false results. See **SPECIFICITY** for lists of substances that will produce positive results, or that do not interfere with test performance.
- If it is suspected that the samples have been mislabeled or tampered with, a new specimen should be collected and the test should be repeated.

QUALITY CONTROL

Good laboratory practice recommends the use of control materials to ensure proper kit performance. Quality control specimens are available from commercial sources. When testing the positive and negative controls, use the same assay procedure as with a urine specimen.

PERFORMANCE CHARACTERISTICS

Accuracy

The accuracy of Drug Screen PCP was evaluated in comparison to a commercially available immunoassay at a cut-off of 25ng/mL for phencyclidine. One hundred twenty two (122) urine samples, collected from presumed non-user volunteers, have been tested by both procedures with 100% agreement.

In a separate study, fifty (50) urine samples, obtained from a clinical laboratory where they were screened and confirmed as positives by the commercially available immunoassay and GC/MS, were tested with The One-Step Phencyclidine MonlabTest®. Of the 50 samples with phencyclidine concentration from 34.6ng/mL to 1846ng/mL, all were found positive by both methods (100% agreement).

Reproducibility

The reproducibility of The One-Step Phencyclidine MonlabTest® was evaluated at four different sites using blind controls. Of the 60 samples without phencyclidine all were determined negative. Of the (60) samples with phencyclidine concentration of 50ng/mL, all were determined positive.

Precision

The precision of The One-Step Phencyclidine MonlabTest® was determined by conducting the test with spiked controls. The control at the 50ng/mL should give a positive result.

Concentration (ng/mL)	Number Tested	Correct Result	Correct Result %
0	50	50	100
50	50	50	100

Specificity

The specificity for The One-Step Phencyclidine MonlabTest® was tested by adding various drugs, drug metabolites, and other compounds that are likely to be present in urine. All compounds were prepared in drug-free normal human urine.

The following structurally related compounds produced positive results when tested at levels equal to or greater than the concentrations listed below.

Compound	Concentration (ng/mL)
Phencyclidine	25
Tenocyclidine	2,000

The following compounds were found not to cross-react when tested at concentrations up to 100µg/mL.

Acentaminophen	Erythromycin
Acetone	Ethanol
Albumin	Furosemide
Amitriptyline	Glucose
D-Amphetamine	Guaiacal glyceryl ether
L-Amphetamine	Hemoglobin
Ampillicin	Hydrocodone
Aspartame	Hydromorphone
Aspirin	(+/-)-Isoproterenol
Benzocaine	Lidocaine
Benzoylcegonine	Meperidine
Bilirubin	Methamphetamine
(+)-Brompheniramine	Methaqualone
Caffeine	(1R, 2S)-(-)-N-Methyl-ephedrine
Chloroquine	Methylphenidate
(+/-)-Chlorpheniramine	Morphine
Chlorpromazine	Naloxone
Cocaine	Natrexone

Codeine	(+)-Naproxen
Creatine	(+/-)-Norephedrine
(-)-Deoxyephedrine	Notriptyline
Dextromethorphan	Oxalic Acid
Diazepam	Oxazepam
4-Dimethylaminoantipyrine	Oxycodone
Dopamine	Penicillin-G
Doxylamine	Pentemine
Ecgonine	Pentobarbital
Ecgonine methyl ester	Pheniramine
(+/-)-Ephedrine	Pheniramine
(-)-Ephedrine	Phenobarbital
(+)-Epinephrine	Phenothiazine
L-Phenylephrine	Sulindac
β-Phenylethylamine	11-nor-Δ ⁹ -tetrahydrocannabinol-9
	- carboxylic acid
Procaine	Thioridazine
Promethazine	Trifluoperazine
d-Propoxyhene	Tyramine
Ranitidine	Secobarbital
Sodium Chloride	Vitamin C

Manufacturer: Shanghai Chemtron Biotech Co., Ltd.

Address: No. 518, Qingdai Rd., International

Medical Park, Pudong, 201318

Shanghai, PEOPLE'S REPUBLIC OF CHINA

European Representative:

Shanghai International Holding Corp. GmbH (Europe)

Address: Eiffestrasse 80, 20537 Hamburg Germany

Tel: 0049-40-2513175 / Fax: 0049-40-255726

MONLAB

MonlabTest®

Importer /

Distributor:

MONLAB SL

Cobalto, 74

08940 Cornellà de Llobregat (Spain)

Tel +34 934 335 860 / Fax +34 934 363 894

REFERENCES

1. Aniline O, Pittes, F, N. Phencyclidine (PCP). A review and perspectives .CRC Crit. Rex. Toxicol. 1982. 10. 145-177.
2. Baselt. R, C. Disposition of Toxic Drugs and Chemicals in Man. Biomedical Publications. Davis. CA 1982.
3. Urine Testing for Drugs of Abuse. National Institute on Drug Abuse (NIDA). Research Monograph 73, 1986.
4. Ellenhorn. M. J. And Barceloux, D. G. Medical Toxicology. Elsevier Science Publishing Company. Inc New York. 1988
5. Fed. Register. Department of Health and Human Services. Mandatory Guidelines for Federal Workplace Drug Testing Programs, 53, 69, 11970-11979, 1988
6. Gilman, A, G, and Goodman, L, S, The Pharmacological Basis of Therapeutics. eds. MacMillan Publishing. New York NY, 1980.
7. Gorodetzky. C, W, Detection of Drugs of Abuse in Biological Fluids. In Martin WR(ed): Drug Addiction 1, New York, Spring-Verlag, 1977.
8. Harvey, R, A, Champe, P, C, Lippincotts illustrated Reviews. Pharmacology, 91-95,1992.
9. Hofmann F. E, A, Handbook on Drug and Alcohol Abuse: The Biomedical Aspects. New York, Oxford University Press, 1983.
10. McBay, A, J, Clin. Chem, 22, 33B-40B, 1987

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REFERENCIAS

1. Aniline O, Pittes, F, N. Phencyclidine (PCP). A review and perspectives .CRC Crit. Rex. Toxicol. 1982. 10. 145-177.
2. Baselt. R, C. Disposition of Toxic Drugs and Chemicals in Man. Biomedical Publications. Davis. CA 1982.
3. Urine Testing for Drugs of Abuse. National Institute on Drug Abuse (NIDA). Research Monograph 73, 1986.
4. Ellenhorn. M. J. And Barceloux, D. G. Medical Toxicology. Elsevier Science Publishing Company. 1nc New York. 1988
5. Fed. Register. Department of Health and Human Services. Mandatory Guidelines for Federal Workplace Drug Testing Programs, 53, 69, 11970-11979, 1988
6. Gilman, A, G, and Goodman, L, S, The Pharmacological Basis of Therapeutics. eds. MacMillan Publishing. New York NY, 1980.
7. Gorodetzky. C, W, Detection of Drugs of Abuse in Biological Fluids. In Martin WR(ed): Drug Addiction 1, New York, Spring-Verlag, 1977.
8. Harvey, R, A, Champe, P, C, Lippincotts illustrated Reviews. Pharmacology, 91-95,1992.
9. Hofmann F. E, A, Handbook on Drug and Alcohol Abuse: The Bionedical Aspects. New York, Oxford University Press, 1983.
10. McBay, A, J, Clin. Chem, 22, 33B-40B, 1987

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Fabricante: Shanghai Chemtron Biotech Co., Ltd.