One-Step MDMA MonlabTest®

Only for professional *in vitro* diagnostic use. Store at 2 - 30°C.

INTENDED USE

The One-Step MDMA MonlabTest® is a one step, lateral flow chromatographic immunoassay for the qualitative detection of MDMA and its metabolites in human urine. The test cut off is 500ng/mL.

The rapid test provides only a preliminary test result. A more specific alternative chemical method, such as Gas Chromatography/Mass Spectrometry (GC/MS), must be used in order to obtain a confirmed analytical result. The National Institute on Drug Abuse (NIDA) has recommended GC/MS method as the preferred confirmation method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary test result is positive.

SUMMARY

MDMA, 3, 4-methylenedioxy-N-methamphetamine, was developed and patented in the early 1900's as a chemical precursor in the synthesis of pharmaceuticals. Chemically, MDMA is similar to the stimulant amphetamine and the hallucinogen mescaline. MDMA can produce both stimulant and psychedelic effects. MDMA didn't gain notoriety as an illicit drug until the 1980's. By 1985, the United States Drug Enforcement Agency placed MDMA on the schedule I list (Drug with a high potential for abuse with no current medical use). Despite the schedule I status, the recent popularity of MDMA as an illicit drug warrants the need for an accurate drug screen specifically designed to detect MDMA.

MDMA is taken orally, usually in a tablet or a capsule with 80-150mg. MDMA's effects last approximately 3 to 6 hours following oral administration, hyperthermia, though confusion, depression, sleep problems, anxiety, and paranoia have been reported to occur even weeks after the drug is taken. MDMA can produce a significant increase in heart rate and blood pressure and a sense of alertness like that associated with amphetamine use. Following a typical dose, 65% of MDMA is excreted unchanged in the urine and up to 7% is demethylated and eliminated in the urine as methylenedioxyamphetamine (MDA). Others metabolites including conjugated mono and di-hydroxy derivatives of both MDMA and MDA. MDMA is detectable in urine for up to 3 days after use.

One step, urine based drug of abuse screening tests are a rapid, visual, sensitive immunoassay test. They are the most widely accepted methods for screening urine for drug of abuse. The One-Step MDMA MonlabTest® is based on the principle of the high affinity and specific antibody and antigen reaction. It can be performed without the use of an instrument. The test utilizes a monoclonal antibody to selectively detect elevated levels of MDMA in urine. The test will show a positive result when the MDMA(Ecstasy) in urine reaches 500ng/mL.

PRINCIPLE

The One-Step MDMA MonlabTest® is an immunoassay based on the principle of competitive binding. The drug or the metabolites, which may be present in the urine specimen competes against limited antibody binding sites. During testing, a urine specimen migrates upward by capillary action. MDMA, if present in the urine specimen below 500ng/mL, will not saturate the antibody binding sites. The colored antibody-particles conjugate will then be captured by immobilized MDMA-BTG on membrane and a visible colored line will show up in the test line region. The colored line will not form in the test line region if the MDMA level is at or above 500ng/mL, because it will saturate all antibodies's binding sites.

A drug-positive urine specimen will not generate a colored line in the test line region because of drug competition; while a drug negative urine specimen will generate a line in the test line region because of the absence of drug competition. To serve as a procedural control, a colored line will always appear at the control line region if the test has been performed properly.

MATERIALS SUPPLIED

- 40 individually wrapped test which include one disposable pipette each.
- One instruction sheet.

MATERIAL REQUIRED BUT NOT PROVIDED

- Specimen collection container
- Timer

STORAGE AND STABILITY

- Store as packaged in the sealed pouch at 2-30°C.
- The kit must remain in the sealed pouch until use.
- DO NOT FREEZE.
- Do not use it after the expiration date.

PRECAUTIONS

- FOR PROFESSIONAL AND IN VITRO DIAGNOSTIC USE ONLY.
- The test device should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The test device should be discarded in a proper biohazard container after testing
- Avoid cross-contamination of urine samples by changing a new specimen collection container and specimen pipette for each sample.


SPECIMEN COLLECTION AND HANDLING

- The urine specimen must be collected in a clean and dry container. Urine samples should be collected such that testing can be performed as soon as possible after the specimen collection. Urine specimens exhibiting visible precipitates should be centrifuged, filtered, or allowed to settle to obtain a clear supernatant for testing.
- Urine specimens may be stored at 2-8°C for up to 48 hours prior to testing. For prolonged storage, specimens may be frozen and stored below -20°C. Refrigerated samples must be equilibrated to room temperature before testing. Frozen specimens should be thawed and mixed before testing.

TEST PROCEDURE


- Read the instruction sheet carefully before testing. Allow the test device, urine specimen and/or controls to equilibrate to room temperature (10-30°C) prior to testing. Do not open pouches until ready to perform the assay.
- Bring the pouched devices to room temperature before opening to avoid condensation of moisture on the membrane. Remove the test devices from the sealed pouch and label the device with patient or control number.
- Place the test device on a clean and level surface. Hold the dropper vertically and transfer 3 full drops of urine (approx. 0.15mL) to the specimen well (S) of the test device, and then start the timer. Avoid trapping air bubbles in the specimen well (S). See the illustration below.
- Wait for the red line(s) to appear. **The result should be read between 3-8 minutes after the addition of sample.** Do not interpret the result after 8 minutes.

INTERPRETATION OF RESULTS


Negative

Two colored lines should be observed in the viewing window. The line in the test region (T) is the test line; the line in the control region (C) is the control line, which is used to indicate proper performance of the device. The color intensity of the test line may be weaker or stronger than that of the control line.

Positive

Only one colored line appears in the control region (C). The absence of a test line indicates a positive result.

Invalid

The test is invalid if the control line is not visible at five minutes. The test failed, or the test procedure was not followed properly. Verify the test procedure and repeat the test with a new testing device.

NOTE: A very faint line on the test region indicates that the MDMA concentration in the sample is near the cut-off level. These samples should be re-tested or confirmed with a more specific method, such as GC/MS.

LIMITATIONS OF THE PROCEDURE

- The assay is designed for human urine only.
- A positive result with any of the test indicates the presence of a drug/metabolite only and does not indicate or measure intoxication. The test result provides only a preliminary analytical result. A secondary analytical method must be used to obtain a confirmed result. Gas chromatography and mass spectrometry (GC/MS) are the preferred confirmatory methods.
- It is possible that technical or procedural errors, as well as other interfering substances in the urine specimen may cause erroneous results. See **SPECIFICITY** for list of substances that will produce positive results, or that do not interfere with test performance.
- Adulterants, such as bleach and/or alum, in urine specimens may produce erroneous results regardless of the analytical method used. If adulteration is suspected, the test should be repeated with another urine specimen.
- Certain medications containing opiates or opiate derivatives may produce a positive result. Additionally, foods and tea containing poppy products may also produce a positive result.

QUALITY CONTROL

- A procedural control is included in the test. A red line appearing in the control region (C) is the internal procedural control. It confirms sufficient specimen volume and correct procedural technique. A clear background is required.
- Control standards are not supplied with this kit; however it is recommended that positive and negative controls be tested as good laboratory testing practice to confirm the test procedure and to verify proper test performance. When testing the **positive** and **negative** controls, use the same assay procedure as with a urine specimen.

PERFORMANCE CHARACTERISTICS

Accuracy

The accuracy of One-Step MDMA MonlabTest® was evaluated in comparison to a commercially available one-step immunoassay. One hundred (100) urine samples collected from presumed non-user health people have been tested by both procedures with 100% agreement. In a separate study, fifty-three (53) urine samples, obtained from a clinical laboratory where they were confirmed as positives by GC/MS, were tested by both One-Step MDMA MonlabTest® and a commercially available one-step immunoassay. Of the forty-five (45) samples with MDMA concentration between 600-2600ng/mL, all were found to be positives by both methods (100% agreement). Of the eight (8) samples with MDMA concentration between 200-400ng/mL, all were determined negatives by both methods (100% agreement).

Reproducibility

The reproducibility of One-Step MDMA MonlabTest® was evaluated at three different site using blind controls. Of the fifty samples with MDMA concentration of 250ng/mL all were determined negative. Of the fifty samples with MDMA concentration of 1,000ng/mL all were determined positive.

Precision

The precision of the One-Step MDMA MonlabTest® was determined by conducting the test with spiked controls. The control at the 250ng/mL should give a negative result and the control at the 1,000ng/mL should give a positive result.

Concentration (ng/mL)	Test number	Correct result	% Correct result
250	50	50	100
1,000	50	50	100

Specificity

The specificity for One-Step MDMA MonlabTest® was tested by preparing various drug, drug metabolites, and other compounds that are likely to be present in urine. All compounds were prepared in drug-free normal human urine.

The following table lists compounds that are positively detected in urine by One-Step MDMA MonlabTest® at levels equal to or greater than the concentration list below.

Compound	Concentration (ng/mL)
MDMA	500
MDA	2,000

A study was conducted to determine the cross-reactivity of the test with compounds in urine not associated with MDMA. **The following compounds were found not to cross-react when tested at concentrations up to 100µg/mL.**

Acetamidophen	Glucuronide
Acetylsalicylate	Glutethimide
Aminopyrine	Guaifenesin
Amitriptyline	Hippuric acid
Amobarbital	Hyalazine
Amoxapine	Hydrochlorothiazide
Amoxicillin	Hydrocodone
L-Amphetamine	Hydrocortisone
Apomorphine	Hydromorphone
Ascorbic acid	O-Hydroxyhippuric acid
Aspartame	Ibuprofen
Atropine	Imipramine
Benzocaine	Iproniazid
Benzoylcegonine	(-) Isoproterenol
Benzphetamine	Isosuprine

Butabarbital	Ketamine
Cannabidiol	Ketoprofen
Chloralhydrate	Labetalol
Chloramphenicol	Levorphanol
Chlordiazepoxide	Lidocaine
Chlorothiazide	Loperamide
Chlorpromazine	Loxapine succinate
Chlorquine	Maprotiline
Cholesterol	Meperidine
Clomipramine	Mephentermine
Clonidine	Meprobamate
Cocaine	Methadone
Cortisone	Methaqualone
(-) Cotinine	Methoxyphenamine
Creatinine	Methylphenidat
Deoxycorticosterone	Methypylon
Dextromethorphan	Morphine-3-β-D-glucuronide
Diazepam	Nalidixic acid
Diethylpropion	Nalorphine
Diflunisal	Naloxone
Digoxin	Naltrexone
Diphenhydramine	Naproxen
Doxylamine	Noroxymorphone
Ecgonine	Norpseudoephedrine
(-) ψ Ephedrine	Noscapine
β-Estradiol	Nylidrin
Estrone-3-sulfate	D,L-Octopamine
Ethyl-p-aminobenzoate	Oxalic acid
Fenoprofen	Oxazepam
Furoximide	Quinine
Gentisic acid	Ranitidine
Oxalic acid	Salicylic acid
Oxazepam	Secobarbital
Oxolinic acid	Serotonin
Oxycodone	Sulfamethazine
Oxymetazoline	Sulindac
Oxymorphone	Temazepam
p-Hydroxymethamphetamine	Tetracycline
Papaverine	Tetrahydrocortisone
Penicillin-G	Tetrahydrozoline
Pentazocaine	Δ ⁹ -THC
Pentobarbital	Thebaine
Perphenazine	Thiamine
Phencyclidine	Thioridazine
Phenelzine	D,L-Thyroxine
Phendimetrazine	Tolbutamide
Phenobarbital	Tranylcypromine
Phenytoin	Triamterene
L-Phenylethylamine	Trifluoperazine
L-Phenylpropanolamine	Trimethoprim
Prednisolone	Trimipramine
Prednisone	D,L-Tryptophan
Promazine	D,L-Tyrosine
Promethazine	Uric acid
D,L-Propranolol	Verapamil
Propiomazine	Zomepirac
D-Propoxyphene	Promazine
D-Pseudoephedrine	Promethazine

Quinidine	D,L-Propranolol
Ecgonine methylester	Propiomazine
Erythromycin	D-Propoxyphene
D-Norpropoxyp	

REFERENCES


- Baselt, R, C. Disposition of Toxic Drugs and Chemicals in Man. Biomedical Publications. Davis. CA 1982.
- Urine Testing for Drugs of Abuse. National Institute on Drug Abuse (NIDA). Research Monograph 73, 1986.
- Ellenhorn, M. J. And Barceloux, D. G. Medical Toxicology. Elsevier Science Publishing Company. 1nc New York. 1988
- Fed. Register. Department of Health and Human Services. Mandatory Guidelines for Federal Workplace Drug Testing Programs, 53, 69, 11970-11979, 1988
- Gilman, A, G, and Goodman, L, S, The Pharmacological Basis of Therapeutics. eds. MacMillan Publishing. New York NY, 1980.
- Gorodetzky, C, W, Detection of Drugs of Abuse in Biological Fluids. In Martin WR(ed): Drug Addiction 1, New York, Spring-Verlag, 1977.
- Harvey, R, A, Champe, P, C, Lippincotts illustrated Reviews. Pharmacology, 91-95,1992.
- Hofmann F. E, A, Handbook on Drug and Alcohol Abuse: The Biomedical Aspects. New York, Oxford University Press, 1983.
- McBay, A, J, Clin. Chem, 22, 33B-40B,1987

PACKAGING

CHEM-7002C-40

40 MDMA Tests

GRAPHICAL SYMBOLS USED

	Manufacturer		For <i>in vitro</i> diagnostic use only
	Don't re-use		Consult instructions for use
	Contains sufficient for <n> tests		Keep dry
	Catalogue Code		Temperature limitation
	Lot Number		Use by

 **Manufacturer:** Shanghai Chemtron Biotech Co., Ltd.

Address: No. 518, Qingdai Rd., International

Medical Park, Pudong, 201318

Shanghai, PEOPLE'S REPUBLIC OF CHINA

 **European Representative:**

Shanghai International Holding Corp. GmbH (Europe)

Address: Eiffestrasse 80, 20537 Hamburg Germany

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REFERENCIAS

- Baselt, R, C. Disposition of Toxic Drugs and Chemicals in Man. Biomedical Publications. Davis. CA 1982.
- Urine Testing for Drugs of Abuse. National Institute on Drug Abuse (NIDA). Research Monograph 73, 1986.
- Ellenhorn, M. J. And Barceloux, D. G. Medical Toxicology. Elsevier Science Publishing Company, Inc New York. 1988
- Fed. Register. Department of Health and Human Services. Mandatory Guidelines for Federal Workplace Drug Testing Programs, 53, 69, 11970-11979, 1988
- Gilman, A, G, and Goodman, L, S, The Pharmacological Basis of Therapeutics. eds. MacMillan Publishing. New York NY, 1980.
- Gorodetzky, C, W, Detection of Drugs of Abuse in Biological Fluids. In Martin WR(ed): Drug Addiction 1, New York, Spring-Verlag, 1977.
- Harvey, R, A, Champe, P, C, Lippincotts illustrated Reviews. Pharmacology, 91-95,1992.
- Hofmann F, E, A, Handbook on Drug and Alcohol Abuse: The Biomedical Aspects. New York, Oxford University Press, 1983.
- McBay, A, J, Clin. Chem, 22, 33B-40B,1987

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Fabricante: Shanghai Chemtron Biotech Co., Ltd.