One-Step THC MonlabTest®

Only for professional *in vitro* diagnostic use. Store at 2 - 30°C.

INTENDED USE

The One-Step THC MonlabTest® is a lateral flow, one-step immunoassay for the qualitative detection of 11-nor-tetrahydrocannabinol-9-THC-carboxylic acid in human urine at a cut-off of 50ng/mL. This product is used to obtain a visual, qualitative result and is intended for professional use. The assay should not be used without proper supervision and is not intended for over the counter sale to lay persons.

This assay provides only a preliminary analytical test result. A more specific alternative chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) has been established as the preferred confirmatory method by the National Institute on Drug Abuse (NIDA). Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are indicated.

SUMMARY

Marijuana is a hallucinogenic agent derived from the flowering portion of the hemp plant. Smoking is the primary method of use of marijuana/cannabis; higher doses used by abusers produce central nervous system effects, altered mood and sensory perceptions, loss of coordination, impaired short-term memory, anxiety, paranoia, depression confusion, hallucinations and increased heart rate. A tolerance to the cardiac and psychotropic effects can occur, and withdrawal syndrome produces restlessness, insomnia, anorexia and nausea.

When marijuana is ingested, the drug is metabolized by the liver. The primary urinary metabolite of marijuana is 11-nor- Δ^9 -THC-9-carboxylic acid, and its glucuronide. This means that the presence of detected cannabinoids, including the primary carboxyl metabolite, in the urine indicate marijuana/cannabis use. Urine based screening tests for drugs of abuse range from simple immunoassay tests to complex analytical procedures. The speed and sensitivity of immunoassays have made them the most widely accepted method for screening urine for drugs of abuse. The One-Step THC MonlabTest® is based on the principle of highly specific immunochemical reactions of antigens and antibodies which are used for the analysis of specific compounds in biological fluids. This test is a rapid, visual, competitive immunoassay that can be used for the qualitative detection of 11-nor- Δ^9 -THC-9-carboxylic acid in human urine at a 50ng/mL cut-off concentration.

PRINCIPLE

The One-Step THC MonlabTest® is a one-step immunoassay in which a chemically labeled drug (drug conjugate) competes with the drug which may be present in urine for limited antibody binding sites. The test device contains a membrane strip which was pre-coated with drug conjugate on the test band. A colored anti-THC monoclonal antibody-colloidal gold conjugate pad is placed at the right end of the membrane.

In the absence of drug in the urine, the solution of colored antibody-colloidal gold conjugate and urine moves upward,

chromatographically by capillary action, across the membrane. This solution then migrates to the immobilized drug conjugate zone on the test band region. The colored antibody-colloidal gold conjugate then attaches to the drug conjugate to form a visible line as the antibody complexes with the drug conjugate. Therefore, the formation of a visible precipitant in the test zone occurs when the test urine is negative for the drug. When the drug is present in the urine, the drug/metabolite antigen competes with the drug conjugate on the test band region for limited antibody sites on the anti-THC monoclonal antibody-colloidal gold conjugate. When a sufficient concentration of drug is present, it will fill the limited antibody binding sites. This will prevent attachment of the colored antibody-colloidal gold conjugate to the drug-conjugate zone on the test band region. Therefore, absence of the color band on the test region indicates a positive result.

A control band with a different antigen/antibody reaction is also added to the immunochromatographic membrane strip at the control region (C) to indicate that the test has performed properly. This control line should always appear, regardless of the presence of drug or metabolite. This means that **negative** urine will produce **two** colored bands, and **positive** urine will produce only **one** band. The presence of this colored band in the control region also serves as 1) verification that sufficient volume has been added, and 2) that proper flow was obtained.

REAGENTS AND MATERIALS SUPPLIED

- 40 individually wrapped test which include one disposable pipette each.
- One instruction sheet.

MATERIALS REQUIRED BUT NOT PROVIDED

- Specimen collection container
- Timer

STORAGE AND STABILITY

The test kit should be stored at temperature 2-30°C in the sealed pouch for the duration of the shelf-life.

PRECAUTIONS

- FOR IN VITRO DIAGNOSTIC USE.
- For professional use only.
- Urine specimens may be potentially infectious. Proper handling and disposal methods should be established.
- Avoid cross-contamination of urine samples by using a new specimen collection container and specimen pipette for each urine sample.
- Do not use it after the expiration date.

SPECIMEN COLLECTION AND HANDLING

The One-Step THC MonlabTest® is formulated for use with urine specimens. Fresh urine does not require any special handling or pretreatment. Urine samples should be collected such that testing can be performed as soon as possible after the specimen collection, preferably during the same day. The specimen may be refrigerated at 2-8 °C for 2 days, or frozen at -20 °C for a longer period of time. Specimens that have been refrigerated must be equilibrated to room temperature prior to testing.


Specimens previously frozen must be thawed, equilibrated to room temperature, and mixed thoroughly prior to testing.

Note: Urine specimens and all materials coming in contact with them should be handled and disposed of as if capable of transmitting infection. Avoid contact with skin by wearing gloves and proper laboratory attire.

TEST PROCEDURE


- Review "Specimen collection" instructions. Test device, patient's samples, and controls should be brought to room temperature (10-30°C) prior to testing. Do not open pouches until ready to perform the assay.
- Remove the test device from its protective pouch (bring the device to room temperature before opening the pouch to avoid condensation of moisture on the membrane). Label the device with patient or control number.
- Hold the dropper vertically and transfer 3 full drops of urine (approx. 0.15ml) to the specimen well (S) of the test device, and then start the timer. Use a separate pipette and device for each sample or control.
- **Read result between 3 to 8 minutes after the addition of samples.** Do not read result after 8 minutes.

INTERPRETATION OF RESULTS


Negative:

Two colored lines appear in the viewing window. The line in the test region (**T**) is the drug probe line; the line on the control region (**C**) is the control line, which is used to indicate proper performance of the device. **The color intensity of the test line may be weaker or stronger than that of the control line.**

Positive

Only **one** colored line appears in the control region (**C**). The **absence** of a test line indicates a positive result.

Invalid

The test is invalid if the control line is not visible at five minutes. The test failed, or the test procedure was not followed properly. Verify the test procedure and repeat the test with a new testing device.

Note: A very faint line in the test region indicates that the THC in the samples is near to the cut-off level for the test. These should be re-tested or confirmed with a more specific method before a positive determination is made.

LIMITATIONS OF PROCEDURE

- The assay is designed for use with human urine only.
- A positive result with the tests indicates the presence of a drug/metabolite only and does not indicate or measure intoxication.
- There is a possibility that technical or procedural errors as well as other substances and factors not listed may interfere with the test and cause false results. See SPECIFICITY for lists of substances that will produce positive results, or that do not interfere with test performance.

- If it is suspected that the samples have been mislabeled or tampered with, a new specimen should be collected and the test should be repeated.

QUALITY CONTROL

Good laboratory practice recommends the use of control materials to ensure proper kit performance. Quality control specimens are available from commercial sources. When testing the positive and negative controls, use the same assay procedure as with a urine specimen.

PERFORMANCE CHARACTERISTICS

Accuracy

The One-Step THC MonlabTest® was evaluated in comparison to a commercially available one-step immunoassay at a cut-off 50ng/mL. One hundred (100) urine samples, collected from presumed non-user volunteers, have been tested by both procedures with 100% agreement.

In a separate study, sixty two (62) urine samples, obtained from a clinical laboratory where they were screened and confirmed as positive by the commercially available immunoassay and GC/MS, were tested by One-Step THC MonlabTest®

Test. The concentration of 11-nor- Δ^9 -THC-9-carboxylic acid in the urine samples were ranged from 30 to 300ng/mL. Of the fifteen (15) samples with 11-nor- Δ^9 -THC-9-carboxylic acid concentrations from 60 to 80ng/mL, all were found to be positives by both procedures (100% agreement). Of the thirty two (32) samples with 11-nor- Δ^9 -THC-9-carboxylic acid concentration ≥ 80 ng/mL, all were found to be positives by both procedures (100% agreement). Of the fifteen (15) samples with 11-nor- Δ^9 -THC-9-carboxylic acid concentration from 33 to 42ng/mL, all were found to be negatives by both procedures (100% agreement).

Reproducibility

The reproducibility of the One-Step THC MonlabTest® was evaluated at four different sites using blind controls. Of the sixty samples with a carboxylic acid concentration of 25ng/mL, all were determined negatives. Of the sixty samples with carboxylic acid concentration of 100ng/mL, all were determined positive.

Precision

The precision of the One-Step THC MonlabTest® was determined by conducting the test with spiked controls. The control at the 25ng/mL should give a negative result and the control at the 75ng/mL should give a positive result.

11-nor Δ^9 -THC-9-carboxylic acid Conc. (ng/mL)	Number Tested	Control Result	Correct Result
25	60	60	100
75	60	60	100

Specificity

The specificity for the One-Step THC MonlabTest® was tested by adding various drugs, drug metabolites, and other compounds that are likely to be present in urine. All compounds were prepared in drug-free normal human urine.

The flowing structurally related compounds produced positive results when tested at levels equal to or greater than the concentrations listed below.

Compound	Concentration (ng/mL)
11-nor- Δ^8 -THC-9-carboxylic acid	50
11-nor- Δ^9 -THC-9-carboxylic acid	50
11-hydroxy- Δ^9 -tetrahydrocannabinol	2,500
Δ^8 -tetrahydrocannabinol	7,500
Δ^9 -tetrahydrocannabinol	10,000
Cannabinol	10,000

The following compounds were found not to cross-react when tested at concentrations up to 100 μ g/mL.

Acetaminophen	Guaiacol
Acetone	Glycerol
Albumin	Ether
Amitriptyline	Hemoglobin
Ampicillin	Imipramine
Aspartame	(+/-)-Isoproterenol
Aspirin	Lidocaine
Atropine	(1R,S)-(-)-N-Methyl-Ephedrine
Benzocaine	(+)-Naproxen
Bilirubin	(+/-)-Norephedrine
Caffeine	Oxalic Acid
Chloroquine	Penicillin-G
(+)-Chlorpheniramine	Pheniramine
(+/-)-Chlorpheniramine	Phenothiazine
Creatine	I-Phenylephrine
Dexbrompheniramine	Phenylethylamine
Dextromethorphan	Procaine
4-Dimethylaminoantipyrine	Quinidine
Dopamine	Ranitidine
(+/-)-Ephedrine	Riboflavin
(-)-Ephedrine	Sodium Chloride
(+)-Epinephrine	Sulindac
Erythromycin	Tyramine
Ethanol	Vitamin C
Furosemide	
Glucose	

REFERENCES

- Baselt, R, C. Disposition of Toxic Drugs and Chemicals in Man. Biomedical Publications. Davis. CA 1982.
- Urine Testing for Drugs of Abuse. National Institute on Drug Abuse (NIDA). Research Monograph 73, 1986.
- Ellenhorn, M. J. And Barce1oux, D. G. Medical Toxicology. Elsevier Science Publishing Company. 1nc New York. 1988
- Fed. Register. Department of Health and Human Services. Mandatory Guidelines for Federal Workplace Drug Testing Programs, 53, 69, 11970-11979, 1988
- Gilman, A, G, and Goodman, L, S, The Pharmacological Basis of Therapeutics. eds. MacMillan Publishing. New York NY, 1980.
- Gorodetzky, C, W, Detection of Drugs of Abuse in Biological Fluids. In Martin WR(ed): Drug Addiction 1, New York, Spring-Verlag, 1977.
- Harvey, R, A, Champe, P, C, Lippincotts illustrated Reviews. Pharmacology, 91-95,1992.



- Hofmann F. E, A, Handbook on Drug and Alcohol Abuse: The Biomedical Aspects. New York, Oxford University Press, 1983.
- McBay, A, J, Clin. Chem, 22, 33B-40B,1987

PACKAGING

CHEM-7005C-40

40 THC Tests

GRAPHICAL SYMBOLS USED

	Manufacturer		For <i>in vitro</i> diagnostic use only
	Don't re-use		Consult instructions for use
	Contains sufficient for <n> tests		Keep dry
	Catalogue Code		Temperature limitation
	Lot Number		Use by

 **Manufacturer:** Shanghai Chemtron Biotech Co., Ltd.

Address: No. 518, Qingdai Rd., International

Medical Park, Pudong, 201318

Shanghai, PEOPLE'S REPUBLIC OF CHINA

 **European Representative:**

Shanghai International Holding Corp. GmbH (Europe)

Address: Eiffestrasse 80, 20537 Hamburg Germany

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REFERENCIAS

- Baselt, R, C. Disposition of Toxic Drugs and Chemicals in Man. Biomedical Publications. Davis. CA 1982.
- Urine Testing for Drugs of Abuse. National Institute on Drug Abuse (NIDA). Research Monograph 73, 1986.
- Ellenhorn, M. J. And Barce1oux, D. G. Medical Toxicology. Elsevier Science Publishing Company. Inc New York. 1988
- Fed. Register. Department of Health and Human Services. Mandatory Guidelines for Federal Workplace Drug Testing Programs, 53, 69, 11970-11979, 1988
- Gilman, A, G, and Goodman, L, S, The Pharmacological Basis of Therapeutics. eds. MacMillan Publishing. New York NY, 1980.
- Gorodetzky, C, W, Detection of Drugs of Abuse in Biological Fluids. In Martin WR(ed): Drug Addiction 1, New York, Spring-Verlag, 1977.
- Harvey, R, A, Champe, P, C, Lippincotts illustrated Reviews. Pharmacology, 91-95,1992.
- Hofmann F. E, A, Handbook on Drug and Alcohol Abuse: The Biomedical Aspects. New York, Oxford University Press, 1983.
- McBay, A, J, Clin. Chem, 22, 33B-40B, 1987

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**Fabricante:** Shanghai Chemtron Biotech Co., Ltd.