
One-Step Methamphetamine 
MonlabTest® 

Only for professional *in vitro* diagnostic use. Store at 2 - 30°C.

INTENDED USE

The One-Step Methamphetamine MonlabTest® is a lateral flow, one-step immunoassay for the qualitative detection of methamphetamine in human urine at a cut-off of 1,000ng/mL. This product is used to obtain a visual, qualitative result and is intended for professional use. The assay should not be used without proper supervision and is not intended for over the counter sale to lay persons.

This assay provides only a preliminary analytical test result. A more specific alternative chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) has been established as the preferred confirmatory method by the National Institute on Drug Abuse (NIDA). Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are indicated.

SUMMARY

Methamphetamine is a potent sympathomimetic agent with therapeutic applications. Acute higher doses lead to enhanced stimulation of the central nervous system and induce euphoria, alertness, and sense of increased energy and power. More acute responses produce anxiety, paranoia, psychotic behavior, and cardiac dysrhythmias. The pattern of psychosis which may appear at high doses may be indistinguishable from schizophrenia.

Methamphetamine is excreted in urine as amphetamine and oxidized as deaminated derivatives. However, 40% of methamphetamine is excreted unchanged. Thus the presence of the parent compound in the urine indicates methamphetamine use.

PRINCIPLE

The One-Step Methamphetamine MonlabTest® is a one-step immunoassay in which a chemically labeled drug (drug-protein conjugate) competes with the drug which may be present in urine for limited antibody binding sites. The test device contains a membrane strip which has been pre-coated with drug-protein conjugate on the test band region. A colored anti-methamphetamine antibody-colloidal gold conjugate pad is placed at the end of the membrane. The colored antibody colloidal gold conjugate moves along with urine, chromatographically by capillary action, across the membrane. In the absence of drug in the urine, the solution of colored antibody-colloidal gold conjugate and urine moves upward chromatographically by capillary action across the membrane. This solution then migrates to the immobilized drug-protein conjugate zone to form a visible line as the antibody complexes with the drug-protein conjugate. Therefore, the formation of a visible precipitant in the test band region occurs when the test urine is negative for the drug.

When the drug is present in the urine, the drug/metabolite antigen competes with the drug-protein conjugate on the test band region for the limited antibody sites on the anti-methamphetamine antibody-colloidal gold conjugate. When a sufficient concentration of drug is present, it will fill the limited

antibody binding sites. This will prevent attachment of the colored antibody-colloidal gold conjugate to the drug-protein conjugate zone on the test band region (T). Therefore, absence of the color band on the test region indicates a positive result.

A control band that has a different antigen/antibody reaction is added to the immunochromatographic membrane strip at the control region (C) to indicate that the test has performed properly. This control line should always appear regardless of the presence of drug or metabolite. This means that negative urine will produce two colored bands, and positive urine will produce only one band. The presence of this colored band in the control region also serves as verification that 1) sufficient volume has been added, and 2) that proper flow was obtained.

MATERIALS SUPPLIED

- 40 individually wrapped test which include one disposable pipette each.
- One instruction sheet.

MATERIAL REQUIRED BUT NOT PROVIDED

- Specimen collection container.
- Timer

STORAGE AND STABILITY

- The test kit should be stored refrigerated or at temperature 2-30°C.
- Each device should remain in its sealed pouch for the duration of the shelf-life.

PRECAUTIONS

- FOR IN-VITRO DIAGNOSTIC USE.
- For professional use only.
- Urine specimens may be potentially infectious. Proper handling and disposal methods should be established.
- Avoid cross-contamination of urine samples by using a new specimen collection container and specimen pipette for each urine sample.
- Do not use it after the expiration date.

SPECIMEN COLLECTION AND HANDLING

The One-Step Methamphetamine MonlabTest® is formulated for use with urine specimens. Fresh urine does not require any special handling or pretreatment. Urine samples should be collected such that testing may be performed as soon as possible after the specimen collection, preferably during the same day. The specimen may be refrigerated at 2-8°C for 2 days or frozen at -20°C for a longer period of time. Specimens that have been refrigerated must be equilibrated to room temperature prior to testing. Specimens previously frozen must be thawed, equilibrated to room temperature, and mixed thoroughly prior to testing.


Note: Urine specimens and all materials coming in contact with them should be handled and disposed of as if capable of transmitting infection. Avoid contact with skin by wearing gloves and proper laboratory attire.

TEST PROCEDURE

- Review "Specimen collection" instructions. Test device, patient's samples, and controls should be brought to room temperature (10-30°C) prior to testing. Do not open pouches until ready to perform the assay.
- Remove the test device from its protective pouch (bring the device to room temperature before opening the pouch to avoid condensation of moisture on the membrane). Label the device with patient or control number.


- Hold the dropper vertically and transfer 3 full drops of urine (approx. 0.15ml) to the specimen well (S) of the test device, and then start the timer. Use a separate pipette and device for each sample or control.
- **Read result between 3 to 8 minutes after the addition of samples.** Do not read result after 8 minutes.

INTERPRETATION OF RESULTS


Positive

The control line appears in the test window, but the test line is not visible.

Negative

Two colored lines should be observed in the viewing window. The line in the test region (T) is the drug probe line; the line in the control region (C) is the control line, which is used to indicate proper performance of the device. The color intensity of the test line may be weaker or stronger than that of the control line.

Invalid

The test is invalid if the control line is not visible at five minutes. The test failed, or the test procedure was not followed properly. Verify the test procedure and repeat the test with a new testing device.

Note: A very faint line on the test region indicates that the methamphetamine in the sample is near the cut-off level for the test. These samples should be re-tested or confirmed with a more specific method before a positive determination is made.

LIMITATIONS OF PROCEDURE

- The assay is designed for use with human urine only.
- A positive result with any of the tests indicates the presence of a drug/metabolite only and does not indicate or measure intoxication.
- There is a possibility that technical or procedural errors as well other substances or factors not listed may interfere with the test and cause false results. See **SPECIFICITY** for lists of substances that will produce positive results, or that do not interfere with test performance.
- If it is suspected that the samples have been mislabeled or tampered with, a new specimen should be collected and the test should be repeated.

QUALITY CONTROL

Good laboratory practice recommends the use of control materials to ensure proper kit performance. Before using a new kit with patient specimens, positive and negative controls should be tested. Quality control specimens are available from commercial sources. When testing the positive and negative controls, use the same assay procedure as with a urine specimen.

PERFORMANCE CHARACTERISTICS

Accuracy

The accuracy of the One-Step Methamphetamine MonlabTest® was evaluated in comparison to a commercially available immunoassay at a cut-off of 1,000ng/mL for methamphetamine. One hundred (100) urine samples, collected from presumed non-user volunteers, have been tested by both procedures with 100% agreement.

In a separate study, fifty one (51) urine samples, (obtained from a clinical laboratory where they were screened and confirmed as positives by GC/MS) were tested with the One-Step Methamphetamine MonlabTest® and the commercially available immunoassay. The concentration of methamphetamine in the urine samples were all above 1,300ng/mL. All the samples were found positive by both methods (100% agreement).

Reproducibility

The reproducibility of the One-Step Methamphetamine MonlabTest® was evaluated at four different sites using blind controls. Of the sixty (60) samples with methamphetamine concentration of 500ng/mL, all were determined negative. Of the sixty samples with methamphetamine concentration of 2,000ng/mL, all were determining positive.

Precision

The precision of the One-Step Methamphetamine MonlabTest® was determined by conducting the test with spiked controls. The control at 500ng/mL gave negative results and the control at 1,500ng/mL gives positive results.

Methamphet-amineConc.	Tested	Positive (+)	Borderline (+/-)	Negative (-)	Correct Results(%)
500ng/mL	40	0	0	40	100
1500ng/mL	40	40	0	0	100

Specificity

The specificity for the One-Step Methamphetamine MonlabTest® was tested by adding various drugs, drug metabolites, and other compounds that are likely to be present in urine. All compounds were prepared in drug-free normal human urine.

The following structurally related compounds produced positive results when tested at levels equal to or greater than the concentrations listed below.

Compound	Concentration (ng/mL)
(+)-Methamphetamine	1,000
D-Amphetamine	50,000
Chloroquine	50,000
(+/-)-Ephedrine	50,000
(-)-Methamphetamine	25,000
(+/-)-3,4-methylenedioxymethamphetamine	2,000
Procaine	10,000
β-Phenylethylamine	50,000
Ranitidine	50,000

The following compounds were found not to cross-react when tested at concentrations up to 100µg/mL.

Acetaminophen	Dexbrompheniramine
Acetone	Dextromethorphan
Albumin	Diazepam
L-Amphetamine	4-Dimethylaminoantipyrine
Ampicillin	Dopamine
Aspartame	Doxylamine

Aspirin	Ecgonine
Atropine	Ecgonine methyl ester
Benzocaine	(+)-Epinephrine
Benzoyllecgonine	Erythromycin
Bilirubin	Ethanol
(+)-Brompheniramine	Furosemide
Caffeine	Glucose
(+)-Chlorpheniramine	Guaiacol Glyceryl Ether
(+/-)-Chlorpheniramine	Hemoglobin
Chlorpromazine	Hydrocodone
Creatine	Hydromaphone
Imipramine	Pentamine
(+/-)-Isoproterenol	Phencyclidine
Lidocaine	Pheniramine
Meperidine	Phenobarbital
Methaqualone	Phenothiazine
(1R,2S)-(-)-N-Methyl-Ephedrine	L-Phenylephrine
Methylphenidate	Promethazine
Morphine	Quinidine
Naltrexone	Riboflavin
(+)-Naproxen	Secobarbital
(+/-)-Norephedrine	Sodium Chloride
Nortriptyline	Sulindac
Oxalic Acid	Thioridazine
Oxazepam	Trifluoperazine
Oxycodone	Tyramine
Penicillin-G	Vitamin C
Oxycodone	

REFERENCES


1. Baselt. R, C. Disposition of Toxic Drugs and Chemicals in Man. Biomedical Publications. Davis. CA 1982.
2. Urine Testing for Drugs of Abuse. National Institute on Drug Abuse (NIDA). Research Monograph 73, 1986.
3. Ellenhorn. M. J. And Barceloux, D. G. Medical Toxicology. Elsevier Science Publishing Company. 1nc New York. 1988
4. Fed. Register. Department of Health and Human Services. Mandatory Guidelines for Federal Workplace Drug Testing Programs, 53, 69, 11970-11979, 1988
5. Gilman, A, G, and Goodman, L, S, The Pharmacological Basis of Therapeutics. eds. MacMillan Publishing. New York NY, 1980.
6. Gorodetzky. C, W, Detection of Drugs of Abuse in Biological Fluids. In Martin WR (ed): Drug Addiction 1, New York, Spring-Verlag, 1977.
7. Harvey, R, A, Champe, P, C, Lippincotts illustrated Reviews. Pharmacology, 91-95,1992.
8. Hofmann F. E, A, Handbook on Drug and Alcohol Abuse: The Biomedical Aspects. New York, Oxford University Press, 1983.
9. McBay, A, J, Clin. Chem, 22, 33B-40B, 1987.

PACKAGING

CHEM-7004C-40

40 Methamphetamine Tests

GRAPHICAL SYMBOLS USED

	Manufacturer		For <i>in vitro</i> diagnostic use only
	Don't re-use		Consult instructions for use
	Contains sufficient for <n> tests		Keep dry
	Catalogue Code		Temperature limitation
	Lot Number		Use by

 **Manufacturer:** Shanghai Chemtron Biotech Co., Ltd.

Address: No. 518, Qingdai Rd., International

Medical Park, Pudong, 201318

Shanghai, PEOPLE'S REPUBLIC OF CHINA

 **European Representative:**

Shanghai International Holding Corp. GmbH (Europe)

Address: Eiffestrasse 80, 20537 Hamburg Germany

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REFERENCIAS

- Baselt, R, C. Disposition of Toxic Drugs and Chemicals in Man. Biomedical Publications. Davis. CA 1982.
- Urine Testing for Drugs of Abuse. National Institute on Drug Abuse (NIDA). Research Monograph 73, 1986.
- Ellenhorn, M. J. And Barceloux, D. G. Medical Toxicology. Elsevier Science Publishing Company. Inc New York. 1988
- Fed. Register. Department of Health and Human Services. Mandatory Guidelines for Federal Workplace Drug Testing Programs, 53, 69, 11970-11979, 1988
- Gilman, A, G, and Goodman, L, S, The Pharmacological Basis of Therapeutics. eds. MacMillan Publishing. New York NY, 1980.
- Gorodetzky, C, W, Detection of Drugs of Abuse in Biological Fluids. In Martin WR (ed): Drug Addiction 1, New York, Spring-Verlag, 1977.
- Harvey, R, A, Champe, P, C, Lippincotts illustrated Reviews . Pharmacology, 91-95,1992.
- Hofmann F. E, A, Handbook on Drug and Alcohol Abuse: The Biomedical Aspects. New York, Oxford University Press, 1983.
- McBay, A, J, Clin. Chem, 22, 33B-40B, 1987.

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**Fabricante:** Shanghai Chemtron Biotech Co., Ltd.