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
Quantitative determination of anti-streptolysin O (ASO)

Only for professional *in vitro* diagnostic use.
 Store at 2-8°C.

PRODUCT CHARACTERISTICS

The calibrator is a lyophilised serum from human origin containing an ASO concentration suitable to calibrate ASO determinations by turbidimetry method.

REAGENTS

ASO Calibrator	Calibrator. Human serum. ASO concentration is stated on the vial label.
-----------------------	---

PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

CALIBRATION

The concentration value of the Calibrator has been standardized against the ASO International Standard from NIBSC ASO.

PREPARATION

Reconstitute (→) with 1.0 mL of distilled water. Mix thoroughly, avoiding foam forming. Bring to room temperature for about 10 minutes before use.

STORAGE AND STABILITY

The calibrator is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

After reconstitution, stable for 1 month at 2-8°C or 3 months at -20°C.

PROCEDURE

To be used in turbidimetric assays.

PRESENTATION

MO-165046	1 x 1 mL
-----------	----------