Calibrador Lp(a) MonlabTest®

Turbidimetría

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Lp(a) Calibrator MonlabTest®

Turbidimetry

Quantitative determination of Lipoprotein (a) [Lp(a)]

Only for professional *in vitro* diagnostic use.
Store at 2-8°C.

PRODUCT CHARACTERISTICS

The calibrator is a lyophilised serum from human origin containing a Lp(a) concentration suitable to calibrate Lp(a) determinations by turbidimetry method.

REAGENTS

Calibrator Lp(a)	Calibrator. Human serum. Sodium azide 0.95 g/L. Lp(a) concentration is stated on the vial label.
-------------------------	--

PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However, handle cautiously as potentially infectious.

CALIBRATION

The concentration value of the Calibrator has been standardized against an Internal Reference Material.

PREPARATION

Reconstitute (→) with 1.0 mL of distilled water.
Mix gently, avoiding foam forming. Bring to room temperature for 15 minutes before use.

STORAGE AND STABILITY

The calibrator is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

After reconstitution, stable for 2 weeks at 2-8°C.

PROCEDURE

To be used in Lp(a) turbidimetric assays.

PRESENTATION

MO-165051

1 x 1 mL