PCR Calibrador

MonlabTest®

Turbidimetría

IVD

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Turbidimetry

IVD

Quantitative determination of C-Reactive Protein (CRP)

Only for professional *in vitro* diagnostic use.
Store at 2-8°C.

PRODUCT CHARACTERISTICS

The calibrator is a lyophilised serum from human origin containing a CRP concentration suitable to calibrate CRP determinations by turbidimetry method.

REAGENTS

CRP-Calibrator	Calibrator. Human serum. CRP concentration is stated on the vial label.
-----------------------	--

PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However, handle cautiously as potentially infectious.

CALIBRATION

The concentration value of the Calibrator has been standardized against the Reference Material ERM-DA 474/IFCC.

PREPARATION

Reconstitute (→) with 1.0 mL of distilled water. Mix thoroughly, avoiding foam forming. Bring to room temperature for about 10 minutes before use.

STORAGE AND STABILITY

The calibrator is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

After reconstitution, stable for 1 month at 2-8°C or 3 months at -20°C.

PROCEDURE

To be used in turbidimetric assays.

PRESENTATION

MO-165048

1 x 1 mL