ASO/PCR/FR Control L
MonlabTest®

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MO-165050 1 x 1 mL

ASO/PCR/FR Control L
MonlabTest®

Quantitative determination of ASO/CRP/RF

Only for professional *in vitro* diagnostic use.
 Store at 2-8°C.

PRODUCT CHARACTERISTICS

The control is a lyophilised human serum used to evaluate the precision and accuracy of ASO, CRP and RF determinations in human serum by turbidimetry method.

REAGENTS

ASO/CRP/RF Control L	Human serum. With a low concentration of ASO, CRP and RF.
----------------------	---

PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However, handle cautiously as potentially infectious.

CALIBRATION

The concentration values of the Control have been standardized against the Reference Material ERM-DA 474/IFCC for CRP, International Reference Standard from NIBSC 64/002 for RF and International Standard from NIBSC ASO for ASO.

PREPARATION

Reconstitute (→) accurately with 1.0 mL of distilled water. Mix thoroughly, avoiding foam forming. Bring at room temperature for about 10 minutes before use.

STORAGE AND STABILITY

The control is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

After reconstitution, stable for 10 days at 2-8°C or 6 months at -20°C.

PROCEDURE

To be used in turbidimetric assays.

COMPONENT	METHOD	VALUE	RANGE	UNITS
Anti streptolysin O (ASO)	Turbidimetry			IU/mL
C- Reactive Protein (CRP)				mg/L
Rheumatoid factor (RF)				IU/mL

Lot:

PRESENTATION

MO-165050 1 x 1 mL