Control Lipídico MonlabTest®

Suero Control Humano

Control Lipídico

Para uso profesional de diagnóstico *in vitro*.
Conservar a 2-8°C.

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Human Control Serum

Lipidic Control

Only for professional *in vitro* diagnostic use.
Store at 2-8°C.

PRODUCT CHARACTERISTICS

Lipidic Control MonlabTest is a human lyophilised control serum. Is intended for control of accuracy for use with manual and automated analytical procedures.

REAGENTS

Human serum. Biological additives. Bacteriostatics agents. The concentration of the components are lot-specific. See chart below for the values of ingredients.

PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

CALIBRATION

Information available on request.

PREPARATION

- Reconstitute (→) with 1.0 mL of distilled water.
- Mix thoroughly, avoiding foam forming.
- Bring to room temperature for about 30 min. before use.

Improper handling and/or storage can affect results. Inaccurate reconstitution and errors in assay technique can cause erroneous results.

STORAGE AND STABILITY

The Control serum is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date or if there is visible evidence of microbial growth.

Store tightly capped when not use.

- After reconstitution is stable for:
- At 15°C to 25°C 30 hours
- At 2°C to 8°C 14 days
- At -25°C to -15°C 3 months

PRESENTATION

MO-165112

1 x 1 mL

Component	Method	Value	Range	Units
Cholesterol	CHOD-PAP			mmol/l
HDL cholesterol	Direct- Liquid			mmol/l
LDL cholesterol	Direct- Liquid			mmol/l
Phosfolipids	CHOLIN-PAP			mmol/l
Triglycerides	GPO-PAP			mmol/l