Multicontrol de Proteínas MonlabTest®

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Lot:

Multicontrol Protein MonlabTest®

Quantitative determination of Plasma Proteins

Only for professional *in vitro* diagnostic use.
Store at 2-8°C.

PRODUCT CHARACTERISTICS

The Multicontrol Protein Serum is a liquid serum from human origin containing protein concentrations suitable to evaluate the accuracy and precision of plasma protein determinations by immunoturbidimetry method.

REAGENTS

CONTROL PROTEIN SERUM	Control. Human serum. Sodium azide 0.95 g/L. Protein concentrations are stated on the list below.
-----------------------	---

PRECAUTIONS

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However, handle cautiously as potentially infectious.

CALIBRATION

The concentration values of the Control are traceable to Standard Reference Material ERM-DA470k/IFCC (all constituents except antithrombin III).

STORAGE AND STABILITY

The control is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.
Opened vials, stable for 30 days at 2-8°C.

PROCEDURE

To be used in immunoturbidimetric assays.

PRESENTATION

MO-165045

1 x 2 mL

PROTEINS	METHOD	CONCENTRATION	RANGE
α ₁ -antitrypsin	Inmunoturbidimetría	62.1	49.7 – 74.5
α ₁ -ac. glycoprotein		39.8	31.8 – 47.8
Antithrombin III		13.8	11.0 – 16.6
C3		70.5	56.4 – 84.6
C4		14.3	11.4 – 17.2
Ceruloplasmin		18.1	14.5 – 21.7
Haptoglobin		53.5	42.8 – 64.2
IgA		120	96.0 – 144
IgG		612	490 – 734
IgM		92.5	74.0 – 111
Prealbumin		18.6	14.9 – 22.3
Transferrin		126	101 – 151

Lot: