## VDRL Estabilizado MonlabTest®

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1. George P. Schmid. Current Opinion in Infectious Diseases 1994; 7: 34-40.
2. Sandra A Larsen et al. Clinical Microbiology Reviews 1995; 8 (1): 1-21.
3. Sandra Larsen et al. A manual of Test for Syphilis American Public Health 4. Association 1990: 1-192.
4. Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACCC Press, 1995.

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Ref.: MO-165023 250 tests	5 mL Antígeno VDRL estabilizado 1 mL Control + 1 mL Control -
---------------------------	---

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## VDRL stabilized MonlabTest®

### Qualitative determination of plasma reagins

Only for professional *in vitro* diagnostic use. Store at 2 - 8°C.

### PRINCIPLE OF THE METHOD

The VDRL MonlabTest is a non-treponemal slide agglutination test for the qualitative and semi-quantitative detection of plasma reagins. The antigen suspension, a lipid complex, is agglutinated when mixed with samples containing reagins of patient affected by syphilis.

### CLINICAL SIGNIFICANCE

Reagins are a group of antibodies against some components produced in the damage tissues from patients infected by *Treponema palladium*, the agent which causes the syphilis. This microorganism produces some damage to the liver and heart, releasing some tissue fragments. Immunological patient system reacts producing reagins, antibodies against these fragments. The assay is useful to follow the antibiotic therapy answer.

### REAGENTS

<b>VDRL antigen stabilized</b>	Solution containing cardiolipin 0.3 g/L, lecithin 2.1g/L and cholesterol 9 g/L in phosphate buffer 1.5 mmol/L. Preservative, pH 7.0.
<b>Control + Red cap</b>	Artificial serum with a reagin titer $\geq 1/8$ .
<b>Control - Blue cap</b>	Animal serum. Preservative

### PRECAUTIONS

Control +/- : H317-May cause an allergic skin reaction. Contains 2-Methylisothiazol-3(2H)-one (Proclin 950). Follow the precautionary advice indicated on the SDS and product label. Components from human origin have been tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However, handle cautiously as potentially infectious.

### CALIBRATION

The sensitivity is calibrated against the International Reference WHO (1<sup>st</sup> International Standard for Syphilitic Human Serum, ref. 05/132).

### PREPARATION

The reagents are ready to use.

### STORAGE AND STABILITY

All the kit components will remain stable until the expiration date printed on the label, when stored tightly closed at 2-8°C and contaminations are prevented during their use. Mix reagents gently before use.

**Do not freeze.** The freezing of VDRL antigen may cause a loss of its functionality.

### ADDITIONAL EQUIPMENT

- Mechanical rotator with adjustable speed at 180 r.p.m.
- Glass slides
- Light microscope (10x objective lens)
- Pipettes 50  $\mu$ L.

### SAMPLES

Fresh serum or plasma. Stable 7 days at 2-8°C or three months at -20°C.

The samples with presence of fibrin should be centrifuged before use. Do not use highly hemolyzed or lipemic samples.

### PROCEDURE

#### Qualitative method

1. Allow the reagents and samples to reach room temperature. The sensitivity of the test may be reduced at low temperatures.
2. Place 50  $\mu$ L of the sample and one drop of each Positive and Negative controls into separate circles on the slide test.
3. Swirl the VDRL suspension gently before using and add 20  $\mu$ L of this reagent onto each sample.
4. Place the slide on a mechanical rotator at 160-180 r.p.m. for 4 minutes. False positive results could appear if the test is read later than 4 minutes.

#### Semi-quantitative method

1. Make serial two-fold dilutions of the sample in 9 g/L saline solution.
2. Proceed for each dilution as in the qualitative method.

### READING AND INTERPRETATION

Examine the presence or absence of agglutination immediately after rotation using the light microscope (10x objective lens).

#### Interpretation

Agglutination	Reading	Report
Medium or large clumps	R	Reactive
Small clumps	W	Weakly reactive
No clumping or very slight "roughness"	N	Non-Reactive

In the semi-quantitative method, the titer is defined as the highest dilution showing a positive result.

### QUALITY CONTROL

Positive and Negative controls are recommended to monitor the performance of the procedure, as well as a comparative pattern for a better result interpretation.

All result different from the negative control result, will be considered as a positive.

### PERFORMANCE CHARACTERISTICS

1. **Analytical sensitivity:** Accurate titer determination of the reference material, under the described assay conditions (see, Calibration).
2. **Prozone effect:** No prozone effect was detected up to titers  $>1/64$ .
3. **Diagnostic sensitivity:** 100%
4. **Diagnostic specificity:** 100%

### INTERFERENCES

Bilirubin (20 mg/dL), hemoglobin (2 g/dL) and lipids (1000 mg/dL), do not interfere. Rheumatoid factor (300 IU/mL) interferes. Other substances may interfere<sup>4</sup>.

### LIMITATIONS OF THE PROCEDURE

- VDRL test is non-specific for syphilis. All Reactive samples should be retested with treponemic methods such as TPHA and FTA-Abs to confirm the results.
- A Non-Reactive result by itself does not exclude a diagnosis of syphilis.
- False positive results have been reported in diseases such as infectious mononucleosis, viral pneumonia, toxoplasmosis, pregnancy and autoimmune diseases.

### BIBLIOGRAPHY

1. George P. Schimid. Current Opinion in Infectious Diseases 1994; 7: 34-40.
2. Sandra A Larsen et al. Clinical Microbiology Reviews 1995; 8 (1): 1-21.
3. Sandra Larsen et al. A manual of Test for Syphilis American Public Health Association 1990: 1-192.
4. Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACC Press, 1995.

### PACKAGING

	5 mL VDRL stabilized antigen
Ref.: MO-165023 250 tests	1 mL Control +
	1 mL Control -

### SYMBOLS FOR IVD COMPONENTS AND REAGENTS

	Manufacturer		For <i>in vitro</i> diagnostic use only
	Don't re-use		Consult instructions for use
	Contains sufficient for <n> tests		Keep dry
	Catalogue Code		Temperature limitation
	Lot Number		Use by