One-Step Cocaine MonlabTest®

Only for professional *in vitro* diagnostic use. Store at 2 - 30°C.

INTENDED USE

The One-Step Cocaine MonlabTest® is a lateral flow, one-step immunoassay for the qualitative detection of benzoylcegonine (cocaine metabolite) in human urine at a cut-off of 300ng/mL. This product is used to obtain a visual, qualitative result and is intended for professional use. The assay should not be used without proper supervision and is not intended for over the counter sale to lay persons.

This assay provides only a preliminary analytical test result. A more specific alternative chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) has been established as the preferred confirmatory method by the National Institute on Drug Abuse (NIDA). Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are indicated.

SUMMARY

Derived from leaves of coca plant, cocaine is a potent central nervous system stimulant and a local anesthetic. Among the psychological effects induced by using cocaine are euphoria, confidence and a sense of increased energy, accompanied by increased heart rate, dilation of the pupils, fever, tremors and sweating. Cocaine is excreted in urine primarily as benzoylcegonine in a short period of time. Benzoylcegonine has a biological half-life of 5 to 8 hours, which is much longer than that of cocaine (0.5 to 1.5 hour) and can be generally detected for 24 to 60 hours after cocaine use.

Urine based screening tests for drugs of abuse range from simple immunoassay tests to complex analytical procedures. The speed and sensitivity of immunoassays have made them the most widely accepted method for screening urine for drugs of abuse. The One-Step Cocaine MonlabTest® is based on the principle of the highly specific immunochemical reactions of antigens and antibodies which are used for the analysis of specific compounds in biological fluids. This test is a rapid, visual, competitive immunoassay that can be used for the qualitative detection of benzoylcegonine, the major metabolite of cocaine, in human urine.

PRINCIPLE

The One-Step Cocaine MonlabTest® is a one-step immunoassay in which a chemically labeled drug (drug conjugate) competes with the drug which may be present in urine for limited antibody binding sites. The test device contains a membrane strip which was pre-coated with drug conjugate on the test band. A colored anti-benzoylcegonine monoclonal antibody-colloidal gold conjugate pad is placed at the right end of the membrane. In the absence of drug in the urine, the solution of colored antibody-colloidal gold conjugate and urine moves upward, chromatographically by capillary action, across the membrane. This solution then migrates to the immobilized drug conjugate zone on the test band region. The colored antibody-colloidal gold conjugate then attaches to the drug conjugate to form a visible line as the antibody complexes with the drug conjugate. Therefore, the formation of a visible precipitant in the test zone occurs when the test urine is **negative** for the drug. When the drug is present in the urine, the drug/metabolite antigen competes with the drug conjugate on the test band region for the limited antibody sites on the antibody-colloidal gold conjugate. When a sufficient concentration of drug is present, it will fill the limited antibody concentration of drug is present. It will fill the limited antibody binding sites. This will prevent attachment of the colored

antibody colloidal gold conjugate to the drug conjugate zone on the test band region. Therefore, absence of the color band on the test region indicates a **positive** result.

A control band that has a different antigen/antibody reaction is added to the immunochromatographic membrane strip at the control region (**C**) to indicate that the test has performed properly. This control line should always appear, regardless of the presence of drug or metabolite. This means that **negative** urine will produce **two** colored bands, and **positive** urine will produce only **one** band. The presence of this colored band in the control region also serves as 1) verification that sufficient volume of sample has been added, and 2) that proper flow was obtained.

MATERIALS SUPPLIED

- 40 individually wrapped test which include one disposable pipette each. The test contains a membrane coated with drug conjugate and a colloidal gold conjugate pad coated with anti-benzoylcegonine monoclonal antibody.
- One instruction sheet.

MATERIAL REQUIRED BUT NOT PROVIDED

- Specimen collection container.
- Timer

STORAGE AND STABILITY

The test cassette should be stored or refrigerated at temperature 2-30°C in the sealed pouch for the duration of the shelf-life.

PRECAUTIONS

- FOR IN VITRO DIAGNOSTIC USE.
- For professional use only.
- Urine specimens may be potentially infectious. Proper handling and disposal methods should be established.
- Avoid cross-contamination of urine samples by using a new specimen collection container and specimen pipette for each urine sample.
- Do not use it after the expiration date.

SPECIMEN COLLECTION AND HANDLING


The One-Step Cocaine MonlabTest® is formulated for use with urine specimens. Fresh urine does not require any special handling or pretreatment. Urine samples should be collected such that testing can be performed as soon as possible after the specimen collection, preferably during the same day. The specimen may be refrigerated at 2-8°C for 2 days, or frozen at -20°C for a longer period of time. Specimens that have been refrigerated must be equilibrated to room temperature prior to testing. Specimens previously frozen must be thawed, equilibrated to room temperature, and mixed thoroughly prior to testing.

Note: Urine specimens and all materials coming in contact with them should be handled and disposed of as if capable of transmitting infection. Avoid contact with skin by wearing gloves and proper laboratory attire.

TEST PROCEDURE


- Review "Specimen collection" instructions. Test device, samples of patient, and controls should be brought to room temperature (10-30°C) prior to testing. Do not open pouches until ready to perform the assay.
- Remove the test device from its protective pouch (bring the device to room temperature before opening the pouch to avoid condensation of moisture on the membrane). Label the device with patient or control number.
- Hold the dropper vertically and transfer 3 full drops of urine (approx. 0.15mL) to the specimen well (S) of the test device, and then start the timer. Use a separate pipette and device for each sample or control.
- **Read result between 3 to 8 minutes after the addition of samples.** Do not read result after 8 minutes.

INTERPRETATION OF RESULTS


Negative

Two colored lines appear in the viewing window. The line in the test region (T) is the drug probe line; the line in the control region (C) is the control line, which is used to indicate proper performance of the device. **The color intensity of the test line may be weaker or stronger than that of the control line.**

Positive

Only **one** colored line appears in the control region (**C**). The **absence** of a test line indicates a positive result.

Invalid

The test is invalid if the control line is not visible at five minutes. The test failed, or the test procedure was not followed properly. Verify the test procedure and repeat the test with a new testing device.

Note: A very faint line in the test region indicates that the Morphine in the samples is near to the cut-off level for the test. These should be re-tested or confirmed with a more specific method before a positive determination is made.

LIMITATIONS OF PROCEDURE

- The assay is designed for use with human urine only.
- A positive result with any of the tests indicates the presence of a drug/metabolite only and does not indicate or measure intoxication.
- There is a possibility that technical or procedural errors, as well as other substances or factors not listed, may interfere with the test and cause false results. See **SPECIFICITY** for lists of substances that will produce positive results, or that do not interfere with test performance.
- If it is suspected that the samples have been mislabeled or tampered with, a new specimen should be collected, and the test should be repeated.

QUALITY CONTROL

Good laboratory practice recommends the use of control materials to ensure proper kit performance. Quality control specimens are available from commercial sources. Before using a new kit with patient specimens, positive and negative controls should be tested. When testing the positive and negative controls, use the same assay procedure as with a urine specimen.

PERFORMANCE CHARACTERISTICS

Accuracy

The accuracy of the One-Step Cocaine MonlabTest® was evaluated in comparison to a commercially available one-step immunoassay at a cut-off of 300ng/mL. One hundred (100) urine samples, collected from presumed non-user volunteers, were tested by both procedures with 100% agreement.

In a separate study, fifty-six (56) urine samples, obtained from a clinical laboratory where they were confirmed as positives by GC/MS, were tested by both One-Step Cocaine MonlabTest® and a commercially available one-step immunoassay. Of the forty-six (46) samples with benzoylecgonine concentration ≥ 350 ng/mL, all were found to be positives by both methods (100% agreement). Of the ten (10) samples with benzoylecgonine concentration from 180 to 260ng/mL, all were determined negatives by both methods (100% agreement).

Reproducibility

The reproducibility of the One-Step Cocaine MonlabTest® was evaluated at four different sites using blind controls. Of the sixty samples with benzoylecgonine concentration of 150ng/mL, all were determined negatives. Of the sixty samples with benzoylecgonine concentration of 600ng/mL, all were determined positives.

Precision

The precision of the One-Step Cocaine MonlabTest® was determined by conducting the test with blind spiked controls. The control at the 150ng/mL of benzoylecgonine should give a negative result and the control at the 450 ng/mL of benzoylecgonine should give a positive result.

Concentration (ng/mL)	Number Tested	Correct Result	% Correct Result
150	50	50	100
450	50	50	100

Specificity

The specificity for the One-Step Cocaine MonlabTest® was tested by adding various drugs, drug metabolites, and other compounds that are likely to be present in urine. All compounds were prepared in drug-free normal human urine.

The following compounds were tested for cross-reactivity in the One-Step Cocaine MonlabTest® and found to be positives if the levels were greater than 300ng/mL:

Benzoylecgonine

Cocaine

The following compounds were found not to cross-react tested at concentrations up to 100µg/mL.

Acetaminophen	Bilirubin
Amitriptyline	Caffeine
Albumin	Chloroquine
Aspirin	(+)-Chlorpheniramine
Acetone	(+/-)-Chlorpheniramine
Ampicillin	Creatine
Aspartame	Desoxyephedrine
Atropine	Dextromethorphan
Dexbrompheniramine	Dopamine
4-Dimethylaminoantipyrine	Furosemide
Ecgonine	Glucose
Ecgonine Methyl Ester	Guaicol Glyceryl Ether
(+/-)-Ephedrine	Hemoglobin
(-)-Ephedrine	Imipramine
(+)-Epinephrine	(+/-)-Isoproterenol
Erythromycin	(1R,2S)-(-)-N-Methyl-Ephedrine
Ethanol	(+)-Naproxen
Lidocaine	(+/-)-Norephedrine
Oxalic Acid	Pheniramine
Penicillin-G	I-Phenylephrine
Phenothiazine	Sodium Chloride
Phenylethylamine	Thioridazine
Ranitidine	Trimethobenzamide
Sulindac	Vitamin C
Tyramine	

REFERENCES

1. Baselt, R, C. Disposition of Toxic Drugs and Chemicals in Man. Biomedical Publications. Davis. CA 1982.
2. Urine Testing for Drugs of Abuse. National Institute on Drug Abuse (NIDA). Research Monograph 73, 1986.
3. Ellenhorn, M. J. And Barce1oux, D. G. Medical Toxicology. Elsevier Science Publishing Company. 1nc New York. 1988
4. Fed. Register. Department of Health and Human Services. Mandatory Guidelines for Federal Workplace Drug Testing Programs, 53, 69, 11970-11979, 1988
5. Gilman, A, G, and Goodman, L, S, The Pharmacological Basis of Therapeutics. eds. MacMillan Publishing. New York NY, 1980.
6. Gorodetzky, C, W, Detection of Drugs of Abuse in Biological Fluids. In Martin WR(ed): Drug Addiction 1, New York, Spring-Verlag, 1977.
7. Harvey, R, A, Champe, P, C, Lippincotts illustrated Reviews. Pharmacology, 91-95,1992.
8. Hofmann F, E, A, Handbook on Drug and Alcohol Abuse: The Biomedical Aspects. New York, Oxford University Press, 1983.
9. McBay, A, J, Clin. Chem, 22, 33B-40B,1987

PACKAGING

CHEM-7006C-40

40 Cocaine Cassettes

GRAPHICAL SYMBOLS USED

Manufacturer

For *in vitro* diagnostic use only

Don't re-use

Consult instructions for use

Contains sufficient for <n> tests

Keep dry

Catalogue Code

Temperature limitation

Lot Number

Use by

Manufacturer: Shanghai Chemtron Biotech Co., Ltd.

Address: No. 518, Qingdai Rd., International

Medical Park, Pudong, 201318

Shanghai, PEOPLE'S REPUBLIC OF CHINA

EC REP European Representative:

Shanghai International Holding Corp. GmbH (Europe)

Address: Eiffestrasse 80, 20537 Hamburg Germany

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
REFERENCIAS

- Baselt. R, C. Disposition of Toxic Drugs and Chemicals in Man. Biomedical Publications. Davis. CA 1982.
- Urine Testing for Drugs of Abuse. National Institute on Drug Abuse (NIDA). Research Monograph 73, 1986.
- Ellenhorn. M. J. And BarceLoux, D. G. Medical Toxicology. Elsevier Science Publishing Company. Inc New York. 1988
- Fed. Register. Department of Health and Human Services. Mandatory Guidelines for Federal Workplace Drug Testing Programs, 53, 69, 11970-11979, 1988
- Gilman, A, G, and Goodman, L, S, The Pharmacological Basis of Therapeutics. eds. MacMillan Publishing. New York NY, 1980.
- Gorodetzky. C, W, Detection of Drugs of Abuse in Biological Fluids. In Martin WR(ed): Drug Addiction 1, New York, Spring-Verlag, 1977.
- Harvey, R, A, Champe, P, C, Lippincotts illustrated Reviews. Pharmacology, 91-95,1992.
- Hofmann F. E, A, Handbook on Drug and Alcohol Abuse: The Biomedical Aspects. New York, Oxford University Press, 1983.
- McBay, A, J, Clin. Chem, 22, 33B-40B,1987

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**Fabricante:** Shanghai Chemtron Biotech Co., Ltd.