CRYPTO-GIARDIA-ENTAMOEBIA MonlabTest® CE

MO-076005 20 TESTS IVD

Immunochromatographic test for the detection of and *Cryptosporidium parvum*, *Giardia lamblia* and *Entamoeba histolytica* antigens in human stool samples

For *in vitro* use only. Store at 2 - 30°C.

INTRODUCTION

Intended use

A single-step qualitative immunochromatographic test with separate lines to detect and help diagnose infection caused by the parasites, *Cryptosporidium spp.*, *Giardia lamblia* and *Entamoeba histolytica* based on their antigens in human stools.

The Crypto-Giardia-Entamoeba MonlabTest results should be evaluated by a doctor in combination with other clinical information about the patient. They should not be used as the only criterion for diagnosing the infection.

This test is exclusively for *in vitro* diagnostic use.

This test is for professional use. This is not a diagnostic test for use at the patient care location. This is not a self-diagnostic test.

General information

The Crypto-Giardia-Entamoeba MonlabTest qualitative immunochromatographic assay is a procedure for detecting *Cryptosporidium spp.*, *Giardia lamblia* and *Entamoeba histolytica* in human stools on separate lines. A positive signal on any of the positive test lines gives a good indication that an infection caused by *Cryptosporidium spp.*, *Giardia lamblia* and/or *Entamoeba histolytica* may be present.

Intended test population

The Crypto-Giardia-Entamoeba MonlabTest is intended for use on the general population. However, in the case of *Cryptosporidium* and *Giardia*, a significant number of samples can come from children, given that *Cryptosporidium* is considered the second main cause of gastroenteritis among children, behind Rotavirus, and *Giardia lamblia* is highly prevalent among children.

Disease incidence among the population

Cryptosporidiosis is an endemic disease worldwide, occurring much more frequently in developing countries. Its prevalence among patients with gastrointestinal symptoms is 1-4% in Europe and North America, and 3-20% in Africa, Asia, Australia, Central America, and South America. According to the report published by the ECDC (European Centre for Disease Prevention and Control) based on 2017 data, the reporting rate was 3.2 confirmed cases per 100,000 persons in the European Union and countries of the European Economic Area. The reporting rate was higher for children aged 0 to 4 years old with 12.5 cases per 100,000 persons.

When it comes to *Giardia lamblia*, giardiasis is a highly prevalent parasitism occurring worldwide, particularly among children. Worldwide, the estimated frequency is 200,000,000 infected people, 500,000 of which become ill. It is a cause of diarrhoea for up to 20% of cases in developing countries, but only for 3-7% in developed countries. The report published by the ECDC based on 2017 data indicates that there were 19,437 confirmed giardiasis cases and the reporting rate was 5.5 confirmed cases per 100,000 persons. The reporting rate was higher for children aged 0 to 4 years old, with 17.6 cases per 100,000 persons for boys and 14.9 cases per 100,000 persons for girls. The reporting rate reduced with age and was lowest among people ≥ 65 years old.

When it comes to *Entamoeba histolytica*, approximately 50 million cases of invasive disease from *E. histolytica* occur each year worldwide, resulting in up to 100,000 deaths. This is only the tip of the iceberg as only 10%-20% of infected people are symptomatic. The incidence of amebiasis is greater in developing countries due to deficient sanitation and socio-economic conditions, but with the increase in travel and emigration to developed countries, the infection is becoming more common in non-endemic regions.

The overall prevalence of amebiasis in the United States of America is approximately 4%. Infection by the non-pathogenic and morphologically identical species, *Entamoeba dispar*, is 10 times more common than infection by *E. histolytica*. Furthermore, only 10% of infections by *E. histolytica* cause invasive disease. Therefore, only 1% of people with positive results for Entamoeba via microscopy develop symptomatic amebiasis.

Characteristics of *Cryptosporidium spp.* and its infection

Cryptosporidium spp. is a gastrointestinal protozoan parasite that infects humans and animals, causing a disease known as intestinal, tracheal and pulmonary cryptosporidiosis. *Cryptosporidium parvum* belongs to the *Cryptosporidium* genus, and along with *Cryptosporidium hominis*, causes over 90% of human cryptosporidiosis infections. Its most common route of transmission is the ingestion of cysts present in water or food contaminated with human faecal matter.

Cryptosporidium has a very short life cycle and each generation can develop and mature in 12-24 hours. Along with its self-infecting nature, this means the intestinal tract becomes infected by a large number of parasites in just a few days. It is usually a short, acute infection, although it can also be asymptomatic. Diarrhoea, weight loss, and abdominal cramps are the clinical symptoms of the disease, which can be more serious among children or immunosuppressed patients, as the diarrhoea episodes are worse and can lead to electrolyte imbalance.

Pulmonary and tracheal cryptosporidiosis is associated with coughing and fever, often accompanied by intestinal changes. These symptoms generally develop 4-6 days after infection, although they can appear at any moment within 2-10 days.

Characteristics of *Giardia lamblia* and its infection

Giardia lamblia is the gastrointestinal parasite that causes giardiasis in humans, with a clinical condition that varies from asymptomatic patient carriers to cases of acute or long-lasting diarrhoea, particularly among children.

The disease's most common route of transmission involves the ingestion of cysts present in water or food contaminated with faecal matter.

Symptoms appear 1-3 weeks after contagion, following an incubation period of approximately 8 days. The most frequent presentation is characterised by diarrhoea, weight loss, colic-type abdominal pain, and arrested growth and development. Symptom onset may be sudden or gradual; the disease can be self-limited or produce severe diarrhoea with poor intestinal absorption. Gastrointestinal function changes are common in children with persistent symptoms.

Characteristics of *Entamoeba histolytica* and its infection

Entamoeba histolytica is the agent that causes dysentery and amebiasis. This parasitic disease is the third cause of mortality among humans after malaria and schistosomiasis, and is estimated to be responsible for 50,000-100,000 deaths a year.

It is usually transmitted by the ingestion of cysts present in food or water contaminated with infected faecal material. The disease can manifest as an acute, chronic, or asymptomatic infection. The symptoms of amebiasis are both intermittent and mild among most infected human beings (different gastrointestinal disorders, including diarrhoea, colitis). In more severe cases, dysentery is the most common manifestation of invasive intestinal amebiasis. Extra-intestinal amebiasis occurs in a small percentage of infected adults, with liver abscesses the most common location.

THE BASIS OR BASIC PRINCIPLES OF THE TEST

The Crypto-Giardia-Entamoeba MonlabTest is an immunochromatographic test that uses a combination of:

- Blue latex particles conjugated to a specific antibody that targets *Cryptosporidium spp.* antigens and works in conjunction with a specific antibody for *Cryptosporidium spp.* located in the membrane (position 1 on the housing).
- Red latex particles conjugated to specific antibodies that target *Giardia lamblia* antigens and work in conjunction with specific antibodies for *Giardia lamblia* located in the membrane (position 2 on the housing).
- Green latex particles conjugated to specific antibodies that target *Entamoeba histolytica* antigens and work in conjunction with specific antibodies for *Entamoeba histolytica* located in the membrane (position 3 on the housing).
- Purple latex particles conjugated to an antigen recognised by a specific antibody for said antigen bound to the membrane, comprising the so-called test control line (position C on the housing).

This test first involves treating the sample with the sample dilution buffer (included in the kit) to extract the parasites. Once extracted, it simply requires adding a specific volume of the supernatant to the reagent strip and waiting for 10 minutes.

The coloured particles migrate when the extracted sample flows along the strip. In the event of a positive sample, the antibodies present in the particles will react with the parasite's antigens, forming a particle-antigen complex, which will also be retained in the membrane by other parasite-specific antibodies.

Different coloured lines will appear depending on the parasite content in the sample. These lines are used to interpret the result after 10 minutes of incubation at room temperature (see Fig. 1).

MATERIALS INCLUDED IN THE KIT

- 20 Reaction devices,
- 1 Sample dilution buffer (30 mL),
- Instructions for use.

MATERIALS NOT INCLUDED IN THE KIT

- Suitable devices for sampling stools,
- Vortex,
- Test tubes/microtubes,
- Graduated pipettes,
- Chronometer.

PRECAUTIONS

- The patient samples (stools) must be handled with care as they may contain infectious agents. Disposable gloves must be worn throughout the process.
- The sample dilution buffer contains 0.095% sodium azide as an anti-microbial agent. Avoid direct contact with skin and mucous membranes. Dispose of appropriately. Do not use the buffer if there are any signs of contamination or precipitation. The sample dilution buffer contains animal-sourced material (albumins). All animal-sourced material used are trade references and come with their relevant certificate stating they are not infectious.
- Do not eat, drink, smoke or prepare food in the area where reagents and samples are handled.
- On completing the task, clean all work surfaces with soap and water, and disinfect them with a suitable solution. Lastly, dispose of the protective gloves correctly and wash your hands with soap and water, rubbing them well.
- Do not exchange components from kits with different batch numbers.
- If the test is stored refrigerated, allow all the kit's components and the stool samples to reach room temperature, as cold reagents and/or samples can reduce the test's functionality. Allowing 20 to 30 minutes to reach room temperature is recommended.
- Only use the reagents *in vitro*.
- The immunochromatographic strips are single use. One strip is used for each sample requiring analysis.
- Do not use the kit components after their expiry dates.
- If the outer box is damaged, the product can still be used providing none of the components have been damaged.
- Throw the test away if the primary packaging is either broken or impaired.
- Remove the strip or device from the aluminium bag when the test is going to be performed to avoid unnecessary excessive exposure of the product to any potentially damaging environmental factors.
- Adding the correct volume of extracted sample to the reaction device is very important. If lower than indicated, the chromatography may not develop because insufficient sample reached the reaction zones; excessive volume could spill inside the device and the chromatography may not develop correctly.
- Please dispose of the used product in accordance with current legislation.
- Do not use the test if there are any coloured lines in the results zone prior to use.

- Sampling a suitable amount is very important: around **50 mg** if they are **solid or soft samples** (a small pellet around 3 mm in diameter) and **100 µl** if they are **liquid samples**; these amounts are extracted in 1 mL of the sample dilution buffer included with the kit. Too much sample in relation to the amount of added buffer prevents the chromatography from developing correctly. This is especially critical with solid or soft samples, as sampling the correct amount is not that simple.
Too less sample in relation to the amount of added buffer can affect the test result. This is especially critical for samples with a low concentration of analyte.
- Do not discard the kit's outer box until all the content has been used. The box contains essential information about the product's CE marking and the batch number.
- Any serious incident involving the product must be reported to the manufacturer and the relevant authority in the member state where the user and/or patient are located.

STORAGE

The Crypto-Giardia-Entamoeba MonlabTest product can be stored at any temperature from +2 to +30 °C.
The use-by date of each component is printed on the aluminium wrappers.
The reagent strips must be used once the protective packaging has been opened after reaching room temperature (if stored refrigerated).
The product maintains the claimed stability once open.

SAMPLES

- This test is designed to analyse human stool samples. Collecting the stool sample as soon as symptoms appear is recommended.
- Do not use samples that have been collected in transport media, or those with added preservative agents (such as formalin, SAF, PVA or similar) or enrichment media, as their presence could interfere with correct performance of the test.
- Analysing **fresh, untreated samples** is recommended. If they need to be kept for a certain amount of time, they can be stored in a refrigerator (+2 to 8°C) for 1 or 2 days. For longer time periods, they should be frozen to -20°C, bearing in mind that some samples may lose their immunoreactivity after freezing.
- Ensure any frozen samples have fully thawed at room temperature before analysing them.
- Avoid freeze-thaw cycles, as this could alter the immunological recognition of the parasites.

PREPARING STOOL SAMPLES

General note: use disposable gloves throughout the test as infectious samples are being handled. Once the work is concluded, do not forget to comply with the hygiene procedures detailed in point 4 of the "Precautions" section.

- Homogenise the sample beforehand to make it as representative as possible.
- For **liquid or semi-liquid stools**, use a pipette to add **100 µl** of sample to a labelled 1.5 mL microtube. If the **sample is solid**, sample approximately **50 mg** (a small pellet around 3 mm in diameter) and add it to a labelled 1.5 mL microtube.
- Add **1 mL of the sample diluent** to the aforementioned 1.5 mL microtube (or the volume required to ensure a proportion of 100 µl (~50 mg) of sample in 1 mL of dilution buffer).
- Agitate in the vortex for 30 seconds or the time required to ensure the sample suspends completely in the buffer. If a vortex is not available, agitate vigorously by hand to ensure the stool sample suspends in the buffer.
- Centrifuge the 1.5 mL microtubes for 5 minutes at 700 xg (approximately 3,000 rpm) in order to sediment the solid particles, using a small centrifuge adapted to these microtubes. If a suitable centrifuge is unavailable, please wait 3-5 minutes until the solid particles have settled to the bottom of the tube. However, the test works best when using a clear sample solution extracted using a centrifuge.
- Once extracted, stools should be analysed immediately or within 6 hours of extraction using the Crypto-Giardia-Entamoeba MonlabTest (see the "Procedure" section).

PROCEDURE

Once the samples have been prepared as indicated above, proceed as follows:

- Remove the reaction device from the aluminium bag. Discard the small desiccant sachet as its sole purpose is to keep the test dry and it is not used to perform the test.
- Use a graduated pipette to sample **125 µl of the supernatant** of the sample and add it to the sample area of the reaction device (the circular window indicated by an arrow). Do not add solid particles with the liquid.
- Wait **10 minutes**, then read and interpret the result.

READING THE RESULTS

Five images in Fig. 1 show some of the possible results with the Crypto-Giardia-Entamoeba MonlabTest, with combinations possible should the patient present co-infections with the detected parasites.

Four different coloured lines can be seen:

- An upper purple line: this is the control line that indicates that the test functions correctly.
- A middle green line: this indicates that *Entamoeba histolytica* is present in the sample.
- A middle red line: this indicates *Giardia lamblia* is present in the sample.
- A lower blue line: this indicates that *Cryptosporidium spp.* is present in the sample.

The purple control line must always appear. If any of the other three coloured lines also appear, this indicates that *Entamoeba histolytica*, *Giardia lamblia* and/or *Cryptosporidium spp.* is/are present in the sample.


Fig. 1: Examples of possible results using the Crypto-Giardia-Entamoeba (C-G-E) MonlabTest.

- Strip 1: NEGATIVE result.**

Only a horizontal **PURPLE** line appears on the reaction device (aligned with the letter "C" marked on the housing). This is the control line, which must always appear, as it indicates that the chromatography worked correctly.

- Strips 2-5: POSITIVE results:**

- Strip 2: *Cryptosporidium positivo***

A **PURPLE** (control) line appears aligned with the letter "C" and a **BLUE** line appears below the control line (aligned with the number "1" marked on the housing).

- Strip 3: *Giardia positivo***

A **PURPLE** (control) line appears aligned with the letter "C" and a **RED** line appears below the control line (aligned with the number "2" marked on the housing).

- Strip 4: *Entamoeba positivo***

A **PURPLE** (control) line appears aligned with the letter "C" and a **GREEN** line appears below the control line (aligned with the number "3" marked on the housing).

- Strip 5: *Cryptosporidium positivo + Giardia positivo + Entamoeba positivo***

A **PURPLE** (control) line appears aligned with the letter "C". Directly below, a **GREEN** line appears (aligned with the number "3" marked on the housing), below which a **RED** line appears (aligned with the number "2" marked on the housing), below which a **BLUE** line appears (aligned with the number "1" marked on the housing).

- Strip 6: INVALID results**

The **PURPLE** control line does not appear, the control line colour is clearly altered, or different colours appear on the positive lines (compared to those previously indicated). This is an indication that the test has not worked correctly. This may be due to any of the following reasons:

- some of the reagents have deteriorated or the test has expired.
- the sample was not prepared in accordance with the instructions for use.
- the sample has a high blood content.

Repeating the test using a new strip in the event of an invalid result is recommended, strictly adhering to the instructions for use described in this manual. In the case of samples with blood, using an alternative technique is advised, as the instability problem does not usually depend on the strip used, but on the matrix of the sample itself.

Any line that appears after 10 minutes will have no diagnostic value due to the nature of the sample.

NOTE: a doctor should make the final and definitive diagnosis. This test only detects *Cryptosporidium spp.*, *Giardia lamblia*, and *Entamoeba histolytica* in a sample. It does not constitute an argument to state that a person is infected by any of these parasites.

QUALITY CONTROL

The test is invalid if a **PURPLE** line does not appear, either because the test was not performed properly or because the reagents have deteriorated. In this case, repeat the test, strictly following the operating protocol detailed in these instructions for use.

WARNING: including our controls with a known result is recommended to ensure that correct data is obtained.

LIMITATIONS OF THE PROCEDURE

- The Crypto-Giardia-Entamoeba MonlabTest is for differentiating and identifying *Cryptosporidium spp.*, *Giardia lamblia*, and *Entamoeba histolytica*, detecting their presence in human stools as long as the parasites' load is equal to or above the assay's limit of detection.
Detecting other complications such as respiratory cryptosporidiosis described among immunosuppressed patients, most of whom have advanced AIDS/HIV disease, or extraintestinal amebiasis in the case of *Entamoeba histolytica*, is not the intended use of the test.
- The test is qualitative, not quantitative, although the intensity of the positive lines does relate to the amount of parasites detectable in the stool sample.
- The test has shown good correlation with other techniques such as microscopy and PCR (in the case of *E. histolytica* detection) after analysing a large amount of stool samples.

- Using the correct amount of sample is important. A defective sample may lead to very weak positive results. In this case, the test should be repeated with more sample. On the other hand, excess sample could cause the test to develop very slowly and even prevent it from developing correctly (control lines not visible). In this case, the test should be repeated with a smaller amount of sample (see the "Preparing stool samples" section).
- A negative result does not rule out the possibility of infection by *Cryptosporidium spp.*, *Giardia lamblia*, and *Entamoeba histolytica*. Not detecting these parasites may result from factors such as sampling during an unsuitable phase of the disease (parasite excretion in stools is usually intermittent); incorrect storage of the sample; unsuitable transportation of the sample, and so on.
- A positive result does not rule out the presence of other pathogenic agents, which may cause a co-infection of any of these three parasites with other micro-organisms. In any case, co-infections can only be clarified by differential diagnosis.
- The test results must be interpreted alongside the available data from epidemiological studies, the patient's clinical evaluation, and other diagnostic procedures.
- Analysing some samples can produce undefined coloured lines, caused by negative samples in most cases. The test should be repeated if undefined coloured lines appear. If the same result is obtained again, performing the analysis using another analytical method is suggested.
- Given the similarity between *Entamoeba histolytica* and *Entamoeba dispar*, some cross-reactivity has been found to occur with the non-pathogenic species, *E. dispar* when it is present in the sample at concentrations equal to or higher than 1.56×10^5 cysts/mL.
- Test results should be interpreted in conjunction with information available from epidemiological studies, patient clinical evaluation and other diagnostic procedures.
- Test results to verify the eradication efficacy (test of cure) have not been established.

DETECTION CAPABILITY

When using oocysts of *Cryptosporidium parvum* inactivated by gamma irradiation (batch B971-627) from Biopoint (Biomérieux) as a reference standard, the detection capability of the Crypto-Giardia-Entamoeba MonlabTest for *Cryptosporidium* in faecal matter is established at 2.25×10^5 oocysts/mL. Using the supernatant collected after the first washing of faecal matter from mice experimentally infected with *Cryptosporidium parvum* from R-Biopharm as a reference, the detection capability of the Crypto-Giardia-Entamoeba MonlabTest in the sample dilution buffer for *Cryptosporidium* is established at 5.63×10^4 cells/mL. Using cultures of *Giardia lamblia* (ATCC 30957) in their two possible biological forms (cysts and trophozoites), the detection capability of the Crypto-Giardia-Entamoeba MonlabTest in a faecal matter for *Giardia* is established as 3120 cysts/mL and 6.09×10^4 trophozoites/mL. Using the supernatant of a quantified culture of *Entamoeba histolytica* HM-1:IMSS from R-Biopharm, the detection capability of the Crypto-Giardia-Entamoeba MonlabTest for *E. histolytica* in faecal matter is established as 913 trophozoites/mL.

ANALYTICAL SPECIFICITY

The monoclonal antibodies used to detect *Cryptosporidium spp.*, *Giardia lamblia*, and *Entamoeba histolytica* specifically recognise antigens from these parasites. However, antibodies that target *E. histolytica* have been detected to present some cross-reactivity with the species, *Entamoeba dispar*. No cross-reactions with other antigens and/or micro-organisms present in stool samples have been detected. See the cross-reactivity section.

DIAGNOSTIC SENSITIVITY AND SPECIFICITY

Internal evaluation

The study included the analysis of 114 samples for *Cryptosporidium spp.* and 113 samples for *Giardia lamblia* characterised using microscopy as a reference technique. For *Entamoeba histolytica*, a total of 104 samples characterised using PCR as a reference technique were analysed, given that when using microscopy, *E. histolytica* is morphologically indistinguishable from the non-pathogenic species, *E. dispar* and *E. moshkovskii*.

Microscopy (*Cryptosporidium y Giardia*)

<i>Cryptosporidium spp.</i>	Microscopy	
	Positive	Negative
Crypto-Giardia-Entamoeba MonlabTest	26	0
	0	88

Sensitivity Concordance = $26/(26+0) = >99,9\%$

Specificity Concordance = $88/(88+0) = >99,9\%$

<i>Giardia lamblia</i>	Microscopy	
	Positive	Negative
Crypto-Giardia-Entamoeba MonlabTest	54	0
	1	58

Sensitivity Concordance = $54/(54+1) = 98,2\%$

Specificity Concordance = $58/(58+0) = >99,9\%$

The false negative divergent sample was analysed using the R-Biopharm ELISA RIDASCREEN® Giardia, which also provided a negative result, suggesting that it may be a positive sample with a very low concentration in the parasite and it is below the limit of sensitivity for the evaluated immunoassays.

PCR (*Entamoeba histolytica*)

<i>Entamoeba histolytica</i>	PCR	
	Positive	Negative
Crypto-Giardia-Entamoeba MonlabTest	9	4
	0	91

Sensitivity Concordance = $9/(9+0) = >99,9\%$

Specificity Concordance = $91/(91+4) = 95,8\%$

The 4 false positive divergent samples were analysed using PCR on *E. dispar*, giving a negative result, ruling out that the false positive was caused by a cross reaction with this species.

REPEATABILITY

INTRA-ASSAY PRECISION:

Five replicates of a series of 1/2 dilutions of internal standards (sensitivity curve), along with real samples established as NC ("negative sample"), LPC ("weak positive sample"), and PC ("positive sample") for each analyte detected by the test, were measured on the same day by the same person.

No differences between the replicates of each analyte were detected, or with either the standard or the samples, so the repetitiveness was 100%.

REPRODUCIBILITY

INTER-DAY PRECISION:

The same batch of Crypto-Giardia-Entamoeba MonlabTest was used to measure a sensitivity curve and real samples established as NC ("negative sample"), LPC ("weak positive sample"), and PC ("positive sample") for each analyte in duplicate over a five day period spaced out over time. Reproducibility with all analytes was 100% using the standard. For some LPC samples, differences around 1/2 half dilution were observed for Crypto and Entamoeba analytes, while the reproducibility for Giardia was 100%. The reproducibility with the three analytes was 100% for the other samples.

INTER-OPERATOR PRECISION:

Three persons measured a sensitivity curve for each analyte and real samples established as NC ("negative sample"), LPC ("weak positive sample") and PC ("positive sample") for each analyte in triplicate. The observed differences were no superior to a 1/2 dilution with the standard in any of the cases, and the differences were equal to or less than 1/2 half dilution with the real samples.

INTER-BATCH PRECISION:

Three different batches of the Crypto-Giardia-Entamoeba MonlabTest were evaluated in parallel by analysing the usual sensitivity curves for this test and real samples established as NC ("negative control"), LPC ("low positive control") y PC ("positive control") for each analyte. The analysis was performed by one person on the same day. The biggest observed differences were equal to less than 1/2 half dilution to the standard, and the same results were obtained using real samples, providing evidence for high inter-batch precision for the test.

THE HOOK OR PROZONE EFFECT

There is evidence that the Crypto-Giardia-Entamoeba MonlabTest:

- Does not present a Hook effect for Crypto until *Cryptosporidium parvum* concentrations of 7.2×10^6 oocysts/mL in faecal matter when using oocysts of *Cryptosporidium parvum* inactivated by gamma irradiation (batch B971-627) from Biopoint (Biomérieux) as a reference standard.
- A Hook effect was not detected for Giardia until *Giardia lamblia* concentrations of 1.6×10^6 cysts/mL in faecal matter using a quantified culture grown from the strain ATCC 30957 as standard.
- Does not present a Hook effect for Entamoeba until *Entamoeba histolytica* concentrations of 2.3×10^6 trophozoites/mL in faecal matter, when using the supernatant of a quantified culture of *Entamoeba histolytica* HM-1:IMSS from R-Biopharm as standard, albeit that very high concentrations of 9.3×10^6 trophozoites/mL slightly lower the signal intensity.

INTERFERING SUBSTANCES

At the concentrations indicated in the table, the following substances did not affect the test results when added to extracted stool samples (negative and weak positive in the analyte):

Interfering substances	Concentration	Interfering substances	Concentration
Loperamide (Fortasec®)	6.25 µg/mL	Omeprazole	7.8 µg/mL
Racecadotril (Tiorfan®)	0.12 mg/mL	Lansoprazole	23.5 µg/mL
Gelatin tannate (Tanagel)	1.17 mg/mL	Ibuprofen	0.94 mg/mL
Diosmectite (Dialiv)	7 mg/mL	Acetylsalicylic acid	1.6 mg/mL
Ultra levura (<i>Saccharomyces boulardii</i> CNM I-745)	0.2 mg/mL	Paracetamol	1.6 mg/mL
Infloran (<i>Lactobacillus acidophilus</i> and <i>Lactobacillus bifidus</i>)	1.17×10^6 cell/mL	Palmitic/ stearic acid	40 % (w/w)
Metronidazole	0.78 mg/mL	Mucin	3.5 % (w/w)
Vancomycin	2.34 mg/mL	Barium sulphate	2.9 mg/mL
Bismuth subsalicylate (Peptobismol®)	5 % (v/v)	Blood	40 % (v/v)
Almagate (Almax®)	3.1 mg/mL		

CROSS-REACTIVITY

The Crypto-Giardia-Entamoeba MonlabTest was tested against different organisms with no cross-reactivity occurring:

Staphylococcus aureus, *Staphylococcus epidermis*, *Staphylococcus haemolyticus*, *Staphylococcus saprophyticus*, *E.coli* 0157:H7, *E.coli* ST131-025b, *E.coli* ST131, *E.coli* O6, *Salmonella typhimorium*, *Yersinia enterocolitica*, *Campylobacter coli*, *Campylobacter jejuni*, *Listeria monocytogenes*, *Shigella flexneri* and *Salmonella enteritidis* present at a concentration of 1×10^7 cells/mL; *Entamoeba dispar* and *H. pylori* CPY6271 of unknown concentration and sourced from BEI Resources diluted ten times in faecal matter; and the viruses Rotavirus A and Adenovirus 51 at $TCID_{50} \geq 5 \times 10^2$ /mL, and Astrovirus Type 1 at $TCID_{50} = 2.8 \times 10^6$ /mL and Norovirus 1 (clon W3) a $TCID_{50} = 3,3 \times 10^6$ /mL.

Cross-reactivity with *Entamoeba dispar* at concentrations equal to or above 1.56×10^5 cysts/mL was observed on the *Entamoeba histolytica* test line.

MODIFICATIONS COMPARED TO THE PREVIOUS VERSION OF THIS DOCUMENT

The structure of the general information provided at the beginning of the document is improved, breaking down the different sections in order to be coherent with the IVDR requirements.



The components of the kit are better identified.

Regarding the test performance parameters, modifications have been introduced in the following sections: "Precautions", "Storage", "Samples", "Preparing stool samples", "Reading the results", "Limitations of the procedure", "Detection capability", "Analytical specificity", "Diagnostic sensitivity and specificity", "Repeatability", "Reproducibility", "Hook or prozone effect", formerly "Hook effect", "Interfering substances" and "Cross-reactivity".

REFERENCES

- Goñi, P. *et al.*, Evaluation of an immunochromatographic dip strip test for simultaneous detection of *Cryptosporidium* spp, *Giardia duodenalis*, and *Entamoeba histolytica* antigens in human faecal samples. *Eur J Clin Microbiol Infect Dis.* 31(8):2077-82 (2012).
- Roka, M. *et al.*, Prevalence of intestinal parasites in HIV-positive patients on the island of Bioko, Equatorial Guinea: Its relation to sanitary conditions and socioeconomic factors. *Sci Total Environ.* 432:404-11 (2012).
- Weitzel, T. *et al.*, Evaluation of seven commercial antigen detection tests for *Giardia* and *Cryptosporidium* in stool samples. *Clin Microbiol Infect.* 12:656-659 (2006).
- Lake, I.R. *et al.*, Case-control study of environmental and social factors influencing cryptosporidiosis. *Eur J Epidemiol.* 22:805-811 (2007).
- Khalakina, A. *et al.*, Is drinking water a risk factor for endemic cryptosporidiosis? A case-control study in the immunocompetent general population of the San Francisco Bay Area. *BMC Public Health.* 3: 1471-2458 (2003).
- Cacciò, S.M. *et al.*, Unravelling *Cryptosporidium* and *Giardia* epidemiology. *Trends Parasitol.* 21: 430-437 (2005).
- Chalmers RM, *et al.*, Looking for *Cryptosporidium*: the application of advances in detection and diagnosis. *Trends Parasitol;* 29(5):237-251(2013).
- Nygard, K. *et al.*, A large community outbreak of waterborne giardiasis-delayed detection in a non- endemic urban area. *BMC Public Health.* 6: 141 (2006).
- Monis, P.T., *et al.*, *Cryptosporidium* and *Giardia*- zoonoses: fact or fiction? *Infect Genet. Evol.* 3: 233- 244 (2003).
- Ekdahl, K., *et al.*, Imported giardiasis: impact of international travel, immigration, and adoption. *Am. J. Trop. Med.* 72: 825-830 (2005).
- Hooshyar H., *et al.*, *Giardia lamblia* infection: review of current diagnostic strategies, *Gastroenterol Hepatol Bed Bench;*12(1):3-12 (2019).
- Rosoff JD, *et al.*, Isolation and identification of a *Giardia lamblia*-specific stool antigen (GSA 65) useful in coprodiagnosis of giardiasis. *J Clin Microbiol;* 23(5):905-910 (1986).
- Van den Bossche D, *et al.*, Comparison of four rapid diagnostic tests, ELISA, microscopy and PCR for the detection of *Giardia lamblia*, *Cryptosporidium* spp, and *Entamoeba histolytica* in feces. *J Microbiol Methods;*110:78-84 (2015).
- Kantor M, *et al.*, *Entamoeba histolytica*: Updates in Clinical Manifestation, Pathogenesis, and Vaccine Development. *Can J Gastroenterol Hepatol.* Dec 2 (2018).
- Dolabella SS, *et al.*, Amoebic liver abscess production by *Entamoeba dispar*. *Ann Hepatol;* 11(1):107-117 (2012).
- Ungar, B. L. P. *et al.*, Use of a monoclonal antibody in an enzyme immunoassay for the detection of *Entamoeba histolytica* in fecal specimens. *Am. J. Trop. Med. Hyg.* 34: 465-472 (1985).
- Strachan, W. D. *et al.*, Immunological differentiation of pathogenic and non-pathogenic isolates of *Entamoeba histolytica*. *The Lancet* 12: 561 - 563 (1988).
- Yau, Y. C. W. *et al.*, Development of monoclonal antibodies which specifically recognise *Entamoeba histolytica* in preserved stool samples. 39: 716-719 (2001).
- Braga, L.L. Seropositivity for intestinal colonization with *Entamoeba histolytica* and *Entamoeba dispar* in individuals in Northeastern Brazil. 36: 3044-3045 (1998).
- Petri Jr. W.A., *et al.*, Pathogenic and nonpathogenic strains of *Entamoeba histolytica* can be differentiated by monoclonal antibodies to the galactose-specific adherence lectin. 58: 1802-1806 (1990).

SYMBOLS FOR IVD COMPONENTS AND REAGENTS

	Manufactured by		For <i>in vitro</i> diagnostic use
	Contains sufficient for <n> tests		Please read pack insert
	Catalogue number		Store at
	Lot number		Expiry date
	This product fulfils the requirements of Directive 98/79/EC on <i>in vitro</i> diagnostic medical devices		

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- Goñi, P. *et al.*, Evaluation of an immunochromatographic dip strip test for simultaneous detection of *Cryptosporidium spp*, *Giardia duodenalis*, and *Entamoeba histolytica* antigens in human faecal samples. *Eur J Clin Microbiol Infect Dis.* 31(8):2077- 82 (2012).
- Roka, M. *et al.*, Prevalence of intestinal parasites in HIV-positive patients on the island of Bioko, Equatorial Guinea: Its relation to sanitary conditions and socioeconomic factors. *Sci Total Environ.* 432: 404-11 (2012).
- Weitzel, T. *et al.*, Evaluation of seven commercial antigen detection tests for *Giardia* and *Cryptosporidium* in stool samples. *Clin Microbiol Infect.* 12:656-659 (2006).
- Lake, I.R. *et al.*, Case-control study of environmental and social factors influencing cryptosporidiosis. *Eur J Epidemiol.* 22:805-811 (2007).
- Khalakdina, A. *et al.*, Is drinking water a risk factor for endemic cryptosporidiosis? A case-control study in the immunocompetent general population of the San Francisco Bay Area. *BMC Public Health.* 3: 1471-2458 (2003).
- Cacciò, S.M. *et al.*, Unravelling *Cryptosporidium* and *Giardia* epidemiology. *Trends Parasitol.* 21: 430- 437 (2005).
- Chalmers RM, *et al.*, Looking for *Cryptosporidium*: the application of advances in detection and diagnosis. *Trends Parasitol;* 29(5):237-251(2013).
- Nygard, K. *et al.*, A large community outbreak of waterborne giardiasis-delayed detection in a non- endemic urban area. *BMC Public Health.* 6: 141 (2006).
- Monis, P.T., *et al.*, *Cryptosporidium* and *Giardia*- zoonoses: fact or fiction?. *Infect Genet. Evol.* 3: 233- 244 (2003).
- Ekdahl, K. *et al.*, Imported giardiasis: impact of international travel, immigration, and adoption. *Am. J. Trop. Med.* 72: 825-830 (2005).
- Hooshyar H., *et al.*, *Giardia lamblia* infection: review of current diagnostic strategies, *Gastroenterol Hepatol Bed Bench;*12(1):3-12 (2019).
- Rosoff JD, *et al.*, Isolation and identification of a *Giardia lamblia*-specific stool antigen (GSA 65) useful in coprodiagnosis of giardiasis. *J Clin Microbiol;* 23(5):905-910 (1986).
- Van den Bossche D, *et al.*, Comparison of four rapid diagnostic tests, ELISA, microscopy and PCR for the detection of *Giardia lamblia*, *Cryptosporidium spp*, and *Entamoeba histolytica* in feces. *J Microbiol Methods;*110:78-84 (2015).
- Kantor M, *et al.*, *Entamoeba histolytica*: Updates in Clinical Manifestation, Pathogenesis, and Vaccine Development. *Can J Gastroenterol Hepatol.* Dec 2 (2018).
- Dolabella SS, *et al.*, Amoebic liver abscess production by *Entamoeba dispar*. *Ann Hepatol;* 11(1):107-117 (2012).
- Ungar, B. L. P. *et al.*, Use of a monoclonal antibody in an enzyme immunoassay for the detection of *Entamoeba histolytica* in fecal specimens. *Am. J. Trop. Med. Hyg.* 34: 465-472 (1985).
- Strachan, W. D. *et al.*, Immunological differentiation of pathogenic and non-pathogenic isolates of *Entamoeba histolytica*. *The Lancet* 12: 561 - 563 (1988).
- Yau, Y. C. W. *et al.*, Development of monoclonal antibodies which specifically recognise *Entamoeba histolytica* in preserved stool samples. *39: 716-719 (2001).*
- Braga, L.L. Seropositivity for intestinal colonization with *Entamoeba histolytica* and *Entamoeba dispar* in individuals in Northeastern Brazil. *36: 3044-3045 (1998).*
- Petri Jr. W.A., *et al.*, Pathogenic and nonpathogenic strains of *Entamoeba histolytica* can be differentiated by monoclonal antibodies to the galactose-specific adherence lectin. *58: 1802-1806 (1990).*

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