## Legionella MonlabTest®

MO-804003 20 TESTS

One step test to detect *Legionella pneumophila* antigens

IVD

A rapid, one step test for the qualitative detection of *Legionella pneumophila* serogroup 1 antigens in urine specimens.  
For *in vitro* diagnostic use only.

### INTENDED USE

An *in vitro* rapid immunochromatographic assay for the qualitative detection of *Legionella pneumophila* serogroup 1 antigen (*L. pneumophila* serogroup 1 antigen) in urine specimens from patients with symptoms of pneumonia. It is intended to aid in the presumptive diagnosis of *Legionella* infection (Legionnaires' Disease) caused by *L. pneumophila* serogroup 1 in conjunction with culture and other methods.

### SYNTHESIS

Legionnaires' disease is a serious form of pneumonia that carries with it a mortality rate in the order of 10-15% in otherwise healthy individuals. Symptoms include a flu-like illness, followed by a dry cough and frequently progress to pneumonia. Approximately 33% of people infected may also present with diarrhoea and vomiting and around 50% may show signs of mental confusion. The incubation period normally ranges from 2-10 days with 3-6 days the typical illness onset time after exposure. Legionnaires' disease may present as an outbreak of two or more cases following a limited temporal and spatial exposure to a single source, as a series of independent cases in an area in which it is highly endemic or as sporadic cases without any obvious temporal or geographical grouping. Outbreaks have occurred repeatedly in buildings such as hotels and hospitals.

The *Legionella* MonlabTest® allows for early diagnosis of *Legionella pneumophila* serogroup 1 infection through detection of a specific soluble antigen present in the urine of patients with Legionnaires' Disease. *Legionella pneumophila* serogroup 1 antigen has been detected in urine as early as three days after the onset of symptoms. The test is rapid, giving a result within 15 minutes, and utilizes a urine specimen which is convenient for collection, transport, and subsequent detection of early, as well as later, stages of disease.

### PRINCIPLE

The *Legionella* MonlabTest® is an immunochromatographic membrane assay to detect *Legionella pneumophila* soluble antigen in human urine. Anti-*Legionella pneumophila* antibody, the test line, is adsorbed onto nitrocellulose membrane. Antibodies of the control line were adsorbed onto the same membrane as a second band. Anti-*Legionella pneumophila* antibodies are conjugated to visualizing particles that are dried onto an inert absorbent support.

During testing the sample is allowed to react with the conjugate which was pre-adsorbed on the strip test. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the coloured particles migrate. In the case of a positive result the specific antibodies present on the membrane will capture the coloured conjugate. *L. pneumophila* urinary antigen captured by immobilized anti- *L. pneumophila* antibody reacts to bind conjugated antibody. The other immobilized antibodies also captures visualizing conjugate, forming the control line. A positive test result is read visually in 15 minutes or less depending on the concentration of antigen present in the urine specimen. A negative *Legionella pneumophila* result, read in 15 minutes, indicates that *L. pneumophila* antigen was not detected in the urine sample.

The test is interpreted by the presence or absence of visually detectable red and blue colored lines. A positive result will include the detection of both a test and a control line, while a negative assay will produce only the control line. Failure of the control line to appear, whether the test line is present or not, indicates an invalid assay.

### PRECAUTIONS

- For professional *in vitro* diagnostic use only.
- Do not use components past its expiration date.
- Do not mix components from different kit lots.
- Solutions used to make the control swabs are inactivated using standard methods. However, patient samples, controls, and Device tests should be handled as though they could transmit disease. Observe established precautions against microbial hazards.
- The test should be discarded in a proper biohazard container after testing.

### STORAGE AND STABILITY

Store as packaged in the sealed pouch either at refrigerated or room temperature (2-30°C/36-86°F). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

### MATERIALS PROVIDED

- 20 *Legionella* MonlabTest® Tests
- Reagent (sample and controls diluent)
- Instructions for use
- 20 testing tubes
- 1 Positive Control Swab: Inactivated *L. pneumophila* swab + testing tube + pipette
- 1 Negative Control Swab: *L. pneumophila* negative swab + testing tube + pipette

### MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer

### SPECIMEN COLLECTION AND PREPARATION

Urine specimens should be collected in standard containers. The samples can be stored at room temperature (15-30°C/59-86°F) if assayed within 24 hours of collection. Alternatively, specimens may be stored at 2-8°C for up to 14 days or at -10°C to -20°C for longer periods before testing.

When necessary, urine specimens should be shipped in leak-proof containers at 2-8°C or frozen.

Allow all specimens to equilibrate to room temperature before testing.

### ASSAY PROCEDURE

#### Procedure for Patient Samples. See illustration 1:


Do not remove device test from pack until test sample has reached room temperature. Bring patient urine and/or liquid urine control(s) to room temperature (15-30°C/59-86°F).

1. Use a separate testing tube for each sample. Add 3 drops of urine sample and 1 drop of Reagent into the testing tube and mix with the pipette.
2. Remove the device test from its pack just before use. Use a separate test for each sample. Place the test on a flat surface
3. Dispense exactly 3 drops from the testing tube, into the circular window marked with S. Start the timer. Read the result at **15 minutes**.


#### Procedure for Positive and Negative Swab Controls. See illustration 2:

1. Add slowly 13 drops of Reagent into the testing tube.
2. Immediately remove swab control from the pouch and put the swab into the testing tube with reagent, mix 1 minute and extract as much liquid possible from the swab, squeezing the sides of the tube as the swab is withdrawn. Discard the swab.
3. Remove the device test from its sealed bag just before using and place the test on a flat surface. Dispense exactly 3 drops from the testing tube, into the circular window marked with S. Start the timer. Read the result at **15 minutes**. Repeat the procedure for Negative Swab Control using the same Reagent used for test and for positive control swab.

#### Illustration 1

#### Illustration 2

### INTERPRETATION OF RESULTS

Illustration 3

**POSITIVE** sample or positive control swab: two lines appear across the central window, in the result line region, **red** test line marked with the letter T, and in the control line region, **blue** control line marked with the letter C. This means that antigen was detected.

Recommended report: Presumptive positive for *L. pneumophila* antigen in urine, suggesting current or past infection.

#### CONFIRMATION OF POSITIVE SAMPLES

The positive results must be confirmed following the next protocol:

Exposed the positive sample at 100°C for 5 minutes, then centrifuged it for 15 minutes at 1000 rpms to obtain a supernatant (clear). Use this supernatant as a new sample and repeat the assay with a new test (following ASSAY PROCEDURE).

All the positive results should be confirmed with additional laboratory techniques (culture or serology).

**NEGATIVE** sample or negative control swab will give only the **blue** line across the control line region marked with the letter C (control line). The control line means that the detection part of the test was done correctly, but no *L. pneumophila* antigen was detected.

Recommended report: Presumptive negative for *L. pneumophila* antigen in urine, suggesting no recent or current infection. Infection due to *Legionella* cannot be ruled out since other serogroups and species may cause disease, antigen may not be present in urine in early infection, and the level of antigen present in the urine may be below the detection limit of the test.

**INVALID:** if no lines are seen, or if just the test line is seen, the assay is invalid. Invalid tests should be repeated. If the problem persists, discontinue using the test kit and contact your local distributor.

### NOTES ON THE INTERPREATION OF RESULTS

The intensity of the red coloured band in the result region (I) will vary depending on the concentration in the specimen. However, neither the quantitative value, nor the rate of increase in antigens can be determined by this qualitative test.

### QUALITY CONTROL

#### Internal Quality Control:

The *Legionella* MonlabTest® contains built-in control features. The blue line appearing in the control line region (C). It confirms sufficient specimen volume and correct procedural technique.

**Positive Procedural Control** The blue line at the Control line region can be considered an internal positive procedural control. If capillary flow has occurred, this line will always appear.

**Negative Procedural Control** The clearing of background color in the result window provides a negative background control.

**External Positive and Negative Controls:** Good Laboratory Practice recommends the use of positive and negative controls to assure functionality of reagents and proper performance of assay procedure. Positive and negative control swabs that will monitor the entire assay are provided in the kit. To use liquid urine controls, simply process as you would a patient sample. Positive and negative controls should be tested once for each new test kit opened and as otherwise required by your laboratory's standard quality control procedures.

### LIMITATIONS

- Legionella* MonlabTest® has been validated using urine samples only. Other samples (e.g. plasma, serum or other body fluids) that may contain *Legionella* antigen have not been evaluated. The test cannot be used on environmental samples (i.e. potable water).
- This test detects infections caused by *L. pneumophila* serogroup 1 and may detect other serogroups as well. A negative antigen result does not exclude infection with *L. pneumophila*. Culture is recommended for suspected pneumonia to detect causative agents other than *L. pneumophila* and to recover *L. pneumophila* when antigen is not detected in urine.
- The diagnosis of Legionnaires' disease cannot be based on clinical or radiological evidence alone. There is no single satisfactory laboratory test for Legionnaires's disease. Therefore, culture results, serology and antigen detection methods should be used in conjunction with clinical findings to make an accurate diagnosis.
- Excretion of *Legionella* antigen in urine may vary depending on the individual patient. Antigen excretion may begin as early as 3 days after onset of symptoms and persist for up to 1 year afterwards. A positive *Legionella*

MonlabTest® result can occur due to current or past infection and therefore is not definitive for infection without other supporting evidence.

- The test is compatible with boric acid, used as a preservative, up to a maximum of 2%.
- This test is a protein/biotin-based immunoassay (only included in control line preparation) that is designed to minimize the risk of biotin interference. Biotin supplements in patients do not interfere in this system.

### EXPECTED VALUES

The disease occurs in both epidemic and endemic forms and sporadic cases are not easily differentiated from other respiratory infections by clinical symptoms. An estimated 25,000 to 100,000 cases of *Legionella* infection occur in the United States annually. The resulting mortality rate, ranging from 25% to 40%, can be lowered if the disease is diagnosed rapidly and appropriate antimicrobial therapy is instituted early.

### PERFORMANCE CHARACTERISTICS

#### Sensitivity and specificity

The *Legionella* MonlabTest® was highly specific (>99%) and also sensitive (>99%) compared with the results of the Binax NOW® *Legionella* Urinary Antigen Test. The *Legionella* MonlabTest® was used to evaluate 149 samples for evaluation vs. BinaxNOW® *Legionella* Urinary Antigen Card, Binax.

#### Cross reactivity

It was performed an evaluation to determine the cross reactivity of *Legionella* MonlabTest®. There is no cross reactivity with other pathogens occasionally present in urine:

- *Streptococcus pneumoniae*

### REFERENCES

- Roig, J., X. Aquiler, J. Ruiz, et. al. Comparative study of *Legionella pneumophila* and other nosocomial-acquired pneumonias. Chest. 1991;99:344-50.
- Berdal, B.P., C.E. Farshy, and J.C. Feeley. Detection of *Legionella pneumophila* antigen in urine by enzyme-linked immunospecific assay. J. Clin. Microbiol. 1979;9:575-578.
- White A., et al. Rapad diagnosis of Legionnaires' disease. Trans Am Clin. Climatol. Assoc. 1982;93:50-62
- Bibb, W.F., P.M. Arnow, L. Thacker, and R.M. McKinney. Detection of soluble *Legionella pneumophila* antigens in serum and urine specimens by enzyme-linked immunosorbent assay with monoclonal and polyclonal antibodies. J. Clin. Microbiol. 1984;20:478-482.
- Tang, P.W., and S. Toma. Broad-spectrum enzyme-linked immunosorbent assay for detection of *Legionella* soluble antigens. J. Clin. Microbiol. 1986;24:556-558.
- Kohler, R.B., W.C. Winn, Jr., and L.J. Wheat. Onset and duration of urinary antigen excretion in Legionnaires' disease. J. Clin. Microbiol. 1984;20:605-607.

### SYMBOLS FOR IVD COMPONENTS AND REAGENTS

	Manufacturer		For <i>in vitro</i> diagnostic use only		Positive control swab
	Don't re-use		Consult instructions for use		Negative control swab
	Contains sufficient for <n> tests		Keep dry		
	Catalogue Code		Temperature limitation		
	Lot Number		Use by		

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- *Streptococcus pneumoniae*

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- Roig, J., X. Aquiler, J. Ruiz, et. al. Comparative study of *Legionella pneumophila* and other nosocomial-acquired pneumonias. *Chest*. 1991;99:344-50.
- Berdal, B.P., C.E. Farshy, and J.C. Feeley. Detection of *Legionella pneumophila* antigen in urine by enzyme-linked immunospecific assay. *J. Clin. Microbiol.* 1979;9:575-578.
- White A., et al. Rapad diagnosis of Legionnaires' disease. *Trans Am Clin. Climatol. Assoc.* 1982;93:50-62
- Bibb, W.F., P.M. Arnov, L. Thacker, and R.M. McKinney. Detection of soluble *Legionella pneumophila* antigens in serum and urine specimens by enzyme-linked immunosorbent assay with monoclonal and polyclonal antibodies. *J. Clin. Microbiol.* 1984;20:478-482.
- Tang, P.W., and S. Toma. Broad-spectrum enzyme-linked immunosorbent assay for detection of *Legionella* soluble antigens. *J. Clin. Microbiol.* 1986;24:556-558.
- Kohler, R.B., W.C. Winn, Jr., and L.J. Wheat. Onset and duration of urinary antigen excretion in Legionnaires' disease. *J. Clin. Microbiol.* 1984;20:605-607.

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