**CRYPTO-GIARDIA-ENTAMOEBIA MonlabTest®**  
**MO-076005 20 TESTS**

Immunochromatographic test for the detection of and *Cryptosporidium parvum*,  
*Giardia lamblia* and *Entamoeba histolytica* antigens in human stool samples

For *in vitro* use only. Store at 2 - 30°C.

**INTENDED USE**

Crypto-Giardia-Entamoeba MonlabTest chromatographic immunoassay is a procedure for *in vitro* qualitative detection of *Cryptosporidium parvum*, *Giardia lamblia* and *Entamoeba histolytica* antigens in human stool samples.

Cryptosporidiosis is a protozoal infection caused by *Cryptosporidium parvum* in humans. Typically, is an acute short-term infection, but can become severe and non-resolving in children and immunocompromised individuals such as AIDS patients. It may also be asymptomatic. In tropical developing countries, the parasite is often endemic and causes diarrhoea epidemics among children. With immunocompetent patients, the disease manifests itself as self-healing gastroenteritis.

Giardiasis in humans is caused by the protozoan parasite *Giardia lamblia*. This organism is involved in 25% of the cases of gastrointestinal disease and may be present asymptotically. *Giardia lamblia* has become an important cause of chronic diarrhoeas, particularly regarding travel medicine.

*Entamoeba histolytica* is the protozoan parasite responsible for dysentery and amebiasis. It is the third leading cause of morbidity and mortality due to parasitic disease in humans after malaria and schistosomiasis, and is estimated to be responsible for between 50000 and 100000 deaths every year. The disease may manifest itself as an acute, chronic or as an asymptomatic infection.

All these parasites present a simple life cycle that usually consists of an infective cyst stage and a multiplication trophozoite stage. Transmission of these infections occurs via ingestion of cysts, most often via food or water contaminated with human faecal matter.

Usually, diagnosis is performed by microscopic examination and requires experienced technicians. So it would be useful to develop an alternative method, rapid and easy-to-use, not requiring the presence of intact organisms (cysts or trophozoites) in the stool sample.

In the case of *Entamoeba*, a new understanding of this organism has led to the recognition that two species actually exist within what was previously known as *E. histolytica*. Of these two organisms *E. histolytica* is the pathogenic, causing all invasive disease while the other, *E. dispar* is non pathogenic as it is not capable of invading tissue. The two species are morphologically identical, so although diagnosis is usually performed by microscopic examination, it cannot differentiate them. Other techniques are required to detect specific antigens of each specie for an accurate diagnosis and to prevent unnecessary or inappropriate chemotherapy.

The test is based on the immunological capture of coloured microparticles as they move along a membrane on which the monoclonal antibody has been immobilized.

**PRINCIPLE**

Crypto-Giardia-Entamoeba Monlabtest use specific monoclonal antibodies against *Cryptosporidium parvum*, *Giardia lamblia* and *Entamoeba histolytica* that detect all forms of the parasites during their life cycle.

The tests are based on the use of three types of microspheres: blue microspheres covalently linked to a monoclonal anti-*Cryptosporidium parvum* antibody, red microspheres covalently linked to a monoclonal anti-*Giardia lamblia* antibody, green microspheres covalently linked to a monoclonal anti-*Entamoeba histolytica* antibody as well as purple microspheres used as test control.

The parasites present in stool samples react with the latex particles which are coated with specific monoclonal antibodies against the antigen. This *latex particles/antibodies/parasite* complex migrates through a chromatographic process towards the reaction area. In this area, anti-*Cryptosporidium*, anti-*Giardia* and anti-*Entamoeba* antibodies that react with the *latex particles/antibodies/parasite* complex are present. This reaction leads to the appearance of a blue and/or red and/or green lines. These lines are used to interpret the result, following a ten-minute room-temperature incubation.

**MATERIALS PROVIDED**

- 20 cassettes
- 1 Dilution buffer (30 ml)
- Instructions for use

**MATERIALS NOT INCLUDED IN THE KIT**

- Suitable applicators for sample collection (solid or liquid).
- Vortex mixer.
- Centrifuge adapted to sample extraction tubes.
- Test tubes or 1,5ml capped microtubes.
- Variable micropipettes.
- Pipette tips.
- Timer.
- Disposable gloves.

**PRECAUTIONS**

1. Patient samples (stool) may contain infectious agents and should be treated and disposed of as potentially dangerous biological materials.
2. The buffer contains sodium azide as antimicrobial agent. Avoid direct contact with skin and mucosa. Dispose of appropriately. Do not use the buffer if signs of contamination or precipitation are observed.
3. Do not store or prepare food, eat, drink or smoke in the area where the reagents and samples are handled.
4. Wear disposable gloves when handling the samples. Wash your hands thoroughly once you have finished working.
5. Do not exchange components from kits with different lot numbers.

6. All the reagents are for *in vitro* use exclusively.
7. Before using, let all kit components and samples reach room temperature, because cold reagents and/or samples can reduce test functionality. About 20-30 minutes are usually sufficient for reaching room temperature.
8. Do not use kit components beyond the expiration date.
9. If the package is broken, the product can still be used as long as none of the components have been damaged.
10. It is important to add the correct sample amount. If less amount than the indicated is used, the chromatography may not occur because the sample may not reach the reaction area; if an excess amount is used, brown lines may appear instead of blue, red, green or purple ones.
11. The used product should be disposed of as indicated by current legislation.
12. Do not use the test if a coloured line appears in the result area before its use.
13. It is important to take the appropriate amount of stool sample: 50 mg of solid stool or 100 µl of liquid stool to be extracted in 1 ml of diluent. It is very important to keep the correct sample: buffer ratio to ensure the correct performance of the test. **An excessive sample amount in regard to buffer amount prevents a correct chromatography.**
14. In order to ensure an adequate chromatography and optimum results, it is very important to centrifuge the 1.5 ml microtubes prior to extracting the specific quantity of supernatant. This is particularly true in the case of solid stool samples, as the greater number of suspended particles can interfere with the chromatography. However, if centrifugation is not possible, wait some minutes to settle solid particles.
15. Do not discard the external box of the kit until its content has been totally used. The external box contains essential information regarding the EC marking and component lots.

**STORAGE**

The kit can be stored at temperatures between 2°C and 30°C (it can be stored in a refrigerator). The expiration date is printed on the package.

**SAMPLES**

Stool samples should be collected in an appropriate sterile container as soon as possible after symptom onset. The sample should be representative, whenever possible.

The samples can be stored in the refrigerator (at approximately 4°C) for 1-2 days before being analysed. For a longer storage period, store the samples in a freezer at -20°C, with no further handling. In this case, the sample must be completely thawed, brought to room temperature and homogenized before analysis.

Avoid freezing and thawing the samples several times. Do not use specimens that have been collected or stored in transport medium or preserving agents like 10 % formalin, SAF, PVA, Ecofix, etc., because they interfere with the test.

Pay special attention when analyzing hemorrhagic samples as they often give false positive results when the blood content is high.

**PROCEDURE OF SAMPLE PREPARATION**

1. **Important:** Obtain samples from three different sampling sites at least, in order to get a sample as much representative as possible.
2. Dispense **1.0 ml of extraction buffer** in a properly labelled testing tube.
3. **Add** a portion of solid sample of approximately **50 mg** with a swab, a wooden applicator or a bacteriology loop (homogenize the sample before collecting). For **liquid or semi-solid stools** add **100 µl** of stool using an appropriate pipette.
4. **Shake** the test tube thoroughly by using a vortex mixer to assure proper mixing.
5. **Centrifuge** 5 minutes at 700 xg (approximately 3000 rpm in a benchtop centrifuge) to settle solid particles. If a centrifuge is not available, wait 3-5 minutes for the solid particles to settle at the bottom of the tube. In any case, please note that optimum test performance is achieved when centrifugation takes place (this is especially important with solid samples as is indicated in Precautions section).

**PROCEDURE**

1. Following preparation of the sample, take the reaction device out of the aluminum pouch. Discard the small bag of desiccant as it is not used during the test; its only purpose is to protect the test against moisture.
2. Add **125 µl of the supernatant** prepared in the "Procedure of sample preparation" section into the sample area of the reaction device (round window marked with an arrow). Do not add solid particles along with the liquid.
3. Incubate the test at room temperature for **10 minutes** and read the results.

**READING OF RESULTS**

**NEGATIVE:** Only a single **PURPLE** line appears in the result area, aligned with the letter "C" (Control) marked on the device frame. This line should always appear.

**POSITIVE:**

**Cryptosporidium:** A **PURPLE** line and a **BLUE** line appear in the result area, aligned with the number "1" (Test) marked on the device frame. The intensity may vary according to the antigen concentration present in the sample.

**Giardia lamblia:** A **PURPLE** line and a **RED/PINK** line appear in the result area, aligned with the number "2" (Test) marked on the device frame. The intensity can vary according to the antigen concentration present in the sample.

**Entamoeba histolytica:** A **PURPLE** line and a **GREEN** line appear in the result area, aligned with the number "3" (Test) marked on the device frame. The intensity can vary according to the antigen concentration present in the sample.

**INVALID:** the purple control band does not appear or the purple colour of this band is clearly altered, also nonspecific colours in the positive bands are invalid results. This indicates an anomalous test performance. Possible reasons are:

- Some reagents have got damaged or the test has expired.
- The sample was not prepared according to the instructions of use.
- High blood content in the sample.

In the event of an invalid result it is recommended that another test is run, strictly following the protocol described in this manual. For blood samples, it is advisable the use of an alternative technique because the problem of destabilization does not usually depend on the cassette but of the sample matrix.

Any line that may appear beyond 10 minutes due to the nature of the sample has no diagnostic value.

Final diagnosis should not be based only on the results from a single test. The diagnosis should be established by the correlation of test results with other appropriate data and with clinical symptoms.


Interpretation code for cassette images.

C: *Cryptosporidium parvum*

G: *Giardia lamblia*

E: *Entamoeba histolytica*

N: Negative

### QUALITY CONTROL

The test is invalid if no purple line appears, whether because the test was not properly performed or the reagents have deteriorated. If this happens, repeat the analysis, following the working protocol indicated in these instruction sheets closely.

**WARNING:** Including a control with an established result is recommended to ensure that the data obtained are correct.

### PROCEDURE LIMITATIONS

- The test should be used only for detecting *Cryptosporidium parvum*, *Giardia lamblia* and *Entamoeba histolytica* antigens in human stool.
- This is a qualitative test so no quantitative interpretation of the results should be made with respect to the intensity of the positive line.
- Over 200 samples were evaluated to ensure proper test performance. The correlation of results with other techniques (ELISA) was good. However, interferences in the test performance cannot be excluded.
- If there is an excessive sample amount, brown lines may appear instead of the blue, red, green and purple ones, or test and control lines may be absent because chromatography is not correctly conducted. The brown lines have no diagnostic value. If this happens, repeat the test with a smaller amount of stool sample or dilute the previously prepared extract.
- Due to the homology between *E. histolytica* and *E. dispar*, certain cross-reaction may occur between these two species but it is expected to be less than a 5%. No crossreaction with other substances has been observed during test evaluation.
- A negative result does not totally exclude the possibility of *Cryptosporidium parvum*, *Giardia lamblia* and/or *Entamoeba histolytica* infection. The significance of the results should be evaluated in relation to the patient's clinical symptoms.
- It is observed that faecal samples with high blood content may interfere negatively with the test. In these cases, specificity problems (false positive results) may appear with negative samples.
- The analysis of certain samples can produce lines of imprecise colour, mostly corresponding to negative samples. The test should be repeated if lines of indeterminate colour appear. If the same result is obtained again, the analysis should be performed by another analytic method.

### DIAGNOSTIC SENSITIVITY AND SPECIFICITY

Monlab immunochromatographic tests designed for *Cryptosporidium parvum*, *Giardia lamblia* and *Entamoeba histolytica*: Crypto, Giardia, Crypto/Giardia, Entamoeba, and Crypto/Giardia/Entamoeba were externally and internally evaluated against different reference techniques, depending on the case, (Microscopy, ELISA and/or PCR) several times. The results were as follows:

*Cryptosporidium parvum*:

Crypto MonlabTest		Reference techniques	
		Microscopy	ELISA Crypto Ridascreen
MonlabTest Crypto/ Giardia/Entamoeba	Sensitivity	80.0%	79.3%
	95% Confidence Interval	71.1-87.2	72.0-85.5
	Specificity	96.9%	99.5%
	95% Confidence Interval	95.0-98.2	98.5-99.9

*Giardia lamblia*:

Giardia MonlabTest		Reference techniques	
		Microscopy	ELISA Giardia Ridascreen
MonlabTest Crypto/ Giardia/Entamoeba	Sensitivity	92.4%	93.8%
	95% Confidence Interval	87.3-95.9	89.2-96.9
	Specificity	98.6%	98.9%
	95% Confidence Interval	97.3-99.3	97.6-99.6

*Entamoeba histolytica*:

Entamoeba MonlabTest		Reference techniques	
		ELISA Entamoeba Ridascreen	PCR E. histolytica
MonlabTest Crypto/ Giardia/Entamoeba	Sensitivity	82.4%	75.0%
	95% Confidence Interval	71.2-90.5	53.33-90.2
	Specificity	96.4%	89.8%
	95% Confidence Interval	94.1-97.9	86.4-92.6

### REPEATABILITY

INTRA-ASSAY PRECISION

Triplicates for each concentration of the sensitivity curve were assayed with one lot of the product, obtaining the same results.

### REPRODUCIBILITY

INTER-DAY PRECISION: Using one product lot, triplicates of the sensitivity curve were performed over 10 consecutive days by the same person. Only one difference, of less than one 1:2 dilution, was found, which is acceptable and tolerable for a qualitative immunochromatography technique.

INTER-LABORATORY PRECISION: Five different laboratory technicians assayed the same samples, obtaining high precision and concordance. Only a difference of one 1:2 dilution was found, which is acceptable and tolerable for a qualitative immunochromatography technique.

INTER-LOT PRECISION: A sensitivity curve was performed using three product lots. The analyses were performed by the same person on the same day. Only one difference, of less than one 1:2 dilution, was found, which is acceptable and tolerable for a qualitative immunochromatography technique.

### HOOK EFFECT

In the case of *Crypto* and *Giardia*, a standard concentration was not available. However, analysis of the sensitivity curve at the maximum internal standard concentration available, did not show the existence of any hook effect.

Concerning *Entamoeba histolytica* the prozone or hook effect of the test was studied using a quantified number of cysts. The results shown that there is IC prozone but is not inhibitory, at least, at the concentrations studied.

The analysis shown a lower reactivity around  $9 \times 10^6$  cysts/ml of *E. histolytica* HM-1. This cyst concentration is approximately 640 times the concentration of the sensitivity limit of the test (14000-28000 cysts/ml).

### INTERFERING SUBSTANCES

The substances indicated in the table, at the specified concentration, did not interfere with the results. Only Ibuprofen\* at higher concentration could interfere in the detection of a very low *Giardia* positive sample. A single lot was used to conduct the study.

Metronidazole	1.2 mg/ml
Racecadotril	0.72 mg/ml
Loperamide	24 µg/ml
Braun Atropine	1.68 µg/ml
Cimetidine	1.92 mg/ml
Omeprazole	0.144 mg/ml
Neomycin	6.62 mg/ml
Ampicillin	7.2 mg/ml
Ibuprofen*	5.76 mg/ml
Acetylsalicylic Acid	4.8 mg/ml
Sucrose	8 mg/ml
Blood	10% (v/v) →The sample was prepared by adding fresh blood directly on the matrix.

INTERFERING MICROORGANISMS

The microorganisms indicated below did not cause any interference in the results.

**Bacteria/Parasites/Viruses**


<i>Entamoeba coli</i>	<i>Iodamoeba bütschlii</i>	<i>Hymenolepis nana</i>
<i>Giardia lamblia</i>	<i>Cryptosporidium parvum</i>	<i>Helicobacter pylori</i>
<i>Blastocystis hominis</i>	<i>Endolimax nana</i>	<i>Clostridium difficile</i>
<i>Entamoeba hartmannii</i>	Rotavirus	Adenovirus

**REFERENCES**

1. Goñi, P. et al., Evaluation of an immunochromatographic dip strip test for simultaneous detection of *Cryptosporidium* spp, *Giardia duodenalis*, and *Entamoeba histolytica* antigens in human faecal samples. *Eur J Clin Microbiol Infect Dis*. 31(8):2077-82 (2012).
2. Roka, M. et al., Prevalence of intestinal parasites in HIV-positive patients on the island of Bioko, Equatorial Guinea: Its relation to sanitary conditions and socioeconomic factors. *Sci Total Environ*. 432: 404-11 (2012).
3. Weitzel, T. et al., Evaluation of seven commercial antigen detection tests for *Giardia* and *Cryptosporidium* in stool samples. *Clin Microbiol Infect*. 12:656-659 (2006).
4. Lake, I.R. et al., Case-control study of environmental and social factors influencing cryptosporidiosis. *Eur J Epidemiol*. 22:805-811 (2007).
5. Khalakdina, A. et al., Is drinking water a risk factor for endemic cryptosporidiosis? A case-control study in the immunocompetent general population of the San Francisco Bay Area. *BMC Public Health*. 3: 1471-2458 (2003).
6. Cacciò, S.M. et al., Unravelling *Cryptosporidium* and *Giardia* epidemiology. *Trends Parasitol*. 21: 430-437 (2005).
7. Certad, G. et al., *Cryptosporidium parvum*, a potential cause of colic adenocarcinoma. *Infect. Agent. Cancer* 2: 22 (2007).
8. Nygard, K. et al., A large community outbreak of waterborne giardiasis-delayed detection in a non-endemic urban area. *BMC Public Health*. 6: 141 (2006).
9. Monis, P.T., and Thompson, R.C.A. *Cryptosporidium* and *Giardia*-zoonoses: fact or fiction?. *Infect Genet. Evol*. 3: 233-244 (2003).
10. Ekdahl, K. and Andersson, Y. Imported giardiasis: impact of international travel, immigration, and adoption. *Am. J. Trop. Med*. 72: 825-830 (2005).
11. Ungar, B. L. P. et al., Use of a monoclonal antibody in an enzyme immunoassay for the detection of *Entamoeba histolytica* in fecal specimens. *Am. J. Trop. Med. Hyg*. 34: 465-472 (1985).
12. Strachan, W. D. et al., Immunological differentiation of pathogenic and non-pathogenic isolates of *Entamoeba histolytica*. *The Lancet* 12: 561 - 563 (1988).
13. Yau, Y. C. W. et al., Development of monoclonal antibodies which specifically recognise *Entamoeba histolytica* in preserved stool samples. 39: 716-719.
14. Braga, L.L. Seropositivity for intestinal colonization with *Entamoeba histolytica* and *Entamoeba dispar* in individuals in Northeastern Brazil. 36: 3044-3045.
15. Petri Jr. W.A., et al., Pathogenic and nonpathogenic strains of *Entamoeba histolytica* can be differentiated by monoclonal antibodies to the galactose-specific adherence lectin. 58: 1802-1806 (1990).
16. Paniagua, G.L. et al., Two or more enteropathogens are associated with diarrhoea in Mexican children. 28: 17 (2007).
17. Haque, R. et al., *Entamoeba histolytica* infection in children and protection from subsequent amebiasis. 74: 904-909 (2006)

**SYMBOLS FOR IVD COMPONENTS AND REAGENTS**

	Manufactured by		For in vitro diagnostic use
	Do not re-use		Please read pack insert
	Contains sufficient for <n> tests		Dry storage
	Catalogue number		Store at
	Lot number		Expiry date
	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices		Dilution buffer
	Reagent		

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
- Goñi, P. et al., Evaluation of an immunochromatographic dip strip test for simultaneous detection of *Cryptosporidium* spp, *Giardia duodenalis*, and *Entamoeba histolytica* antigens in human faecal samples. *Eur J Clin Microbiol Infect Dis.* 31(8):2077-82 (2012).
- Roka, M. et al., Prevalence of intestinal parasites in HIV-positive patients on the island of Bioko, Equatorial Guinea: Its relation to sanitary conditions and socioeconomic factors. *Sci Total Environ.* 432: 404-11 (2012).
- Weitzel, T. et al., Evaluation of seven commercial antigen detection tests for *Giardia* and *Cryptosporidium* in stool samples. *Clin Microbiol Infect.* 12:656-659 (2006).
- Lake, I.R. et al., Case-control study of environmental and social factors influencing cryptosporidiosis. *Eur J Epidemiol.* 22:805-811 (2007).
- Khalakdina, A. et al., Is drinking water a risk factor for endemic cryptosporidiosis? A case-control study in the immunocompetent general population of the San Francisco Bay Area. *BMC Public Health.* 3: 1471-2458 (2003).
- Cacciò, S.M. et al., Unravelling *Cryptosporidium* and *Giardia* epidemiology. *Trends Parasitol.* 21: 430-437 (2005).
- Certad, G. et al., *Cryptosporidium parvum*, a potential cause of colic adenocarcinoma. *Infect. Agent. Cancer* 2: 22 (2007).
- Nygard, K. et al., A large community outbreak of waterborne giardiasis-delayed detection in a non-endemic urban area. *BMC Public Health.* 6: 141 (2006).
- Monis, P.T., and Thompson, R.C.A. *Cryptosporidium* and *Giardia*-zoonoses: fact or fiction?. *Infect Genet. Evol.* 3: 233-244 (2003).
- Ekdahl, K. and Andersson, Y. Imported giardiasis: impact of international travel, immigration, and adoption. *Am. J. Trop. Med.* 72: 825-830 (2005).
- Ungar, B. L. P. et al., Use of a monoclonal antibody in an enzyme immunoassay for the detection of *Entamoeba histolytica* in fecal specimens. *Am. J. Trop. Med. Hyg.* 34: 465-472 (1985).
- Strachan, W. D. et al., Immunological differentiation of pathogenic and non-pathogenic isolates of *Entamoeba histolytica*. *The Lancet* 12: 561 - 563 (1988).
- Yau, Y. C. W. et al., Development of monoclonal antibodies which specifically recognise *Entamoeba histolytica* in preserved stool samples. 39: 716-719.
- Braga, L.L. Seropositivity for intestinal colonization with *Entamoeba histolytica* and *Entamoeba dispar* in individuals in Northeastern Brazil. 36: 3044-3045.
- Petri Jr. W.A., et al., Pathogenic and nonpathogenic strains of *Entamoeba histolytica* can be differentiated by monoclonal antibodies to the galactose-specific adherence lectin. 58: 1802-1806 (1990).
- Paniagua, G.L. et al., Two or more enteropathogens are associated with diarrhoea in Mexican children. 28: 17 (2007).
- Haque, R. et al., *Entamoeba histolytica* infection in children and protection from subsequent amebiasis. 74: 904-909 (2006)

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