C. DIFFICILE Ag (GDH) MonlabTest®

MO-076012 20 TESTS

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- Rupnik M, Dupuy B, Fairweather NF, Gerding DN, Johnson S, Just I, Lyrer DM, Popoff MR, Rood JJ, Sonenshein AL, Thelestam M, Wren BW, Wilkins TD and von Eichel-Streiber C. Revised nomenclature of *Clostridium difficile* toxins and associated genes. J Med Microbiol. Feb 2005; Vol: 54(Pt 2) pp: 113-117.
- Vonberg RP, Reichardt C, Behnke M, Schwab F, Zindler S and Gastmeier P. Costs of nosocomial *Clostridium difficile* associated diarrhoea. J Hosp Infect. Sep 2008; Vol: 70(1) pp:15-20.
- Finn E, Andersson FL and Madin-Warburton M. Burden of *Clostridioides difficile* infection (CDI) – a systematic review of the epidemiology of primary and recurrent CDI. BMC Infectious Diseases 2021; Vol: 21 pp: 456.
- Kuehne SA, Cartman ST, Heap JT, Kelly ML, Cockayne A and Minton NP. The role of toxin A and toxin B in *Clostridium difficile* infection. Nature. Oct 2010; Vol: 467(7316) pp: 711-713.
- Shetty N, Wren MW, Coen PG. The role of glutamate dehydrogenase for the detection of *Clostridium difficile* in faecal samples: a meta-analysis. J Hosp Infect. 2011. Jan 77 (1): 1-6.
- Cohen SH, Gerding DN, Johnson S, Kelly CP, Loo VG, McDonald LC, Pepin J, Wilcox MH; Society for Healthcare Epidemiology of America; Infectious Diseases Society of America. Clinical practice guidelines for *Clostridium difficile* infection in adults: 2010 update by the society for healthcare epidemiology of America (SHEA) and the infectious diseases society of America (IDSA). Infect. Control Hosp. Epidemiol. 2010. May; 31 (5): 431-55.
- Wilcox MH, Planche T, Fang FC and Gilligan P. What is the current role of algorithmic approaches for diagnosis of *Clostridium difficile* infection. J Clin Microbiol. 2010 Dec 48 (12): 4347-53
- Wren MW, Sivapalan M, Kinson R and Shetty NR. Laboratory diagnosis of *Clostridium difficile* infection. An evaluation of tests for faecal toxin, glutamate dehydrogenase, lactoferrin and toxigenic culture in the diagnostic laboratory. Br J Biomed Sci. 2009; 66(1): 1-5.

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C. DIFFICILE Ag (GDH) MonlabTest®

MO-076012 20 TESTS

One-step immunochromatographic test for the qualitative detection of Glutamate Dehydrogenase (GDH) from *C. difficile* in human faeces

For *in vitro* use only. Store at 2 - 30°C.

INTRODUCTION

Intended use

A single-step immunochromatographic test to aid in the diagnosis of the gastrointestinal infection caused by *Clostridium difficile*, by means of the qualitative detection of Glutamate Dehydrogenase (GDH) from *Clostridium difficile* in human stool samples.

The results obtained with C. DIFFICILE Ag (GDH) MonlabTest should be evaluated by the clinician in combination with other available clinical information from the patient. They should not be used as the only criterion for diagnosing the infection.

The test is qualitative and its execution is not automated, although the results interpretation can be done visually or by using a strip reader.

This test is exclusively for *in vitro* diagnostic use.

This test is for professional use. This is not a diagnostic test for use at the patient care location. This is not a self- diagnostic test.

General information

The test is based on the immunological capture of coloured microparticles as they pass through a membrane on which specific monoclonal antibodies against GDH from *C. difficile* have been immobilised.

A positive signal in the test indicates the presence of *C. difficile* in the stool sample. Subsequent analysis of the sample with a test for detecting *C. difficile* toxins (such as *C. difficile* toxins A+B MonlabTest) confirms whether the strain is toxigenic and, therefore, pathogenic and disease-causing.

The testing population

C. difficile infection can affect all types of people regardless of age, sex or lifestyle, so a priori no sector of the population susceptible to be analysed with this test is excluded.

In those medical centres in which the incidence of nosocomial diseases is higher, the use of this test should be encouraged, as *C. difficile* is the main agent related to the infections acquired when being at hospital. In the same way, any diarrhoeic sample that comes from individuals under antibiotic treatment should be evaluated with this test.

Incidence of the disease

C. difficile can be present asymptotically in up to 5% of the healthy population¹. This carrier rate can increase up to 30% as a result of hospitalisation.

Besides, *C. difficile* is considered to be responsible for approximately 25% of the diarrhoea associated with the consumption of antibiotics and it is one of the main pathogens related to nosocomial infections².

Incidence worldwide oscillates between 1.1 and 631.8 for every 100,000 persons a year. Furthermore, there is clear evidence for disease recurrence, which is around 15-35% for cases with a *C. difficile* infection³.

Characteristics of *C. difficile* and its infection

Clostridium difficile is a spore-forming gram-positive anaerobic bacterium.

In addition to the diarrhoea symptoms, the disease can lead to pseudomembranous colitis, which requires urgent treatment with antibiotics effective against *C. difficile* (metronidazole, vancomycin) as the patient's life may be compromised. Mortality associated with *C. difficile* infection can be between 6% and 30%, especially if the patient suffers from pseudomembranous colitis.

C. difficile releases two high molecular weight toxins⁴, toxin A and toxin B, which are the main virulence factors responsible for the clinical manifestations of disease caused by *C. difficile* infection.

Glutamate dehydrogenase (GDH) is an enzyme produced in large quantities by toxigenic and non-toxigenic strains of *C. difficile*, thus making it an excellent marker for determining the presence of this microorganism⁵.

The Society for Healthcare Epidemiology of America (SHEA), the Infectious Diseases Society of America (IDSA) and the European Society of Clinical Microbiology and Infectious Diseases (ESCMID) propose the use of a two- step protocol for identifying toxigenic *C. difficile*⁶ involving an initial screening of the sample with a GDH test and, if the result is positive, use of a second test to detect toxins from *C. difficile*, as not all strains of this microorganism produce them. MONLAB offers the possibility of following this protocol recommended by some of the leading healthcare societies using its GDH and Toxins A+B (for toxin presence) tests.

BASIC PRINCIPLES OF THE TEST

The GDH test uses a combination of:

- red latex particles conjugated to specific antibodies against GDH from *C. difficile* that cooperate with other specific anti-GDH antibodies located on the test membrane under the control band,
- blue latex particles conjugated to an antigen recognised by a specific antibody for that antigen bound to the membrane; this is the control band.

In this test the sample is first treated with the sample dilution buffer (included in the kit) to extract GDH from the faecal matrix. After the extraction, a specific volume of the supernatant is added to the reactive strip.

When the extracted sample flows through the test membrane the coloured particles migrate. In the case of a positive sample, the specific antibodies present in the membrane capture the coloured antigen-coated particles. Different coloured lines (blue for control and red for test) will be visible depending on the content of GDH in the sample. These lines are used to interpret the result after incubation at room temperature for 10 minutes (see Fig.1).

MATERIALS PROVIDED WITH THE KIT

- 20 cassettes
- 20 vials with diluent (1,5mL)
- Disposable plastic pipettes
- Instructions for use

MATERIALS NOT INCLUDED IN THE KIT

- Vortex
- Timer

Controls are available for the user to validate the results and can be purchased as a separate trade reference.

EQUIPMENT UNITS FOR THE RESULT INTERPRETATION

The test can be interpreted visually or using iPeak® Plus immunochromatographic test reader. For this diagnosis, there is evidence of a 97.6% agreement between visual interpretation and the iPeak® Plus strip reader

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PRECAUTIONS

1. Patient specimens (faeces) should be handled with caution as they may contain infectious agents. The use of disposable gloves throughout the entire procedure is required.
2. The sample dilution buffer contains sodium azide as anti-microbial agent. Avoid direct contact with skin and mucous membranes. Dispose of appropriately. Do not use the buffer if signs of contamination or precipitation are present. The sample dilution buffers contains material from animal origin (albumins). All the materials from animal origin are commercial references and have the certificate that evidences they are not infectious.
3. Do not eat, drink, smoke or prepare or store food in areas where reagents or samples are handled.
4. Once completed the procedure, clean work surfaces with soap and water and finish disinfecting with a suitable solution. Finally, discard the gloves and then wash hands first with soap and water rubbing them well.
5. Do not exchange components from one kit lot to another.
6. Before use, allow all kit components and stool samples to reach room temperature, as cold reagents and/or samples can reduce test functionality. Twenty to 30 minutes are recommended to reach room temperature.
7. All reagents are for *in vitro* use exclusively.
8. The immunochromatographic strips are single use. A strip is used per sample to be evaluated.
9. Do not use kit components beyond their expiration date.
10. In case of external package damage, the product may be used if none of its components has been damaged.
11. In case of damage or alteration of the primary packaging discard the test.
12. Open the primary packaging of the strip, the aluminium pouch when the assay is going to be performed in order to unnecessarily expose the test in excess to external environmental factors that may work against it.
13. Is very important to add the correct volume of extracted sample to the reaction device. If the amount is lower than required, the chromatography may not be completed as the sample does not reach the reaction area; if higher than indicated, excessive volume could spill inside the device, and the chromatography may not develop correctly.
14. Once it is used, it is advisable to manipulate the strip taking into account the same considerations as far the sample. The strips have to be manipulated with the proper precautions and should be managed as biohazardous material.
15. Used product should be disposed of according to applicable laws.
16. Do not use the test if any coloured lines appear in the result area before using the test.
17. It is critical to collect the correct sample quantity: approximately **110 mg of a solid sample (a small ball of 5 mm in diameter)**. If the sample is **semi-liquid** (unable to take it with a pipette), take a sample amount so that it completely **covers the grooves of the stick** attached to the vial cap and **110 µL of a liquid sample (4 drops)** if the disposable pipettes included in the kit are used). These quantities are extracted in the sample diluent supplied with the vials provided in the kit. An excess of sample in relation to the amount of buffer added prevents the chromatography from running correctly; this is especially critical in the case of solid samples, since it is harder to obtain an appropriate quantity.
18. Do not discard the external box of the kit until its content has been totally used. The external box contains essential information regarding the EC marking and component lots.
19. Any serious incident related to the product must be reported to the manufacturer and the competent authority of the Member State in which the user and / or the patient are established.

STORAGE AND STABILITY

Product may be stored at any temperature between +2 and +30°C.

The expiry date is printed on the aluminium packages. The reagent strips must be used once the protective packaging has been opened after reaching room temperature (if stored refrigerated). The product maintains the claimed stability once open.

SAMPLES

- This test is designed for the analysis of human faecal samples.
- Do not use samples collected in transport media or added with preservatives (such as formalin, SAF, PVA or similar) or enrichment media as their presence could interfere with the proper performance of the test.
- The use of fresh untreated samples is recommended. If they have to be stored for some time, they can be stored in the refrigerator (+2-8°C) for 3 or 4 days. For a longer storage, they should be frozen at -20°C, considering that some samples lose immunoreactivity after being frozen.

- Pay special attention when analysing haemorrhagic samples since they may cause non-specificity problems if the blood content is high. An indication of this instability of the test is usually a change in the blue colour of the control band (tends to show a purple or very dark blue colour) or the appearance of other bands of an undefined colour.
- In case of freezing, thaw the samples completely at room temperature prior to analysis.
- Avoid freezing-thawing cycles for the faecal samples as the immunological recognition of the antigen can be altered.
- Samples must be well homogenised prior to preparation.

PREPARATION OF SAMPLES

General Note: disposable gloves should be used throughout the test procedure due the handling of infectious samples. Once the work is concluded, do not forget to comply with the hygiene procedures detailed in point 4 of the "Precautions" section.

The protocol for faecal samples preparation is:

- Homogenize previously the sample to get an aliquot as much representative as possible.
- Unscrew the cap from the vial with caution in order not to spill the sample diluent buffer. With **solid samples**, take with the stick attached to the vial cap an approximate amount of **110 mg** of faeces (a small portion of **5 mm in diameter**). If the sample is **semi-liquid** (unable to take it with a pipette), take a sample amount so that it completely **covers the grooves of the stick** attached to the vial cap. With **liquid samples**, take a volume of **110 µL (4 drops)** if the disposable pipettes included in the kit are used).
- Carefully add the sample into the vial containing the dilution buffer. Screw the cap well and shake vigorously to ensure a homogeneous mixture.

PROCEDURE

Once the samples have been prepared as described in previous sections, the procedure is as follows:

- Remove the reaction device from the aluminium pouch. Discard the desiccant bag as it only serves to preserve the test from moisture.
- Put the vial upside down and add **3 drops** to the sample area of the reaction device (round window marked with an arrow).
- Wait for **10 minutes** and read the results (see Fig. 1).

READING OF RESULTS

The strips shown in Figure 1 are an example of the different results that can be obtained with *C. DIFFICILE* Ag (GDH) MonlabTest. 2 different coloured bands can be distinguished:

- Upper (blue) line:** is the control line indicating a proper test performance.
- Lower (red) line:** indicates the presence of GDH in the sample.

The control blue line should always appear. The additional presence of the red line indicates the presence of GDH in the sample.

- Strip 1:** NEGATIVE result. A single transverse **BLUE** line appears aligned with the letter "C" marked on the plastic cassette or casing. This is the control line and should always appear indicating that chromatography takes place normally.
- Strip 2:** POSITIVE result. A **BLUE** line (control) and a **RED** line aligned with the label "T" marked on the casing appear. The line intensity depends on the analyte concentration in the sample.


Fig. 1: Possible results of the GDH test

- Strip 3:** INVALID results. The control line does not appear, the **BLUE** colour of this line is clearly altered or non-specific colours (other than red) appear in the positive line. Any combination of colours other than those indicated above indicates abnormal functioning of the test. Possible reasons are:
 - some reagents have got damaged or the test has expired.
 - the sample was not prepared according to the instructions of use.
 - the sample has a high blood content.

In the event of an invalid result, it is recommended to run another test, strictly following the protocol described in this manual.

In the case of samples with blood, it is recommended to use an alternative technique if the problem persists as this usually depends on the sample matrix itself rather than on the strip used.

Any line that appears beyond 10 minutes, has no diagnostic value.

NOTE: the final and definitive diagnosis is established by the clinician. This test only detects GDH from *C. difficile* in a sample but is not an argument to conclude that the person has been infected by *C. difficile*.

QUALITY CONTROL

The test is invalid if no blue line appears, whether because the test was not performed properly, or the reagents have deteriorated. If this happens, repeat the analysis, following the work protocol indicated in these instruction sheets closely.

WARNING: Including controls with an established result is recommended to ensure that the data obtained are correct.

LIMITATIONS OF THE PROCEDURE

- C. DIFFICILE* Ag (GDH) MonlabTest is intended to aid in the diagnosis of gastrointestinal infection caused by *C. difficile* by detecting the presence of Glutamate Dehydrogenase (GDH) in human faeces as long as the protein load is equal or above the detection limit of the test.
- This is a qualitative, not quantitative test, although the intensity of positive bond relates to the amount of GDH in the stool sample.

- The test has shown good correlation with other techniques (ELISA and rapid test) after analysing a large amount of samples. However, this study does not exclude potential interference in test performance when analysing other faecal samples.
- It is important to control reaction time. If reaction time is less than that indicated, samples containing an analyte concentration significantly higher than the sensitivity limit can be seen clearly, but those that are at the limit will not appear. If the reaction time is longer than indicated, the sensitivity of the test will be affected and can give rise to incorrect interpretations.
- A defect in sample amount may lead to very weak positive results. In this case, the test must be repeated with a larger quantity of sample, maintaining the recommended proportion with respect to sample diluent. An excess of sample can cause the test to develop very slowly and even prevent the test from developing correctly (the control bond is not visible). In this case, the test must be repeated with a smaller quantity of sample (see section "Stool sample preparation").
- A negative result does not exclude the possibility of infection by *C. difficile*. It may be the case that the antigen concentration in the sample is lower than the limit of detection for the test or that antigen-binding substances or enzymes in the faeces inactivate the antigen.
- A positive result must be compared with other techniques that analyse the presence of toxins to ensure that the *C. difficile* strain is virulent. A positive result does not exclude the presence of other pathogens.
- The test results should be interpreted alongside available data from epidemiological studies, the patient's clinical evaluation, and other diagnostic procedures.
- It has been observed that some stool samples with a high blood content may interfere negatively with the test, with possible non-specificity problems for samples that are negative for *C. difficile*. This instability of the test is usually accompanied by a change in the colour of the control bond; instead of a light blue colour, a very dark blue or purple colour appears (see the "Reading the Results" Section).
- C. difficile* colonisation rates of up to 50% have been reported in infants. A high rate has also been reported in cystic fibrosis patients. Normally, these two groups of patients remain asymptomatic and do not require specific treatment. These positive results that are of no clinical significance should be treated with caution as affected patients do not require treatment for *C. difficile* infection.
- The GDH screening as a first step for detecting toxigenic *C. difficile* strains presents a sensitivity limitation as it is indicated by various papers^{7,8}, in which false negative rates up to 5% are indicated.
- The performance of the test has not been established to verify the efficacy of the eradication (test of cure).

DETECTION CAPABILITY

A pure preparation of native glutamate dehydrogenase obtained from the supernatant of a *C. difficile* culture from tgcBionics was used. *C. DIFFICILE* Ag (GDH) MonlabTest interpreted visually and with the IUL Lateral Flow Reader, achieves a detection capacity of 1,5 ng/mL.

ANALYTICAL SPECIFICITY

Monoclonal antibodies used in GDH test offer a unique specificity for the detection of the Glutamate dehydrogenase (GDH) from *C. difficile*, without expected cross-reactions with other antigens and/or micro-organisms present in stool samples have been detected. See the cross-reactivity section.

DIAGNOSTIC SENSITIVITY AND SPECIFICITY

1. Internal Evaluation

A comparison between the rapid test of MONLAB and the ELISA *C. DIFF* CHEK™-60 of Techlab was done. The study included a total of 334 samples, of which 105 were positive and 229 were negative by ELISA.

The results obtained by visual interpretation of the test were as follows:

GDH		ELISA <i>C. DIFF</i> CHEK™-60	
		Positive	Negative
<i>C. DIFFICILE</i> Ag (GDH) MonlabTest	Positive	103	0
	Negative	2	229
Sensitivity concordance = $103/(103+2) \times 100 = 98,1\%$			
Specificity concordance = $229/(0+229) \times 100 > 99,9\%$			

2. External evaluation

A comparison between the rapid test MonlabTest and the rapid test *C. DIFF* QUIK CHEK® of Techlab was done in a hospital in Barcelona. The study included a total of 115 samples, of which 42 were positive and 73 were negative by the Techlab rapid test.

The results obtained by visual interpretation of the test were as follows:

GDH		<i>C. DIFF</i> QUIK CHEK®	
		Positive	Negative
<i>C. DIFFICILE</i> Ag (GDH) MonlabTest	Positive	42	3
	Negative	0	70
Sensitivity concordance = $42/(42+0) \times 100 > 99,9\%$			
Specificity concordance = $70/(70+3) \times 100 = 95,9\%$			

The 3 discordant samples were re-analysed by the *C. difficile* anaerobic culture, verifying that two of the false positives of the Monlab test are in fact positive samples.

REPEATABILITY

- INTRA-ASSAY PRECISION: Five replicates from a series of 1/2 dilutions of a characterised internal standard were analysed visually and with the IUL Lateral Flow Reader. The measurements were taken using 1 batch of the *C. DIFFICILE* Ag (GDH) MonlabTest product on the same day by the same person. 100% repeatability was obtained, indicating a high intra-assay precision of the test.

REPRODUCIBILITY

- INTER-DAY PRECISION: Using the same batch of *C. DIFFICILE* Ag (GDH) MonlabTest, a dilution curve of the standard was measured visually over five days spaced in time. The results were very reproducible as differences of less than ½ dilution were observed, acceptable and tolerable for the assay.
- INTER-OPERATOR / INTER-LABORATORY PRECISION: Five persons measured visually a dilution curve of the standard in duplicate. The differences were not observed to be more than one 1/2 dilution in any case, which is acceptable and tolerable for the technique.
- INTER-LOT PRECISION: Using four batches of *C. DIFFICILE* Ag (GDH) MonlabTest, a dilution curve of the standard was measured visually in duplicate. The analysis was performed by one person on the same day. The differences were observed to be equal than one 1/2 dilution, which is acceptable and tolerable for the technique.

THE HOOK OR PROZONE EFFECT

There is evidence that *C. DIFFICILE* Ag (GDH) MonlabTest does not present a Hook effect when evaluating increasing concentrations of the standard up to 4,000 ng/mL.

INTERFERING SUBSTANCES

The following substances did not affect the test visual results when added to stool samples (positive and negative) at the concentrations shown in the following table:

Racecadotril	5% (p/v)	Ibuprofen	20% (p/v)
Cimetidine	10% (p/v)	Acetylsalicylic acid	30% (p/v)
Loperamide	5% (p/v)	Edulcorant	5% (p/v)
Metronidazole	5% (p/v)	Palmitic Acid	40% (p/v)
Omeprazole	3% (p/v)	Barium Sulphate	5% (p/v)
Ampicillin	15% (p/v)	Mucin	5% (p/v)
Human blood	20% (v/v)		

CROSS-REACTIVITY

This study was developed in two different places, with visual interpretation of the results:

- MONLAB:** *C. DIFFICILE* Ag (GDH) MonlabTest was evaluated against strongly positive samples for the following micro-organisms, without finding problems of cross-reactivity: Adenovirus, Rotavirus, Norovirus, Astrovirus, *Helicobacter pylori*, *Campylobacter jejuni*, *Escherichia coli* O157, *Entamoeba histolytica*, *Giardia lamblia*, *Cryptosporidium parvum*.
- EXTERNAL CENTER:** *C. DIFFICILE* Ag (GDH) MonlabTest was evaluated against different micro-organisms likely to be present in the intestinal tract at any time at concentrations high enough. Tests were carried out on bacterial suspensions at concentration of 10⁸ cfu/mL. The GDH test showed no cross reactivity with any of the micro-organisms listed below:
Aeromonas hydrophila, *Bacillus pumillus*, *Campylobacter jejuni*, *Candida albicans*, *Citrobacter koseri*, *Clostridium perfringens*, *Clostridium sordellii*, *Enterobacter aerogenes*, *Enterobacter cloacae*, *Enterococcus faecalis*, *Enterococcus faecium*, *Escherichia coli* 1, *Escherichia coli* 2, *Klebsiella pneumoniae*, *Klebsiella oxytoca*, *Morganella morganii*, *Proteus mirabilis*, *Pseudomonas aeruginosa*, *Salmonella enteritidis*, *Salmonella typhi*, *Salmonella typhimurium*, *Serratia marcescens*, *Shigella dysenteriae*, *Shigella flexnerii*, *Shigella sonnei*, *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Vibrio cholerae*, *Vibrio parahaemolyticus*, *Yersinia enterocolitica*.

MODIFICATIONS REGARDING THE PREVIOUS VERSION OF THIS DOCUMENT

The structure of the general information provided at the beginning of the instructions is improved, itemising the different points in order to be in line with IVDR Regulation requisites.

Regarding the test's performance parameters, the following sections have been amended: "Precautions", "Samples", "Sample preparation", "Limitations of the procedure", "Detection capability", "Analytical specificity", "Diagnostic sensitivity and specificity", "Repeatability", "Reproducibility", "The Hook or prozone effect", previously "Hook effect", "Interfering substances" and "Cross-reactivity".

REFERENCES

- Rupnik M, Dupuy B, Fairweather NF, Gerding DN, Johnson S, Just I, Lyster DM, Popoff MR, Rood JI, Sonenshein AL, Thelestam M, Wren BW, Wilkins TD and von Eichel-Streiber C. Revised nomenclature of *Clostridium difficile* toxins and associated genes. J Med Microbiol. Feb 2005, Vol: 54(Pt 2) pp: 113-117.
- Vonberg RP, Reichardt C, Behnke M, Schwab F, Zindler S and Gastmeier P. Costs of nosocomial *Clostridium difficile* associated diarrhoea. J Hosp Infect. Sep 2008, Vol: 70(1) pp:15-20.
- Finn E, Andersson FL and Madin-Warburton M. Burden of *Clostridioides difficile* infection (CDI) – a systematic review of the epidemiology of primary and recurrent CDI. BMC Infectious Diseases 2021, Vol: 21 pp: 456.
- Kuehne SA, Cartman ST, Heap JT, Kelly ML, Cockayne A and Minton NP. The role of toxin A and toxin B in *Clostridium difficile* infection. Nature. Oct 2010, Vol: 467(7316) pp: 711-713.
- Shetty N, Wren MW, Coen PG. The role of glutamate dehydrogenase for the detection of *Clostridium difficile* in faecal samples: a meta-analysis. J Hosp Infect. 2011. Jan 77 (1): 1-6.

- Cohen SH, Gerding DN, Johnson S, Kelly CP, Loo VG, McDonald LC, Pepin J, Wilcox MH; Society for Healthcare Epidemiology of America; Infectious Diseases Society of America. Clinical practice guidelines for *Clostridium difficile* infection in adults: 2010 update by the society for healthcare epidemiology of America (SHEA) and the infectious diseases society of America (IDSA). Infect. Control Hosp. Epidemiol. 2010. May; 31 (5): 431-55.
- Wilcox MH, Planché T, Fang FC and Gilligan P. What is the current role of algorithmic approaches for diagnosis of *Clostridium difficile* infection. J Clin Microbiol. 2010 Dec 48 (12): 4347-53
- Wren MW, Sivapalan M, Kinson R and Shetty NR. Laboratory diagnosis of *Clostridium difficile* infection. An evaluation of tests for faecal toxin, glutamate dehydrogenase, lactoferrin and toxigenic culture in the diagnostic laboratory. Br J Biomed Sci. 2009; 66(1): 1-5.

SYMBOLS USED FOR IVD COMPONENTS AND REAGENTS

	Manufactured by		For <i>in vitro</i> diagnostic use
	Do not re-use		Please read pack insert
	Contains sufficient for <n> tests		Store at
	Catalogue number		Expiry date
	Lot number		Precaution

This product fulfils the requirements of Directive 98/79/EC on *in vitro* diagnostic medical devices